## SIGMA-ALDRICH®

#### sigma-aldrich.com

3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

# **Product Information**

### Ribonuclease A Solution from bovine pancreas

Solution of 50% glycerol, for molecular biology

Catalog Number **R4642** Storage Temperature –20 °C

#### CAS RN 9001-99-4

EC 3.1.27.5

Synonyms: Ribonuclease I, Pancreatic ribonuclease, Ribonuclease 3'-pyrimidinooligonucleotidohydrolase, RNase A, Endoribonuclease I

#### **Product Description**

Ribonuclease A (RNase A) is an endoribonuclease that attacks at the 3' phosphate of a pyrimidine nucleotide. The sequence of pG-pG-pC-pA-pG will be cleaved to give pG-pG-pCp and A-pG. The highest activity is exhibited with single-stranded RNA.<sup>1</sup>

RNase A is a single chain polypeptide containing 4 disulfide bridges. In contrast to RNase B, RNase A is not a glycoprotein.<sup>2</sup> RNase A can be inhibited by alkylation of His<sup>12</sup> or His<sup>119</sup>, which are present in the active site of the enzyme.<sup>3</sup> Activators of RNase A include potassium and sodium salts.

Molecular mass:<sup>4</sup> 13.7 kDa (amino acid sequence)

Extinction coefficient:<sup>5</sup> E<sup>1%</sup> = 7.1 (280 nm)

Isoelectric point:<sup>6</sup> pl = 9.6

Optimal temperature: 60 °C (activity range of 15-70 °C)

Optimal pH:7 7.6 (activity range of 6-10)

Inhibitors: ribonuclease inhibitor

The product is supplied in a solution containing 50% glycerol and 10 mM Trizma<sup>®</sup>-HCl, pH 8.0.

Activity: ≥70 Kunitz<sup>8</sup> units/mg protein

A major application for RNase A is the removal of RNA from preparations of plasmid DNA. In this application, the presence of DNase activity as an impurity is a concern. The boiling-water bath method<sup>9</sup> used to eliminate contaminating DNase activity has proven unreliable. For this reason, Sigma-Aldrich developed a proprietary chromatographic preparation method for elimination of DNase activity. This product has been specifically used in studies of nuclear envelope isolation,<sup>10</sup> echinoderm embryos,<sup>11</sup> and molecular cytogenics.<sup>12</sup>

#### **Precautions and Disclaimer**

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

<u>Note</u>: RNase A is stable to both heat and detergents. In addition, it adsorbs strongly to glass. Scrupulous precautions are necessary to ensure that residual RNase A does not cause artifacts in processes that require intact RNA.

#### **Preparation Instructions**

<u>Note</u>: Boiling stock solutions of Catalog No. R4642 to inactivate residual DNase is not necessary, and may cause precipitation of RNase and possible loss of enzymatic activity. If an RNase A solution is heated at a neutral pH, precipitation will occur. When heated at a lower pH, some precipitation may occur because of protein impurities that are present.

#### Storage/Stability

This product remains active for at least 2 years when stored properly at -20 °C.

RNase A is a very stable enzyme and solutions have been reported to withstand temperatures up to 100 °C. At 100 °C, an RNase A solution is most stable between pH 2.0 and 4.5.<sup>13</sup>

#### Procedure

For removal of RNA from preparations of plasmid DNA, DNase-free RNase A is used at a final concentration of 10  $\mu$ g/ml.<sup>14</sup>

#### References

- Sambrook, J., and Russell, D.W., *Molecular Cloning, A Laboratory Manual* (3<sup>rd</sup> ed). Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY), Volume 3, A4.39 (2001).
- Burrell, M.M., "RNase A (EC 3.1.27.5)", in *Methods* in *Molecular Biology*, Vol. 16: *Enzymes of Molecular Biology* (M.M. Burrell, ed.). Humana Press (Totowa, NJ), Ch. 13, pp. 263-270 (1993).
- 3. Plummer, T.H., Jr., and Hirs, C.H.W., *J. Biol. Chem.*, **238(4)**, 1396-1401 (1963).
- Heinrikson, R.L. *et al.*, *J. Biol. Chem.*, **240(7)**, 2921-2934 (1965).
- 5. Smyth, D.G. *et al.*, *J. Biol. Chem.*, **238(1)**, 227-234 (1963).
- Keller, P.J. *et al.*, *J. Biol. Chem.*, **233(2)**, 344-349 (1958).
- Tanford, C., and Hauenstein, J. D., *J. Am. Chem.* Soc., **78(20)**, 5287-5291 (1956).
- 8. Schomberg, D., and Salzmann, M., *Enzyme Handbook*, Vol. 3, 1-3, under EC 3.1.27.5 (1990).
- 9. Kunitz, M., J. Biol. Chem., 164(2), 563-568 (1946).
- Ori, A. *et al.*, "The Use of Targeted Proteomics to Determine the Stoichiometry of Large Macromolecular Assemblies", in *Methods in Cell Biology*, Vol. 122: *Nuclear Pore Complexes and Nucleocytoplasmic Transport – Methods* (V. Doye, ed.). Academic Press (San Diego, CA), Ch. 6, pp. 117-146 (2014).

- Strickland, L. *et al.*, "Light Microscopy of Echinoderm Embryos", in *Methods in Cell Biology*, Vol. 74: *Development of Sea Urchins, Ascidians, and Other Invertebrate Deuterostomes: Experimental Approaches* (C.A. Ettenshohn, G.A Wray, and G.M. Wessel, eds). Elsevier Academic Press (San Diego, CA), Ch. 16, pp. 371-409 (2004).
- Harrison, C.J. et al., "Molecular Cytogenics in Childhood Leukemia", in Methods in Molecular Medicine, Vol. 91: Pediatric Hematology: Methods and Protocols (N.J. Goulden and C.G. Steward, eds.). Humana Press (Totowa, NJ), Ch. 9, pp. 123-137 (2004).
- 13. Crestfield, A.M. *et al.*, *J. Biol. Chem.*, **238(2)**, 618-621 (1963).
- Sambrook, J., and Russell, D.W., *Molecular Cloning, A Laboratory Manual* (3<sup>rd</sup> ed). Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY), Volume 1, 1.78-1.79 (2001).

Trizma is a registered trademark of Sigma-Aldrich Co. LLC.

RBG,MAM,KTA,GCY 08/19-2