

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone (800) 325-5832 (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Nystatin

Product Number **N 3503** Storage Temperature -0 °C

Product Description

CAS Number: 1400-61-9

Melting Point: Gradually decomposes above 160 °C

without melting by 250 °C1

 λ_{max} : 290, 307, and 322 nm (ethanol)

Specific Rotation (at 25 °C)¹: -10° (glacial acetic acid)

+25° (pyridine) +12° (DMF)

-7° (0.1 N HCl in methanol)

Molecular Formula of Nystatin A₁⁻¹: C₄₇H₇₅NO₁₇ Molecular Weight of Nystatin A₁⁻¹: 926.1

Nystatin is a fungistatic and fungicidal polyene antibiotic, which increases the permeability of the cell membrane of sensitive fungi by binding to sterols, chiefly ergosterol. Its main action is against *Candida species*. It is also effective against *Aspergillus*, *Coccidioides immitis*, *Cryptococcus neoformans*, *Histoplasma capsulatum*, *Blastomyces dermatidis*, and other yeasts and fungi. Nystatin has been used to enrich mutants by killing yeast cells. Nystatin has no antibacterial activity. The minimum inhibitory concentration for most sensitive fungi has been reported to range from 1.56 to 6.25 μg/ml.

Nystatin is poorly absorbed from the gastrointestinal tract. It is not absorbed through the skin or mucous membranes when applied topically.

Nystatin is a complex mixture produced by *Streptomyces aureus* and other *Streptomyces* sp., and contains three biologically active components named as A_1 , A_2 and A_3 .

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

Nystatin is soluble at 28 °C in the following: methanol (11.2 mg/ml), ethanol (1.2 mg/ml), carbon tetrachloride

(1.23 mg/ml), chloroform (0.48 mg/ml), benzene (0.28 mg/ml), and ethylene glycol (8.75mg/ml).

Nystatin is soluble in DMSO (5 mg/ml), yielding a clear, bright yellow solution. This product is freely soluble in DMF and formamide.⁴ A 3% suspension in water has pH 6.5-8.0.³

Storage/Stability

Solutions and aqueous suspensions begin to lose activity soon after preparation. Heat, light, and oxygen accelerate decomposition. Aqueous suspensions are stable for 10 minutes when heating to 100 °C at pH 7. Nystatin (in tissue culture media) is stable at 37 °C for three days. It is also stable in moderately alkaline media, but labile at pH 9 and 2. Activity is not diminished by blood or serum. 1,2

Nystatin, cell culture tested (Product No. N 1638) is a sterile nystatin suspension in Dulbecco's Phosphate Buffered Saline. This product is stored frozen and shipped on dry ice. It has a shelf life of 24 months when stored frozen.

References

- 1. The Merck Index, 13th ed., Entry# 6770.
- 2. Brezis, M., et al., Polyene toxicity in renal medulla: injury mediated by transport activity. Science, **224**, 66-68 (1984).
- 3. Martindale The Extra Pharmacopoeia, 29th ed., Reynolds, J. E. F., ed., The Pharmaceutical Press (London, England: 1989), pp. 432-433.
- 4. Clarke's Isolation & Identification of Drugs, 2nd ed., Moffat, A.C., et al., eds., The Pharmaceutical Press (London, England: 1986), p. 829.
- Upson's Handbook of Clinical Veterinary Pharmacology, 2nd ed., Udson, D. W., Veterinary Medicine Publishing Co. (Lenexa, KS: 1985), p. 640.
- 6. Fink, G. R., The Biochemical Genetics of Yeast. Meth. Enzymol., **17-A**, 59-78 (1970).

ALF/RXR 6/06