Instructions for Use

Papanicolaou Stains

Procedure No. HT40

Intended Use

The Sigma-Aldrich Papanicolaou Stains are intended for histochemical staining of exfoliative cells in human cytologic specimens. Papanicolaou staining reagents are for “In Vitro Diagnostic Use.” For professional use only. The data obtained from this manual, qualitative procedure when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of vaginal, uterine and cervical conditions.

Papanicolaou staining techniques, reviewed in a concise report by Street,¹ have changed little in the past 40 years. The stain is used for examining exfoliative cells of sputum as well as vaginal, cervical, and other body secretions. In general, cells are fixed to a slide, treated with a hematoxylin nuclear stain and counterstained with a mixture of orange G, eosin Y and fast green FCF (a replacement for light green SF yellowish). These treatments impart characteristic color to nuclei and cytoplasmic components.²

Reagents

- Papanicolaou Stain OG 6** (Cat. No. HT40116; HT40132; HT40180; HT401128)
Certified orange G, 0.3% w/v, phosphotungstic acid, 0.015% w/v, in denatured alcohol.
- Papanicolaou Stain Modified EA 65** (Cat. No. HT40232)
Certified eosin Y, 0.25% w/v, certified fast green FCF, 0.004% w/v, phosphotungstic acid, 0.4% w/v, in denatured alcohol with stabilizer
- Papanicolaou Stain EA 50** (Cat. No. HT40316; HT40332; HT403128)
Certified eosin Y, 0.23% w/v, certified fast green FCF, 0.08% w/v, certified bismarck brown, 0.05%, phosphotungstic acid, 0.2% w/v, in denatured alcohol
- Papanicolaou Stain EA 65** (Cat. No. HT40432)
Certified eosin Y, 0.23% w/v, certified fast green FCF, 0.01% w/v, certified bismarck brown, 0.05%, phosphotungstic acid, 0.2% w/v, in denatured alcohol

Special Materials Required but Not Provided

- Differentiation Solution (Cat. No. A3179 or A3429)
- Reagent Alcohol (Cat. No. R8382) or Ethanol (absolute)
- Hematoxylin Nuclear Stain:
 - Harris Hematoxylin Solution, Modified (Cat. No. HHS16; HHS32; HHS80; HHS128)
 - Hematoxylin Solution, Gill No. 1 (Cat. No. GHS116; GHS132; GHS1128)
 - Hematoxylin Solution, Gill No. 2 (Cat. No. GHS216; GHS232; GHS280; GHS2128)
- Scott's Tap Water Substitute Concentrate (Cat. No. S5134)
Note: Prepare Scott's Tap Water Substitute working solution as indicated on bottle label.
- Xylene or xylene substitute.

Storage and Stability

Store reagents at room temperature (18-26°C). Reagents are stable until expiration dates shown on labels.

Preparation

Differentiation Solution is ready to use.

Scott's Tap Water Substitute working solution is prepared by mixing 1 part of Scott's Tap Water Substitute Concentrate with 9 parts deionized water.

Filter Papanicolaou reagents before use.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Smears are prepared and fixed immediately with a spray fixative or placed in 95% ethanol for a minimum of 15 minutes.

Notes⁴

- Separate staining solutions should be maintained for gynecologic smears and nongynecologic smears.
- All solutions preceding the hematoxylin solutions should be filtered after each staining sequence and all solutions should be filtered daily.
- Cellulose membrane filters should be dehydrated (Procedure 1, Step 14 and 15) in isopropyl alcohol to prevent dissolution.
- Staining times may be varied for individual color preferences.
- Ethanol may be substituted for Reagent Alcohol.
- Appropriate control slides should be included in each run.

Procedure

Procedure 1: Standard Method

- Reagent Alcohol, 95% for 15 minutes.
- Tap water, rinse.
- Hematoxylin solution for 1-3 minutes.
NOTE: Time will vary with selection of hematoxylin solution.
- Tap water, rinse.
- Differentiation Solution for 20-60 seconds.
NOTE: Necessary for regressive hematoxylin stains only.
- Tap water, rinse.
- Scott's Tap Water Substitute working solution, 10 dips.
- Tap water, rinse.
- Reagent Alcohol, 95%, 10 dips.
- Papanicolaou Stain OG 6 for 1.5 minutes.
- Reagent Alcohol, 95%, 10 dips.
- Papanicolaou Stain Modified EA 65, OR Papanicolaou Stain EA 50, OR Papanicolaou Stain EA 65 for 2.5 minutes.
- Reagent Alcohol, 95%, two changes, 10 dips each.
- Reagent Alcohol, 100% for 1 minute.
- Xylene or xylene substitute, two changes for 2 minutes.
- Coverslip and examine microscopically.

Procedure 2: Rapid Papanicolaou Stain

Reagent Preparation: Prepare Papanicolaou Stain OG 6-EA Solution by mixing equal volumes of Papanicolaou Stain OG 6 and any one of the Papanicolaou EA stains.

- Wash smears in running water for 1 minute.
NOTES:
This step removes spray-fixative and the timing is for use with water soluble fixatives. Initial washing is important and must be thorough to expose the cell to staining.
If slides have been fixed in 95% ethanol, timing may be shortened to a few seconds.
- Stain in Hematoxylin Solution, Gill No. 2 for 1-3 minutes.
Wash in running tap water.
- Scott's Tap Water Substitute working solution for 20 seconds.
Wash in running tap water.
- Reagent Alcohol, 95% for 20 seconds.
- Stain in Papanicolaou Stain OG 6-EA (prepared above) for 1-3 minutes.
- Rinse in two changes of Reagent Alcohol, 100%.
- Clear in two changes of xylene substitute, two changes for 2 minutes each.
- Coverslip and examine microscopically.

Performance Characteristics

Nuclei are stained blue while cytoplasm displays varying shades of blue, orange, pink or red.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT401	Papanicolaou Stain, OG-6	Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
HT402	Papanicolaou Stain, Modified EA 65	Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
HT403	Papanicolaou Stain EA 50	Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
HT404	Papanicolaou Stain EA 65	Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT40116, HT40132, HT40180, HT401128, HT40232, HT40316, HT40332, HT403128, HT40432:

- H225: Highly flammable liquid and vapor.
- H319: Causes serious eye irritation.
- H371: May cause damage to organs (Eyes, Central nervous system).
- P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- P233: Keep container tightly closed.
- P240: Ground and bond container and receiving equipment.
- P241: Use explosion-proof electrical/ ventilating/ lighting/ equipment.
- P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P308 + P311: IF exposed or concerned: Call a POISON CENTER/ doctor.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

- Street CM: Papanicolaou Techniques in Exfoliative Cytology. IN Laboratory Technique in Biology and Medicine, 3rd ed. EV Cowdry Editor, Williams & Wilkins, Baltimore, 1952, p 253
- Johnson PL, Klein MN: Application of Papanicolaou stain to paraffin sections. Stain Technol 31:223, 1956
- Liu W: A simplified cytologic staining technic. Am J Clin Pathol 54:767, 1970
- The Manual of Cytotechnology, edited by CM Keebler, TM Somrak, 7th Ed., American Society of Clinical Pathologists Press, Chicago, 1993, p. 442-444

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 7.0	2021
Rev. 8.0	2022
Rev. 9.0	2022
	Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards.
Rev. 10.0	2025
	Revised Intended Use statement. Updated naming convention. Added CH-REP information

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

1. Street CM: Papanicolaou Techniques in Exfoliative Cytology. IN Laboratory Technique in Biology and Medicine, 3rd ed. EV Cowdry Editor, Williams & Wilkins, Baltimore, 1952, p 253
2. Johnson PL, Klein MN: Application of Papanicolaou stain to paraffin sections. Stain Technol 31:223, 1956
3. Liu W: A simplified cytologic staining technic. Am J Clin Pathol 54:767, 1970
4. The Manual of Cytotechnology, edited by CM Keebler, TM Somrak, 7th Ed., American Society of Clinical Pathologists Press, Chicago, 1993, p. 442–444

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Références

- Street C.M.: Papanicolaou Techniques in Exfoliative Cytology. IN Laboratory Technique in Biology and Medicine, 3rd ed. E.V. Cowdry Editor, Williams & Wilkins, Baltimore, 1952, p 253
- Johnson P.L., Klein M.N.: Application of Papanicolaou stain to paraffin sections. Stain Technol 31:223, 1956
- Liu W.: A simplified cytologic staining technic. Am J. Clin Pathol 54:767, 1970
- The Manual of Cytotechnology, edited by C.M. Keebler, T.M. Somrak, 7th Ed., American Society of Clinical Pathologists Press, Chicago, 1993, p. 442-444

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Bibliografia

1. Street CM: Papanicolaou Techniques in Exfoliative Cytology. IN Laboratory Technique in Biology and Medicine, 3rd ed. EV Cowdry Editor, Williams & Wilkins, Baltimore, 1952, p 253
2. Johnson PL, Klein MN: Application of Papanicolaou stain to paraffin sections. Stain Technol 31:223; 1956
3. Liu W: A simplified cytologic staining technic. Am J Clin Pathol 54:767, 1970
4. The Manual of Cytotechnology, edited by CM Keebler, TM Somrak, 7th Ed., American Society of Clinical Pathologists Press, Chicago, 1993, p. 442-444

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Referencias

1. Street CM: Papanicolaou Techniques in Exfoliative Cytology. IN Laboratory Technique in Biology and Medicine, 3rd ed. EV Cowdry Editor, Williams & Wilkins, Baltimore, 1952, p 253
2. Johnson PL, Klein MN: Application of Papanicolaou stain to paraffin sections. Stain Technol 31:223, 1956
3. Liu W: A simplified cytologic staining technic. Am J Clin Pathol 54:767, 1970
4. The Manual of Cytotechnology, edited by CM Keebler, TM Somrak, 7th Ed., American Society of Clinical Pathologists Press, Chicago, 1993, p. 442-444

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Referenser

- Street CM: Papanicolaou Techniques in Exfoliative Cytology. IN Laboratory Technique in Biology and Medicine, 3:e utg. EV Cowdry redaktör, Williams & Wilkins, Baltimore, 1952, s. 253
- Johnson PL, Klein MN: Application of Papanicolaou stain to paraffin sections. Stain Technol 31:223, 1956
- Liu W: A simplified cytologic staining technic. Am J Clin Pathol 54:767, 1970
- The Manual of Cytotechnology, redigerad av CM Keebler, TM Somrak, 7:e utg, American Society of Clinical Pathologists Press, Chicago, 1993, s. 442-444

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- Street CM: Papanicolaou Techniques in Exfoliative Cytology. IN Laboratory Technique in Biology and Medicine, 3rd ed. EV Cowdry Editor, Williams & Wilkins, Baltimore, 1952, σελ. 253
- Johnson PL, Klein MN: Application of Papanicolaou stain to paraffin sections. Stain Technol 31:223, 1956
- Liu W: A simplified cytologic staining technic. Am J Clin Pathol 54:767, 1970
- The Manual of Cytotechnology, edited by CM Keebler, TM Somrak, 7th Ed., American Society of Clinical Pathologists Press, Chicago, 1993, σελ. 442-444

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Referanser

- Street CM: Papanicolaou Techniques in Exfoliative Cytology. IN Laboratory Technique in Biology and Medicine, 3. utg. EV Cowdry Editor, Williams & Wilkins, Baltimore, 1952, s. 253
- Johnson PL, Klein MN: Application of Papanicolaou stain to paraffin sections. Stain Technol 31:223, 1956
- Liu W: A simplified cytologic staining technic. Am J Clin Pathol 54:767, 1970.
- The Manual of Cytotechnology, edited by CM Keebler, TM Somrak, 7. utg., American Society of Clinical Pathologists Press, Chicago, 1993, s. 442-444

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Referenties

- Street CM: Papanicolaou Techniques in Exfoliative Cytology. IN Laboratory Technique in Biology and Medicine, 3e ed. EV Cowdry uitgever, Williams & Wilkins, Baltimore, 1952, pag. 253
- Johnson PL, Klein MN: Application of Papanicolaou stain to paraffin sections. Stain Technol 31:223, 1956
- Liu W: A simplified cytologic staining technic. Am J Clin Pathol 54:767, 1970.
- The Manual of Cytochemistry, bewerkt door CM Keebler, TM Somrak, 7e editie, American Society of Clinical Pathologists Press, Chicago, 1993, pag. 442-444.

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MDSS CH GmbH
MDSS GmbH
Chemie GmbH