SIGMA-ALDRICH®

sigma-aldrich.com

3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Product Information

Butyrylcholinesterase human, recombinant expressed in transgenic goats

Catalog Number **B7061** Storage Temperature –70 °C

CAS RN 9001-08-5 EC 3.1.1.8 Synonyms: rBChE; Acylcholine acylhydrolase; Pseudocholinesterase; cholinesterase

Product Description

Butyrylcholinesterase (BChE) belongs to the same structural class of proteins, the esterase/lipase family, as acetylcholinesterase (AchE, EC 3.1.1.7). They are serine hydrolases that share substantial structural similarities, but differ in substrate specificities and inhibitor sensitivities.¹ BChE can, unlike AChE, efficiently hydrolyze larger esters of choline such as butyrylcholine and benzoylcholine.

 $\begin{array}{rl} & \text{BChE} \\ \text{Butyrylcholine} & \rightarrow & \text{butyric acid + choline} \\ & \text{H}_2\text{O} \end{array}$

Although BChE is found in the serum, hemopoietic cells, liver, lung, heart, and the central nervous system of vertebrates, it has no known physiological function.^{2,3}

Molecular mass:⁴ 170 kDa (dimer)

Extinction coefficient:⁴ $E^{1\%}$ = 19.8 (280 nm)

Isoelectric point (pl):4 5.0

pH Optimal activity range:⁴ 7.0-8.0

pH Activity range:⁴ 6.0–10.0

Inhibitors:⁵ Pyridostigmine, betaine, Tetraisopropyl pyrophosphoramide, butyrylthiocholine, Diisopropyl phosphorofluoridate, donepezil, Eserine, Physostigmine, Parathion

Substrates: Butyrylthiocholine⁵ Butyrylcholine⁶

Butyrylcholine⁶ Acetylcholine⁷ Acetylthiocholine⁸ Cocaine⁹

Selective inhibition of BChE activity can be used in the detection of organophosphates.¹⁰ Its use in the treatment of organophosphate toxicity shows promise and there is a correlation between the level of BChE in human blood and degree of protection against potentially toxic nerve agents.³

There has also been interest in the roles of cholinesterases with regard to Alzheimer's disease. Investigations into selective inhibitors may provide a clearer picture of the physiological role of BChE in both healthy and diseased individuals.²

This recombinant, human product is isolated from the milk of transgenic goats. It is purified by affinity chromatography and is primarily in the dimeric form, with the remainder being monomers. The product is supplied as a solution in 10 mM sodium phosphate, pH 7.2, with 50 mM NaCl, 1 mM EDTA, and 10% glycerol.

Specific activity: ≥500 units/mg protein

Unit definition: One unit will hydrolyze 1.0 μ mole of butyrylcholine to choline and butyrate per minute at pH 8 at 37 °C. The activity obtained using butyrylcholine as a substrate is ~2.5 times that obtained using acetylcholine.

BChE is assayed titrimetrically in a 50.4 ml reaction mixture containing 4 mM butyrylcholine, 40 mM MgCl₂, 100 mM NaCl, and 12–24 units of BChE at pH 8 and 37 $^\circ$ C.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Store the product at -70 °C. When stored at -70 °C, the enzyme retains activity for at least 2 years.⁴

Product loses activity at elevated temperatures, below pH 6.0. 4

References

- Kovarik, Z., *et al.*, Acetylcholinesterase active centre and gorge conformations analysed by combinatorial mutations and enantiomeric phosphonates. Biochem. J., **373**, 33-40 (2003).
- Savini, L., *et al.*, Specific targeting of acetylcholinesterase and butyrylcholinesterase recognition sites. Rational design of novel, selective, and highly potent cholinesterase inhibitors. J. Med. Chem., **46**, 1-4 (2003).

- Blong, R.M., *et al.*, Tetramerization domain of human butyrylcholinesterase is at the C-terminus. Biochem. J., **327**, 747-57 (1997).
- 4. Supplier Data
- 5. Masson, P., et al., Biochim. Biophys. Acta, **1544**, 166-176 (2001).
- Garcia-Ayllon, M.S., et al., Neuroscience., **107**, 199-208 (2001).
- 7. Giacobini, E., Neurochem., 28, 515-522 (2003).
- 8. Platteborze, P.L., Biotechnol. Appl. Biochem., **31**, 225-229 (2000).
- 9. Zhan, C.G., et al., J. Am. Chem. Soc., **125**, 2462-2474 (2003).
- Beattie, B.D., *et al.*, Determination of butyrylcholinesterase inhibition using ion transfer across the interface between two immiscible liquids. Anal. Chem., **66**, 52-57 (1994).

CS,RBG,MAM 05/08-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.