

BACITRACIN Sigma Prod. No. B0125


CAS NUMBER: 1405-87-4

SYNONYMS: Ayfivin; Baciguent; Baci-Jel; Baciliquin;Bacitek Ointment; Fortracin; Parentracin; Penitracin; Topitracin; USAF CB-7; Zutracin

PHYSICAL AND CHEMICAL PROPERTIES:

Appearance: Powder, white to yellow with a tan cast Molecular Formula: $C_{66}H_{103}N_{17}O_{16}S$ (Bacitracin A) Molecular Weight: 1422.7 (Bacitracin A) Melting Point: 221-225°C¹ Isoelectric point: pI = 8.8 (Messing's method); 8.5 (electrophoresis)² pH of 1% solution = 6-7³ $E^{mM}(225nm) = 6.25^{6}$ $E^{mM}(252nm) = 2.50^{6}$ $\lambda(exc) = 292nm^{2}$ $\lambda(emiss) = 325nm^{2}$

ProductInformation



Bacitracin A

FOR SPECIFICATIONS SEE CATALOG

PRODUCT DESCRIPTION:

Commercial Bacitracin is a mixture of at least 9 Bacitracins, of which Bacitracin A is the major component.⁴ Bacitracin is a polypeptide complex produced by *Bacillus subtilis* and *Bacillus licheniformis*, and is used as an antibacterial agent, primarily against gram-positive organisms.^{4,5} Bacitracin inhibits bacterial cell wall synthesis by inhibiting dephosphorylation of lipid pyrophosphate.¹⁵ Excellent review articles pertaining to Bacitracin (including its mechanism of action) have been published.^{2,16,17,18,19}

METHODS OF PREPARATION AND PURIFICATION:

B-0125 is manufactured by fermentation. No bovine tissue or bovine tissue extract is used. Several methods of preparation (not necessarily those of Sigma's supplier) for Bacitracin and its fragments have been published in the literature.^{7,8,9,10,11,12}

STABILITY / STORAGE AS SUPPLIED:

This hygroscopic product should be stored desiccated at 2-8°C. The stability of Bacitracin at various temperatures has been discussed in the literature.²

SOLUBILITY / SOLUTION STABILITY:

Sigma tests the solubility of B-0125 in water at 50 mg/mL and obtains clear to slightly hazy, yellow to yellow-tan solutions.

The solubility of Bacitracin has been discussed in the literature.^{2,13} Bacitracin is very soluble in water and methanol; soluble in ethanol; slightly soluble in acetone, benzene and ether; and practically insoluble in chloroform, ether and acetone.^{4,14} Aqueous solutions degrade rapidly at room temperature.³ Bacitracin is relatively stable in acidic solutions and unstable above pH 9.¹⁴ An observed loss in potency is probably due to transformation of Bacitracin A into Bacitracin F, which has a low antimicrobial activity.⁴ Sigma's supplier has provided the following solution stability data:

% Activity Retained After Number of Hours Indicated				
<u>pH</u>	23 Hours	<u>71 Hours</u>	<u>120 Hours</u>	<u>143 Hours</u>
2.4	58	37	18	N/A
4.4	94	88	70	N/A
5.9	90	73	N/A	45
7.9	88	56	N/A	31
9.0	80	27	27	N/A

REFERENCES:

- 1. Aldrich Chemical Co. Catalog, p. 115 (1990-1991).
- 2. G. A. Brewer, Anal. Profiles Drug Subst., 9, 1, (1980).
- 3. *Martindale The Extra Pharmacopeia*, 29th ed., J. E. F. Reynolds, Ed, p. 128, The Pharmaceutical Press, London (1989).
- 4. Merck Index, 12th ed., S. Budavari, Ed., pp. 160-161, # 965 (1996).
- 5. Supplier data.
- 6. *Biochem.*, 10, 2429 (1971).
- 7. E. Munekata, et al., Bull. Chem. Soc. Japan., 46(10), 3187 (1973).
- 8. E. Munekata, et al., *Bull. Chem. Soc. Japan.*, 46(12), 3835 (1973).
- 9. J. Porath, Acta Chem. Scand., 6, 1237 (1952).
- 10. L. Chaiet and T. J. Cochrane, Jr., U.S. Pat. 2,915,432, Merck & Co.(1959).
- 11. B. A. Johnson and F. L. Meleney, U.S. Pat. 2,498,165, U. S. Secy. of War (1950).
- 12. T. E. Freaney and L. P. Allen, U.S. Pat. 2,828,246, Commercial Solvents (1958).
- 13. P. J. Weiss, et al., *Antibiot. Chemother.*, 7, 374 (1957).
- 14. Data for Biochemical Research, 3rd ed., R. M. C. Dawson, et. al., Eds., pp. 318-319, Clarendon Press, Oxford (1993).
- 15. G. Siewert and J. L. Strominger, *PNAS*, 57, 767 (1967).
- 16. D. R. Storm, N. Y. Acad. Sci., 235, 387 (1974).
- L. C. Craig, et. al., Ciba Foundation Symposium on Amino Acids and Peptides with Antimetabolic Activity, G. E. W. Wolstenholme and C.M. O'Connor, Eds., pp. 226-246, Little & Brown, Boston (1958).
- 18. E. D. Weinberg, *Antibiotics*, D. Gottlieb and P. D. Shaw, Eds., pp. 90-99 and 240-245, Springer Verlag, New York (1967).
- 19. D. R. Storm and W. A. Toscano, Jr., *Antibiotics*, Vol. 5, Part 1, pp. 1-17, F. E. Hahn, Ed., Springer Verlag, New York (1979).

Sigma warrants that its products conform to the information contain in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.