

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Nerve Growth Factor-7S from mouse submaxillary glands

Catalog Number **N0513** Storage Temperature –20 °C

Synonym: NGF-7S

Product Description

Nerve Growth Factor (NGF) was first discovered in 1953 by Levi-Montalcini, Hamburger, and Cohen¹⁻³ in two mouse sarcomas, and described as a diffusable agent which strongly promotes fiber outgrowth of sensory neurons in chick embryos. Cohen purified NGF from snake venom⁴ and from mouse salivary glands.⁵

NGF is a neurotrophic agent thought to be provided by peripheral tissues for the guidance and sustenance of outgrowing embryonic sympathetic and sensory neurons.⁶ NGF induces the formation of neurite-like filaments from chick embryo dorsal root ganglia² and from rat PC12 pheochromocytoma cells.⁷ *In vivo* NGF may be involved in fetal development^{8,9} and nerve regeneration.¹⁰ NGF may also play a physiological role within the central nervous system.^{8,11,12}

Cellular receptors for NGF have been found in a variety of cell lines¹³ and tissues, including cholinergic neurons of the brain^{14,15} and Schwann cells of damaged nerve axons.¹⁰ Two kinetic types of NGF receptors have been identified from peripheral neurons,¹⁶ neuroblastoma cells,¹⁷ and PC12 cells¹⁸ and are designated as type I (high affinity) and type II (low affinity). The signal transduction mechanism of the receptor has not been clearly identified.

This product (NGF-7S, Catalog Number N0513) is purified from mouse submaxillary glands under nondissociative conditions. Nerve Growth Factor isolated under such conditions has a sedimentation coefficient of 7.1 S.^{19,20} It is generally believed that NGF-7S is a 130 kDa protein composed of 5 non-covalently linked subunits (2α , 1 β , 2γ), although there is recent evidence for a different endogenous form of high molecular weight NGF.²¹ After dissociation of purified NGF-7S by acidic or basic pH, only the β subunit of NGF (NGF- β , Catalog Number N2393) has neurotrophic activity.⁶ NGF- β is a 26.5 kDa dimer of identical 118-residue chains held together tightly by noncovalent bonds. NGF-7S is a form of NGF- β , initially isolated under dissociative conditions, but is often slightly different from NGF- β , due to proteolysis during the purification of NGF- β .⁶ Apparently the 7S complex protects the amino- and carboxy-terminals of NGF- β from hydrolytic enzymes present in the submaxillary gland extract.²²

This product is supplied lyophilized from a 0.2 μm filtered solution in 2.5 mM sodium phosphate, pH 7.0.

The Nerve Growth Factor-7S product has been tested for bioactivity using PC-12 cells.²³ EC₅₀: 2-250 ng/mL

Endotoxin: ≤10 EU/vial

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

To prepare a stock solution, reconstitute the contents of the vial in a solution of buffered saline or tissue culture medium containing 0.1–1.0% BSA or 1-10% serum. This may be diluted immediately before use to the final working concentration of NGF-7S, generally 0.1–10 ng/ml. Additional filtration of the stock solution is not recommended and may result in product loss due to adsorption onto the filter membrane.

Storage/Stability

Prior to reconstitution, store vial at -20 °C. After reconstitution, the product may be stored in a plastic vial for up to two weeks at 2–8 °C or may be stored as aliquots at -20 °C. Prolonged storage of product or repeated freezing and thawing is not recommended.

References

- 1. Levi-Montalcini, R., Science, 237, 1154 (1987).
- Levi-Montalcini, R., *et al.*, Cancer Res., **14**, 49 (1954).
- Cohen, S., *et al.*, Proc. Natl. Acad. Sci. USA, **40**, 1014 (1954).
- 4. Cohen, S., J. Biol. Chem., 234, 1129 (1959).
- Cohen, S., Proc. Natl. Aced. Sci. USA, 46, 302 (1960).
- Server, A., and Shooter, E., Adv. Protein Chem., 31, 339 (1977).
- Greene, R., and Tischler, A., Proc. Natl. Acad. Sci. USA, 73, 2424 (1976).
- Ayer-Lelievre, C., *et al.*, Med. Biol., **61**, 296 (1983).
- Taniuchi, M., *et al.*, Proc. Natl. Acad. Sci. USA, **83**, 4094 (1986).

- Thornburn, G., *et al.*, Growth and Maturation Factors, Vol. 3 (G. Guroff, ed.) John Wiley & Sons, (NY, NY: 1985) p.175.
- 11. Ebendal, T., Prog. Growth Factor Res., **1**, 143 (1989).
- 12. Dreyfus, C., Trends Pharmacol. Sci., **10**, 145 (1989).
- 13. Eveleth, D., In Vitro Cell. Dev. Biol., **24**, 1148 (1988).
- 14. Hefti, F., et al., Neurosci. Let., 69, 37 (1986).
- 15. Raivich, G., and Kreutzberg, G., Neuroscience, **20**, 23 (1987).
- 16. Godfrey, E., and Shooter, E., J. Neurosci., **6**, 2543 (1986).
- 17. Marchetti, D., and Perez-Polo, J., J. Neurochem., **49**, 475 (1987).
- Buxser, S., *et al.*, J. Biol. Chem., **265**, 12701 (1990).
- 19. Varon, S., et al., Biochemistry, 6, 2202 (1967).
- 20. Varon, S., et al., Biochemistry, 7, 1296 (1968).
- 21. Young, M., et al., Biochemistry, 27, 6675 (1988).
- 22. Moore, J., et al., Biochemistry, 13, 833 (1974).
- 23. Green, L., J. Cell Biol., 78, 747 (1928).

KAA, MAM 09/06-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.