

69385 Protein Standard Mix 15 – 600 kDa for testing of SEC/GFC columns

Storage Temperature -20°C

Product Description

The Protein Standard Mix is a calibration standard to test and monitor performance of size exclusion chromatography (SEC) columns. It is a lyophilized mixture of molecular weight markers ranging from 15 kDa to 600 kDa. The mixture consists of 4 proteins thyroglobulin, gamma-globulins, albumin, ribonuclease A and a low molecular weight marker (pABA).

The standard is supplied as a single product or as set of 6 x 2 mL vials.

Components

Thyroglobulin bovine	0.5 g/l	MW ~ 670 000
γ-globulins from bovine blood	1.0 g/l	MW ~ 150 000
Albumin chicken egg grade VI	1.0 g/l	MW ~ 44 300
Ribonuclease A type I-A from		
bovine pancreas	1.0 g/l	MW ~ 13 700
<i>p</i> -aminobenzoic acid (pABA)	0.01 g/l	
Buffers and preservatives		


Preparation Instructions

Add 1 mL deionized H2O to the protein standard vial and swirl gently to solve the lyophilizate. Spin down eventually unsolved fine particles. The clear supernatant can then be applied to the HPLC column.

Storage/Stability

The lyophilized standard is stable when stored at -20° C for 2 years. The hydrated protein mixture can be kept for up to one week at 2-8°C.

Figure





Conditions:

Column 7.8 mm ID x 30 cm L

Elution 0.1 M sodium phosphate, 0.1 M sodium sulfate,

0.05 % sodium azide, pH 6.7

Flow rate 1 ml/min
Detection UV @ 280 nm

Recommended injection volume 20 µl

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

The vibrant M, Supelco, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources. © 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

