

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Poly-D-lysine hydrobromide Cell Culture Tested

Product Number **P 6407** Storage Temperature -0 °C

Product Description

CAS Number: 27964-99-4 Molecular weight range: 70 - 150 kDa

This product is γ -irradiated and is tested for use in cell culture.

Poly-D-Lysine is a charge enhancer, therefore, it can be used for coating many surfaces. When used to coat ELISA plates, it tends to improve binding of IgG. High pH will aid binding to the plates.

Certain cells can digest poly-L-lysine, in which case the poly-D-lysine should be used. The lower molecular weight poly-D-lysine (30-70 kDa) is easier to use because it is less viscous in solution, but the higher molecular weight poly-D-lysine (>300 kDa) provides more attachment sites per molecule. The mid-range product has a MW range of 70-150 kDa.

There might be a small amount of product in the β structure, but very little. This is because the HBr interferes with hydrogen bonding between the amino and carboxyl groups, or between the amino group and other N or O containing moieties. The secondary structures of polyamino acids have never been extensively studied, but it is known that the degree of secondary conformational regularity depends on the amino acid in the polymer.

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in water (50 mg/ml).

Storage/Stability

Transition of a high molecular weight poly-D-lysine (degree of polymerization 1500) in aqueous solution from a helical to a random coiled conformation occurs under the influence of decreasing pH or increasing temperature.¹

Procedure

For coating slides with poly-D-lysine, in general, it is recommended to use a 0.01% (w/v) solution. After a 5 minute incubation, the excess solution is removed and the slides are dried at room temperature or in an oven at gentle heat. Coated slides should be stable at room temperature or refrigerated (2-8 °C for one year if protected from dust. There is no need to sterilize slides coated with poly-D-lysine by autoclaving; simply expose the slides to UV light overnight.

If uneven coating occurs, the glass slides may be pretreated with 1 mM magnesium acetate for 2-3 hours and then rinsed well before coating. Alternatively, they may be acid washed (hydrochloric acid or sulfuric acid). This treatment should allow for an even coating with the poly-D-lysine solution.

This product can also be used to coat tissue culture plasticware. A solution of 0.1 mg/ml is prepared and 0.5 ml is aseptically transferred to a non-tissue culture treated 25 cm² flask. The flask is gently rocked to evenly coat the surface. After five minutes, the excess solution is removed and the surface is thoroughly rinsed with sterile water and allowed to dry for several hours.

References

1. The Merck Index, 11th ed., Entry# 7549.

CMH/ARO 9/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.