

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com

sigma-aldrich.com

ProductInformation

STREPTOZOTOCIN MIXED ANOMERS Sigma No. S0130

CAS NUMBER: 18883-66-4 SYNONYM: Streptozocin

PHYSICAL DESCRIPTION:

Appearance: White to yellow powder

Molecular formula: C₈H₁₅N₃O₇ Molecular weight: 265.2

Melting point: Decomposes at 115EC if anhydrous. Although Sigma does not determine a value, water

content should be #3%.² $E^{mM}(228nm) = 6.36$ (ethanol)

Optical rotation: +39E(equilibrium of α, ß anomers in H₂O, 25EC)¹

Please consult the Material Safety Data Sheet about the properties of this material as a potential carcinogen, mutagen and toxic chemical.

STORAGE / STABILITY AS SUPPLIED:

If the product is stored frozen and protected from moisture and air, it is stable for approximately 2 years. (After 12 year, a sample changed from 94.9% α -anomer to 94.7%, as measured by HPLC.)²

SOLUBILITY / SOLUTION STABILITY:

Streptozotocin is soluble in water, the lower alcohols and in ketones. This product dissolves in water at 50 mg/mL to give a light yellow solution, from clear to slightly hazy. Aqueous solutions rapidly undergo mutarotation to an equilibrium mixture of alpha- and beta-anomers.

STREPTOZOTOCIN MIXED ANOMERS Sigma No. S0130

SOLUBILITY / SOLUTION STABILITY:

Maximum solution stability is at pH 4, with stability decreasing rapidly at higher or lower pH. Freshly prepared solutions are clear and have a light straw color. On standing, they take on a yellow to brown color and effervesce, indicating decomposition. ^{2,3} Solutions should be prepared just before use, since the product is unstable.

GENERAL REMARKS:

This product is an antineoplastic antibiotic produced by the growth of a Streptomyces achromogenes variant or by synthesis. It may affect glucose metabolism. It is used mainly in the treatment of pancreatic (isletcell) tumors.⁴ Burcelin et al. used intravenous injection of streptozotocin in rats at a dose of 65 mg/kg body weight to induce diabetes (using cold 0.1 M citrate buffer pH 4.5).⁵ In rats and dogs, diabetes was induced using intravenous dosage of 50 mg/kg (using 1-2% w/v solutions in saline buffered with citrate dextrose solution at pH 5.0).³ It has been used for the treatment of malignant insulinoma; very precise assays for this drug have been developed.⁶

Streptozotocin does not cross the blood-brain barrier, but its metabolites are found in cerebral spinal fluid.⁴ Its biological half-life in cell culture medium was shown to be approximately 19 minutes.⁷

The antileukemic effects of streptozotocin and its analogs have been reported. Streptozotocin has been shown to be a potent methylating agent that reacts with DNA *in vitro* to form methylated purines. A review article addressed a number of antineoplastic antibiotics, including streptozotocin. A useful handbook offered several references for use in animal studies.

REFERENCES:

- 1. *Merck Index*, 12th ed., #8991 (1996).
- 2. Sigma quality control or supplier information.
- 3. Rakieten, N. et al., Cancer Chemother. Reports, No. 29, 91-98 (1963).
- 4. *Martindale: The Extra Pharmacopoeia*, 29th ed., (Pharmaceutical Press, 1989), p. 649.
- 5. Burcelin, R. et al., *Biochem. J.*, 291, 109-113 (1993).
- 6. Oles, P.J., J. Pharmaceutical Sci., 67 (9), 1300 (1978).
- 7. Jensen, E.M. et al., *J. Natl. Cancer Inst.*, 59, 941-944 (1977).
- 8. Bhuyan, B.K. et al., Cancer Chemother. Reports, Part 1, 45 (6), 709-720 (1972).
- 9. Bennett, R.A. and Pegg, A.E., *Cancer Research*, 41, 2756-2790 (1981).
- 10. Cheng, C.C. and K.-Y. Zee-Cheng, *J. Pharmaceutical Sci.*, 61, 485-501 (1972).
- 11. *Drug Dosage in Laboratory Animals: A Handbook*, 3rd ed., R.E. Borchard et al., eds. (CRC Press, 1992).

Sigma warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.