

HYPOXANTHINE-THYMIDINE (HT) MEDIA SUPPLEMENT [50X]

Product Number **H0137** Storage Temperature -20°C

Product Description

The production of monoclonal antibodies involves the fusion of myeloma cells with spleen cells and the selection of hybridomas (cell hybrids) in culture utilizing hypoxanthine-aminopterin-thymidine (HAT) supplemented medium. Aminopterin blocks the synthesis of DNA by inhibiting dihydrofolate reductase ³. Cells that lack the ability to utilize the salvage pathway for nucleotide synthesis are eliminated. Cells that possess hypoxanthine-guanine phosphoribosyl transferase (HPRTase) and thymidine kinase (TK) enzymes can utilize the salvage pathway if supplied with hypoxanthine and thymidine. ^{1,2}

The purpose of HAT medium [Prod. No. H0262] is to: (1) selectively kill unfused myeloma cells that are well adapted to tissue culture and would otherwise outgrow any hybridomas produced and (2) eliminate any myeloma-myeloma hybridomas that lack HPRTase. HPRTase positive spleen-spleen hybridomas, although not sensitive to aminopterin, are normally short-lived in culture.

After selection is complete (approximately 10-14 weeks), aminopterin is diluted from the culture by several passages of the cells in hypoxanthine-thymidine (HT) supplemented medium (approximately 2-3 weeks) before transfer into normal hybridoma growth medium.

Components

Each vial contains: Hypoxanthine 6.8 mg Thymidine 1.94 mg

ProductInformation

Precautions and Disclaimer

REAGENT: For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

HT media supplement is supplied as a lyophilized, γ -irradiated powder for use in aseptic procedures. When reconstituted to 10 ml, each vial contains 5 x 10 3 M hypoxanthine and 8 x 10 $^{-4}$ M thymidine. When 10 ml of 50X concentrate are diluted to 500 ml with sterile tissue culture medium, the final concentrations of hypoxanthine and thymidine are 100 μ M and 16 μ M, respectively.

Storage/Stability

Store at -20 °C. After reconstitution, store at 2-8 °C.

References

- Freshney, R. I., Culture of Animal Cells: A Manual of Basic Technique, 3rd ed. (John Wiley & Sons, Inc., 1994) pp. 389-391.
- 2. Harlow, E. and D. Lane, *Antibodies: A Laboratory Manual*, (Cold Spring Harbor Laboratory, 1988) p. 280.
- 3. Kennett, R. H., *Meth. Enzymol.*, eds. W. B. Jakoby and I. H. Pastan (Academic Press, 1979) Vol. LVII, p. 352.

pab 10/09/98

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.