

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Succinic acid Cell Culture Tested

Product Number **S 9512** Store at Room Temperature

Product Description

Molecular Formula: $C_4H_6O_4$ Molecular Weight: 118.1 CAS Number: 110-15-6 Melting point: 185-187 °C¹ Boiling point: 235 °C (with partial decomposition into the anhydride)¹ pK_a: 4.21, 5.72 (25 °C)² Synonyms: butanedioic acid, dicarboxylic acid C₄, ethylenesuccinic acid¹

This product is cell culture tested (0.075 mg/ml) and insect cell culture tested (0.01 mg/ml). It is appropriate for use in cell culture and insect cell culture applications.

Succinic acid is a dicarboxylic acid that occurs in nature in such organisms as fungi and lichens.¹ In eukaryotes, succiniate, the anion of succinic acid, is an intermediate in the citric acid cycle, being formed from succinyl CoA and being converted to fumarate. Bacteria and plants can produce succinate from acetate or acetyl CoA, in the glyoxylate cycle.³

Succinic acid has been used in large scale applications including the manufacture of lacquers, dyes, and esters for perfumes.¹ In polymer research, succinic acid has been utilized to prepare biocompatible hybrid dendritic-linear polyester-ethers.⁴ A study of the co-crystallization of cis-itraconazole with various 1,4-dicarboxylic acids, including succinic acid, has been reported.⁵

Succinic acid has been used as a matrix in infrared (IR) MALDI analytical methods.^{6,7,8} An analytical study of various low molecular weight organic acids, including succinic acid, using capillary zone electrophoresis-electrospray ionization mass spectrometry has been published.⁹

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in water (100 mg/ml), yielding a clear, colorless solution. It is also soluble in alcohol (54 mg/ml), methanol (158 mg/ml), acetone (27 mg/ml), glycerol (50 mg/ml), and ether (8.8 mg/ml). It is essentially insoluble in benzene, carbon disulfide, and petroleum ether.¹

References

- 1. The Merck Index, 12th ed., Entry# 9038.
- Data for Biochemical Research, 3rd ed., Dawson, R. M. C., et al., Oxford University Press (New York, NY: 1986), p. 50.
- 3. Biochemistry, 3rd ed., Stryer, L., W. H. Freeman (New York, NY: 1988), pp. 376-378, 388.
- 4. Carnahan, M. A., et al., Hybrid dendritic-linear polyester-ethers for in situ photopolymerization. J. Am. Chem. Soc., **124(19)**, 5291-5293 (2002).
- Remenar, J. F., et al., Crystal engineering of novel cocrystals of a triazole drug with 1,4-dicarboxylic acids. J. Am. Chem. Soc., **125(28)**, 8456-8457 (2003).
- Talrose, V. L., et al., Insight into absorption of radiation/energy transfer in infrared matrixassisted laser desorption/ionization: the roles of matrices, water and metal substrates. Rapid Commun. Mass Spectrom., **13(21)**, 2191-2198 (1999).
- Budnik, B. A., et al., Benefits of 2.94 micron infrared matrix-assisted laser desorption/ionization for analysis of labile molecules by Fourier transform mass spectrometry. Rapid Commun. Mass Spectrom., **14(7)**, 578-584 (2000).

- Feldhaus, D., et al., Influence of the laser fluence in infrared matrix-assisted laser desorption/ionization with a 2.94 μm Er:YAG laser and a flat-top beam profile. J. Mass Spectrom., 35(11), 1320-1328 (2000).
- Hagberg, J., Analysis of low-molecular-mass organic acids using capillary zone electrophoresiselectrospray ionization mass spectrometry. J. Chromatogr. A, **988(1)**, 127-133 (2003).

GCY/JRC 8/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.