




Merck Biosciences

Calbiochem | Novabiochem | Novagen

Novabiochem®

Letters: 03/04

Contents


NEW Derivatives for the introduction of

D-phosphoamino acids	1
<i>NEW</i> Derivative for the introduction of sulfotyrosine.	2
<i>NEW</i> Derivative for the introduction of arginine	2
<i>NEW</i> Glycosylated amino acid derivatives	3
<i>NEW</i> Fmoc-amino acid pentafluorophenyl esters	3
<i>NEW</i> Derivatives for enhanced peptide synthesis	3
<i>NEW</i> Fluorescent dyes for peptide labeling.	4
28th European Peptide Symposium	4

Product Focus: Amino acids for Fmoc SPPS

NEW Derivatives for the introduction of D-phosphoamino acids

Fmoc-D-Aaa(PO(OBzl)OH)-OH; Aaa=Ser, Thr, Tyr



Features & Benefits

- Direct synthesis of phosphopeptides by automated methods
- Introduced using HBTU/DIPEA activation methods
- Compatible with Fmoc SPPS
- No need for difficult post-synthetic phosphorylation

Since their commercial introduction by Novabiochem in 1994, Fmoc-Ser(PO(OBzl)OH)-OH [1], Fmoc-Thr(PO(OBzl)OH)-OH [2, 3], and Fmoc-Tyr(PO(OBzl)OH)-OH [2, 3] have become



the reagents of choice for the synthesis of phosphopeptides by Fmoc SPPS. The use of these building blocks enables phosphoamino acids to be directly incorporated into a sequence during solid phase synthesis, thereby eliminating the need for complex post-synthetic phosphorylation.

Coupling of these derivatives is best achieved using TBTU or HBTU activation, as described in Method 1, with a 3-fold excess of base relative to amino acid derivative [4]. With the hindered Thr derivative, the coupling can be sluggish, particularly to β -branched amino acids, so an extended reaction time may be required. Furthermore, the coupling of consecutive phospho residues is also difficult. When synthesizing peptides containing several phospho residues, the excess of activated amino acid may need to be increased to compensate for loss of reagent that results from acylation of the piperidine associated with each acidic phosphate group.

The benzyl side-chain protecting group is normally removed in 3 hours during the course of the standard 95% TFA cleavage reaction. For peptides containing several phospho residues, a longer cleavage time may be required.

Method 1: TBTU activation of Fmoc-Aaa(PO(OBzl)OH)-OH

1. Dissolve Fmoc-Aaa(PO(OBzl)OH)-OH (5 eq.^a), TBTU/HBTU (5 eq.^a) and HOBt (5 eq.^a) in the minimum volume of DMF.
2. Add DIPEA (15 eq.^a) to mixture, mix and add immediately to Fmoc-deblocked peptide resin.
3. Allow to couple for 1 - 2h.
4. Check completeness of reaction with the Kaiser or TNBS test. Wash resin and repeat reaction if necessary.

^aRelative to resin substitution.

Method 2: Cleavage and side-chain deprotection of phosphopeptides prepared by Fmoc SPPS

1. Place dry resin in a flask and add TFA solution containing appropriate scavengers^b (10 - 25 ml/g). Stopper flask and leave to stand at rt for 3 h with occasional swirling.
2. Remove the resin by filtration under reduced pressure. Wash the resin twice with TFA. Combine filtrates and add dropwise 8 - 10 volumes of cold ether.
3. Isolate peptide by filtration or centrifugation.

^bSee Novabiochem 2004/5 catalog, page 3.16 for appropriate selection of scavengers.

04-13-1078 Fmoc-D-Ser(PO(OBzl)OH)-OH
NEW

1 g
5 g

04-13-1079 Fmoc-D-Thr(PO(OBzl)OH)-OH
NEW


1 g
5 g

04-13-1080 Fmoc-D-Tyr(PO(OBzl)OH)-OH
NEW

1 g
5 g

NEW Derivative for introduction of sulfotyrosine

Fmoc-Tyr(SO₃·NnBu₄)-OH



Features & Benefits

- Direct synthesis of sulfopeptides by automated methods
- No need for complex off-instrument chemistry
- Compatible with standard Fmoc peptide assembly protocols
- Introduced using HBTU/DIPEA activation methods
- Increased acid-stability compared with sodium or barium salts
- Cleave with 90% TFA at 0 °C to minimize desulfation

Sulfation of tyrosine is an important post-translational modification [5] which is involved in protein-protein recognition [6] and is found in a number of biologically-active peptides such as gastrin II, cholecystokinin, and caerulein [7].


For the Fmoc SPPS of sulfotyrosine-containing peptides, Novabiochem has introduced Fmoc-Tyr(SO₃·NnBu₄)-OH. In contrast to other commercially available building blocks for the introduction of this amino acid, which use a metal ion as the sulfate counterion, Novabiochem's derivative employs tetrabutylammonium as this is thought to reduce desulfation during TFA cleavage [8, 9]. Fmoc-Tyr(SO₃·NnBu₄)-OH is best coupled using TBTU/DIPEA activation. Cleavage from the resin should be carried out using 90% TFA at 0 °C. Using this approach, desulfation can be reduced to less than 15% [7].

04-12-1251 Fmoc-Tyr(SO₃·NnBu₄)-OH
NEW

1 g
5 g

NEW Derivative for introduction of arginine

Fmoc-Arg(Boc)₂-OH



Features & Benefits

- No sulfonated modifications of Trp or Tyr
- Compatible with standard Fmoc protocols



Fmoc-Arg(Boc)₂-OH is an excellent derivative for the Fmoc SPPS of arginine-containing peptides [10]. Unlike other arginine derivatives such as Arg(Pmc) or Arg(Pbf), deprotection of the side-chain of Arg(Boc)₂ does not release by-products that can cause sulfonation of Trp or Tyr residues [11, 12]. Fmoc-Arg(Boc)₂-OH can be coupled using standard activation methods such as PyBOP or TBTU, although longer reaction times compared with Arg(Pbf) may be required due to the bulkiness of the side-chain protecting groups.

04-12-1249 Fmoc-Arg(Boc)₂-OH
NEW

1 g
5 g

NEW Glycosylated amino acid derivatives

Fmoc-Ser(Ac₃AcNH- α -Gal)-OH Fmoc-Thr(Ac₃AcNH- α -Gal)-OH



Features & Benefits

- Direct synthesis of glycopeptides using automated methods
- Compatible with standard Fmoc protocols

In contrast with N-linked glycosides, which all contain a mannotrioso-di-N-acetyl-chitobiose core linked via a $\beta(1\rightarrow N\beta)$ linkage to Asn, O-linked glycosides are structurally diverse. One of the most important classes of O-linked glycosides are based on 2-acetamido-2-deoxy- α -D-galactopyranosyl (AcNH- α -Gal) attached to serine or threonine. Such O-glycosides are found in a wide range of proteins, such as mucin secreted from epithelial cells, the tumor-associated Tn-antigen, and gp120 from HIV [13]. Thus, synthetic O-linked glycopeptides are important tools for elucidating the role of O-linked glycosides and for drug discovery and design.

To facilitate the synthesis of such peptides, Novabiochem has introduced Fmoc-Ser(Ac₃AcNH- α -Gal)-OH and Fmoc-Thr(Ac₃AcNH- α -Gal)-OH. The O-glycosidic linkage and the O-acetyl protection in these building blocks are stable to both piperidine and TFA [14], making them completely compatible with standard protocols in Fmoc solid phase peptide synthesis. Removal of the acetyl protecting groups is best carried out by treatment of the peptidyl resin with methanolic ammonia prior to cleavage.

04-12-8103 Fmoc-Ser(Ac₃AcNH- α -Gal)-OH 100 mg

NEW



04-12-8104 Fmoc-Thr(Ac₃AcNH- α -Gal)-OH 100 mg

NEW


NEW Fmoc-amino acid pentafluorophenyl esters

Fmoc-Orn(Boc)-OPfp

Fmoc-D-Orn(Boc)-OPfp



Fmoc-Cit-OPfp



Novabiochem is pleased to offer the pentafluorophenyl esters of Fmoc-protected L- and D-ornithine, and Fmoc-citrulline.

04-12-1545 Fmoc-Cit-OPfp

NEW

1 g

5 g

04-12-1546 Fmoc-Orn(Boc)-OPfp

NEW

1 g

5 g

04-13-1507 Fmoc-D-Orn(Boc)-OPfp


NEW

1 g

5 g

NEW Derivative for enhancing peptide synthesis

Fmoc-Gly-(Hmb)Gly-OH



Features & Benefits

- "Pseudoproline effect" for Gly-Gly-containing sequences
- Eliminates Gly deletions in Gly-Gly repeats
- Compatible with standard Fmoc protocols
- Hmb group removed during TFA cleavage

Fmoc-Gly-(Hmb)Gly-OH is a novel building block for improving the synthesis of peptides containing the commonly occurring Gly-Gly motif. The production of such peptides is often problematic, owing to their propensity to aggregate and the difficulties in separating closely eluting desGly by-products. Insertion of Gly-(Hmb)Gly into a peptide sequence has the same benefits on synthetic efficiency as a pseudoproline dipeptide, preventing aggregation and improving acylation and deprotection kinetics. The use of the analogous Fmoc-Gly-(Dmb)Gly-OH was found to be essential for the synthesis of peptides related to nucleolin [15].

04-12-1253 Fmoc-Gly-(Hmb)Gly-OH


NEW

1 g

5 g

New Fluorescent dye for peptide-labeling

5(6)-Carboxytetramethylrhodamine



This high purity fluorescent dye is a useful tool for preparing fluorescently-labeled peptides and fluorescence-quenched peptide substrates. Novabiochem's 5(6)-carboxytetra-methylrhodamine (TAMRA; λ_{ex} 555 nm; λ_{em} 580 nm) is supplied as a mixture of 5- and 6-isomers, making it cost effective for those applications which do not require a single isomer dye.

5(6)-TAMRA is conveniently introduced during solid phase synthesis by coupling to N-terminal or side-chain amino groups. When it is to be located on a side-chain amino group, the simplest approach is to incorporate an orthogonally-protected derivative, such as Lys(Mtt) or Lys(ivDde), which can be later selectively deprotected on the resin immediately prior to coupling of the dye.

01-63-0134	5(6)-Carboxytetramethylrhodamine	100 mg
NEW		500 mg
Novabiochem's single isomer dyes		
01-63-0114	5-Carboxytetramethylrhodamine	10 mg
		50 mg
01-63-0115	6-Carboxytetramethylrhodamine	10 mg
		50 mg

Merck Biosciences AG · Switzerland
Weidenmattweg 4
4448 Läufelfingen
Phone +41 (62) 285 2525
Fax +41 (62) 285 2520

www.novabiochem.com

Calbiochem, Novabiochem, and Novagen are brands of Merck Biosciences, an affiliate of Merck KGaA, Darmstadt, Germany.
NVL304-NP

28th European Peptide Symposium

Visit our stand at the 28th European Peptide Symposium in Prague to find out more about the products in this letter and to pick-up your copies of the 2004 Innovations and 2004/2005 Catalog. Whilst there take the opportunity to talk to our technical support staff, who will be on-hand each day of the conference to answer your questions.

References

1. T. Wakamiya, et al. (1994) *Chem. Lett.*, 1099.
2. P. White & J. Beythien in "Innovations & Perspectives in Solid Phase Synthesis & Combinatorial Libraries, 4th International Symposium", R. Epton (Ed.), Mayflower Scientific Ltd., Birmingham, 1996, pp. 557.
3. T. Vorherr & W. Bannwarth (1995) *Bioorg. Med. Chem. Lett.*, **5**, 2661.
4. P. White in "Peptides, The Wave of the Future, Proc. 17th American Peptide Symposium", M. Lebl & R. A. Houghten (Eds), Kluwer, Norwell, 2001, pp. 97.
5. C. Niehrs, et al. (1994) *Chem.-Biol. Interact.*, **92**, 257.
6. T. Young & L. L. Kiessling (2002) *Angew. Chem. Int. Ed.*, **41**, 3449: ref [2].
7. K. Kitagawa, et al. (2001) *J. Org. Chem.*, **66**, 1.
8. M. Ueki, et al. in "Peptide Science, Proc. Japanese Peptide Symposium", Kyoto, 1999, pp. 117.
9. S. V. Campos, et al. (2002) *J. Chem. Soc., Perkin Trans. 1*, 682.
10. A. S. Verdini, et al. (1992) *Tetrahedron Lett.*, **33**, 6541.
11. B. Riniker & A. Hartmann in "Peptides, Chemistry, Structure & Biology, Proc. 11th American Peptide Symposium", J.E. Rivier & G.R. Marshall (Eds.), Escom, Leiden, 1990, pp. 950.
12. C.G.Fields & G.B.Fields (1993) *Tetrahedron. Lett.*, **34**, 6661.
13. N. Sharon & H. Lis in "Glycosciences: status & perspectives", H.-J. Gabius & S. Gabius (Eds.), Chapman & Hall, Weinheim, 1997, pp. 133.
14. J. Kihlberg in "Fmoc solid phase peptide synthesis: a practical approach", W. C. Chan & P. D. White (Eds.), Oxford University Press, Oxford, 2000, pp. 195.
15. S. Zahariev, et al. in "Peptides 2002, Proc. 27th European Peptide Symposium", E. Benedetti & C. Pedone (Eds.), Edizioni Ziino, Naples, 2002, pp.336.