



RESTRICTION ENDONUCLEASE Ban II

Product No. **R 2630**

Store at 0 to -20°C

Product Summary

Recognition Sequence: 5' GpuGCPy/C '3

Activity: 10,000 units/ml

Cutting: 100%

Ligation: >95%

Recutting: >95%

No degradation detected with >10 units for 16 hrs.

Fold over digestion: 160 (10 units x 16 hrs.)

Package Size: 1000 units

Unit Definition

One unit is the enzyme activity that completely cleaves 1 µg λ DNA in 1 hr. at 37°C in a total volume of 50 µl of Buffer SB for restriction enzymes.

Specificity

Ban II recognizes the sequence GpuGCPy/C and generates fragments with 3'-cohesive termini.¹

Comments

Digestion Buffer SB is supplied as a 10x concentrate.

Heat inactivation information for Ban II is not available.

Ban II Storage and Dilution Buffer Composition

10 mM Tris-HCl

100 mM KCl

0.1 mM EDTA

10 mM 2-mercaptoethanol

500 µg/ml bovine serum albumin

50% (v/v) glycerol

pH 7.5

Product Information

1x Digestion Buffer SB (B8781) Composition for Ban II: 100 % Digestion at 37°C.

10 mM Tris-HCl

100 mM NaCl

5 mM MgCl₂

1 mM 2-mercaptoethanol

pH 8.0

Quality Control Testing

Absence of unspecific endonuclease activities:

1 µg λDNA is incubated for 16 hrs. in 50 µl buffer SB with excess of Ban II.

Ligation and recutting assay

Ban II fragments, obtained by complete digestion of 1 µg λDNA, are adjusted to pH 7.5 at 20°C. The Ban II fragments are then ligated with 1.0 unit T4-DNA ligase at pH 7.5 at 4°C. A 10 µl reaction mixture, incubated for 16 hrs. at pH 7.5 at 4°C contained 1.0 unit T4- DNA ligase, 66 mM Tris-HCl, 5 mM MgCl₂, 1 mM ATP and 1 mM dithioerythritol.

The degree of ligation and subsequent recutting with Ban II to yield the typical pattern of λDNA-Ban II fragments is determined.

References

1. Sugaisaki, H., et al., Nucl. Acids Res., **10**, 5747 1982.

6/99