



Speed meets performance!

Your first choice for ultra-fast HPLC

Purospher STAR RP-18 endcapped UHPLC columns



When speed meets performance

Purospher® STAR RP-18 endcapped UHPLC columns for Ultra fast HPLC use.

Purospher® STAR RP-18 endcapped, 2 µm and 3 µm UHPLC columns are ideal for ultra fast chromatography applications. The most important attribute of these columns is speed with highest resolution, performance and stability at the same time. They are also ideal for applications requiring maximum efficiency and resolution.

Purospher® STAR RP-18 endcapped UHPLC columns speed separation up to 10 times for:

- high throughput analysis
- process monitoring
- QC analyses
- method development
- LC/MS


and save up to 84.5 % solvent



Hibar® HR Purospher® STAR
UHPLC Columns

Surface modification of Purospher® STAR RP-18 endcapped

Purospher® STAR RP-18 endcapped HPLC columns are designed for universal use. It doesn't matter if the samples are basic, neutral, metal chelating or indeed any other format – Purospher® STAR can do it, naturally without peak tailing!






Your benefits


- High purity of sorbent – best peak symmetry for giving correct results
- High pressure stability (600 bar) – ideal for UHPLC use
- Very broad selectivity – one column type for most applications
- High separation efficiency – high sensitive results
- Excellent pH stability (pH 1.5 - 10) – extremely wide application range

High resolution separation of 9 Alkylphenoles

(Speed up separation from 20 min to 2 min)




Column: LiChroCART® 150-4.6
Purospher® STAR RP-18 endcapped, 5 µm
Temperature: 40 °C
Mobile phase: A: Water
B: Acetonitrile
Gradient: 0 min 45 % B,
from 45 to 95 % B in 15 min,
from 15.1 to 20 min
reequilibration with 45 % B
Flow rate: 1.0 ml/min,
Pressure: 105 bar
Detection UV 247 nm
Inj. volume: 10 µl



Column: Hibar® HR 50-2.1
Purospher® STAR RP-18 endcapped, 2 µm
Temperature: 40 °C
Mobile phase: A: Water
B: Acetonitrile
Gradient: 0 min 55 % B,
from 55 to 100 % B in 0.8 min,
from 0.9 to 2 min
reequilibration with 55 % B
Flow rate: 1.1 ml/min,
Pressure: 505 bar
Detection UV 247 nm
Inj. volume: 1 µl


Separation of 7 Phenoles in 2.5 min




Column: Hibar® HR 50-2.1
 Purospher® STAR RP-18 endcapped, 2 µm
Temperature: 25 °C
Eluent A: 0.5 % Phosphoric acid
Eluent B: Acetonitrile
Gradient: 0 min 34 % B, from 34 to 55 % B in 0.8 min,
 from 55 % to 100 % B in 1 min,
 hold 100 % B to 1.5 min, 34 % B in 1.5 min and
 reequilibration up to 2.5 min
Flow rate: 0.9 ml/min,
Pressure: 515 bar
Detection: UV 225 nm
Inj. volume: 1 µl

Sample:


1. Phenol




2. 2-Chlorophenol




3. 2,4-Dinitrophenol




4. 2-Nitrophenol




5. 4-Chlor-3-Methylphenol



6. 2,4,6-Trichlorophenol



7. Pentachlorophenol



Ordering information

Product	Content	Ord. No.
Hibar® HR 30-2.1 Purospher® STAR RP-18 endcapped, 2 µm	1 piece	1.50645.0001
Hibar® HR 50-2.1 Purospher® STAR RP-18 endcapped, 2 µm	1 piece	1.50646.0001
Hibar® HR 30-2.1 Purospher® STAR RP-18 endcapped, 3 µm	1 piece	1.50650.0001
Hibar® HR 50-2.1 Purospher® STAR RP-18 endcapped, 3 µm	1 piece	1.50651.0001

For more information about our columns and solvents please contact:

www.merck-chemicals.com/chromatography



W288120 10/09, 1st issue

Merck KGaA
64271 Darmstadt, Germany
Fax: +49 (0) 61 51 / 72-60 80
E-mail: chromatography@merck.de
www.merck-chemicals.com/chromatography

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose. Purospher® is a registered trademark of Merck KGaA, Darmstadt, Germany.