

α -Cyano-4-hydroxycinnamic Acid


Product No. **C2020**

Store at 2-8°C

CAS RN 28166-41-8

Synonyms: α -CCA; HCCA; 4-HCCA; α -CHCA;
 α -Cyano; ACCA

Product Description



Appearance: Yellow powder
Molecular Formula: $C_{10}H_7NO_3$
Molecular Weight: 189.17

A specific covalent inhibitor of mitochondrial lactate and pyruvate transport.¹⁻⁴ When rat muscle was examined, tracer lactate uptake was inhibited by 5 mM α -CCA.⁴ It has also been reported to inhibit beta-cell apical anion exchange ($IC_{50} = 2.4$ mM).⁵

Cinnamic acid derivatives are commonly used in Matrix-Assisted Laser Desorption Ionization (MALDI) mass spectrometry. Of the cinnamic acid derivatives, α -cyano-4-hydroxycinnamic acid (α -CCA), sinapinic acid (3,5-dimethoxy-4-hydroxycinnamic acid, Product No. D7927), and ferulic acid (3-methoxy-4-hydroxycinnamic acid, Product No. F3500) are the most popular. The most important properties of MALDI matrices are absorption wavelength, solubility, and crystal formation. α -CCA absorbs well at 337 nm, the output of a nitrogen laser. It has a relatively high solubility in organic solvents and a moderate solubility in water (see Preparation Instructions), and forms fine crystals upon drying from solution.⁶ In addition, α -CCA has a lower tendency to form adducts with the analyte, simplifying the mass spectra.⁶ Finally, peaks arising from the matrix are not observed above a mass to charge ratio (m/z) of 400, making α -CCA particularly amenable to the analysis of low molecular weight peptides.⁶

α -CCA was introduced as a matrix for MALDI time-of-flight (TOF) mass spectrometry in 1992 and is most often used in the analysis of peptides and proteins of less than 10 kDa.⁶ Recently, α -CCA was applied to the analysis of DNA by MALDI-TOF mass spectrometry⁷ and the analysis of peptides and proteins in Fourier-transform mass spectrometry.⁸ Antibiotics,⁹ peptide nucleic acids (a new class of DNA mimics),¹⁰ peptides¹¹ and proteins with masses as high as 66,000 Da have also been successfully analyzed by using this as a matrix.⁸

Preparation Instructions

α -CCA is soluble in methanol at 50 mg/ml, in acetonitrile at ~35 mg/ml, and in water at ~6 mg/ml.

Storage/Stability

Sigma has not tested the solution stability of this chemical.

References

1. Williams, R.J., *et al.*, *Neuroscience*, **74**, 461-468 (1996).
2. Schneider, U., *et al.*, *Neuroscience*, **53**, 1153-1162 (1993).
3. Barnard, J.P., *et al.*, *J. Biol. Chem.*, **268**, 3654-3661 (1993).
4. Watt, P.W., *et al.*, *Biochim. Biophys. Acta*, **944**, 213-222 (1988).
5. Emmons, C., *Am. J. Physiol.*, **276(4 Pt 2)**, F635-F643 (1999).
6. Beavis, R.C. *et al.*, *Organic Mass Spectrometry*, **27**, 156-158 (1992).
7. Koomen, J.M. *et al.*, *Anal. Chem.*, **72**, 3860-3866 (2000).
8. Sze, E.T., *et al.*, *J. Am. Soc. Mass. Spectrom.*, **9**, 166-174 (1998).
9. Ling, Y.C., *et al.*, *Rapid Commun. Mass Spectrom.*, **12**, 317-327 (1998).
10. Butler, J.M., *et al.*, *Anal. Chem.*, **68**, 3283-3287 (1996).
11. Patterson, S.D., *et al.*, *Anal. Chem.*, **66**, 3727-3732 (1994).

MAM,ALC 02/06-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.