Sigma-Aldrich.

Product Information

SIGMA*FAST*[™] Fast Red TR/Naphthol AS-MX Tablets

Tablet, to prepare 1 mL

F4648

Product Description

Fast Red TR/Naphthol AS-MX is the immunohistology substrate of choice for antibodies conjugated to alkaline phosphatase, as it produces an intense red stain. Slides stained with Fast Red TR/Naphthol AS-MX must be cover-slipped using aqueous mounting media, as the reaction product is alcohol-soluble.

SIGMAFAST[™] Fast Red TR/Naphthol AS-MX Phosphate (4-Chloro-2-methylbenzenediazonium/ 3-Hydroxy-2-naphthoic acid 2,4-dimethylanilide phosphate) tablets have been developed for use in immunohistology and blotting, as a precipitating substrate for the detection of alkaline phosphatase activity. Levamisole has been added to a concentration of 0.15 mg/mL to block endogenous alkaline phosphatase activity.

SIGMAFAST[™] Fast Red TR/Naphthol AS-MX tablets require no additional buffers or steps to prepare an active substrate solution. One Fast Red TR/Naphthol AS-MX tablet and one Trizma[®] buffer tablet, dissolved in 1 mL of deionized or distilled water, provides 1 mL of ready-to-use substrate. Each SIGMAFAST[™] Fast Red TR/Naphthol AS-MX tablet set contains the following when dissolved in 1 mL H₂O:

- Fast Red TR: 1.0 mg/mL
- Naphthol AS-MX: 0.4 mg/mL
- Levamisole: 0.15 mg/mL
- Trizma[®] Buffer: 0.1 M

This product has been used to study various systems such as animal models, 1,2 zebrafish tissue samples, $^{3-5}$ and endothelial cells. $^{6-8}$ Several theses 9,10 and dissertations $^{11-26}$ have cited use of product F4648 in their research protocols.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Store the tablets at -20 °C.

Components

SIGMAFAST[™] Fast Red TR/Naphthol AS-MX Phosphate Tablets (Component Number F0775): 5 tablets (for 5SET) or 50 tablets (for 50SET)

Trizma[®] Buffer Tablets (Component Number T9043): 5 tablets (for 5SET) or 50 tablets (for 5OSET)

Reagents and Equipment Required but Not Provided

- Distilled or deionized water
- Pipettes
- Test tubes
- 0.2 µm filter (such as Cat. No. WHA10462701)

Preparation Instructions

- 1. Remove the required number of Fast Red TR/ Naphthol AS-MX and Trizma[®] tablet packages from the freezer.
- 2. Allow the tablets to reach room temperature.
- Open the Trizma[®] tablet package (gold foil) and drop the tablet into an appropriate container.
 Do not touch the tablet with your fingers.
- 4. Add 1 mL of distilled or deionized water.
- 5. Vortex until dissolved.
- Open the Fast Red TR/Naphthol AS-MX tablet package (silver foil). Do not touch the tablet with your fingers.
- 7. Drop one Fast Red TR/Naphthol AS-MX tablet into the Trizma[®] buffer. Vortex until dissolved.

The SIGMAFAST[™] Fast Red TR/Naphthol AS-MX substrate is now ready for use. For best results, the solution should be used within one hour.



Procedure

- 1. Cover the tissue section with 0.1-0.2 mL of Fast Red TR/Naphthol AS-MX solution.
- Fast Red TR/Naphthol AS-MX is a fast-reacting substrate. It should be carefully monitored during the reaction to prevent overdevelopment and high background. Reactions may be stopped by gently washing the slide in water.
- Occasionally, the Fast Red TR/Naphthol AS-MX solution may be hazy. The haziness may be removed by filtering the Fast Red TR/Naphthol AS-MX solution through a 0.2 µm filter.
- When finished, dispose of any remaining substrate solution in a manner consistent with proper hazardous material handling protocols for your institution.

Troubleshooting

Background is too high

- Use a blocking step prior to the application of the primary antibody. Diluted normal serum (10% v/v) from the same species as the secondary antibody generally produces the best results.
- 2. Decrease the staining time.
- 3. Titer the conjugate to optimize the working dilution.

No color develops or color is too faint

- 1. Adjust the concentration of the primary antibody.
- Adjust the concentration of the secondary antibody.
- 3. Determine if the enzyme conjugate is active.
- 4. Consider using an amplifying system such as avidin-biotin.
- 5. Increase the staining time.
- Determine if enzymatic treatment (unmasking) of the antigen is required prior to application of the primary antibody.

References

- Bany, B.M., "Pseudopregnant Bead-Induced Mouse Deciduoma Model", in *The Guide to Investigation of Mouse Pregnancy* (B.A. Croy *et al.*, eds.). Academic Press/Elsevier (London, UK), Chapter 42, pp. 499-504 (2014).
- Benten, D. et al., Methods Mol. Biol., 326, 189-201 (2006).
- Jowett, T., *Methods Mol. Biol.*, **97**, 461-486 (1999).

- 4. Webb, A.E., and Kimelman, D., *Methods Mol. Biol.*, **289**, 137-146 (2005).
- van Boxtel, A.L. *et al.*, *Dev. Cell*, **44(2)**, 179-191 (2018).
- Clancy, R. et al., Arthritis Rheum., 44(5), 1203-1208 (2001).
- Lacorre, D.-A. *et al.*, *Blood*, **103(11)**, 4164-4172 (2004).
- Dath, C. et al., Hum. Reprod., 26(6), 1431-1439 (2011).
- Swanson, Annika J., "Identification of Potential New Biomarkers of Oxygenated Polycyclic Aromatic Hydrocarbon (OPAH) Exposures in Zebrafish". Oregon State University, B.S. Honors Scholar thesis, p. 9 (2014).
- Tatari, Nazarin, "The Effect of Semaphorin 3E on Angiogenesis in Murine Model of Allergic Asthma". University of Manitoba, M.Sc. thesis, p. 58 (2015).
- Ebensen, Thomas, "Das SCIDHu-Tiermodell als Beispiel für eine diskordante xenogene Transplantationssituation: Untersuchungen zur Rolle humaner T-Lymphozyten in der Abstoßungsreaktion" ("The SCIDHu animal model as an example of a discordant xenogeneic transplantation situation: studies on the role of human T lymphocytes in the rejection reaction"). Universität Hannover, Dr. rer. nat. dissertation, p. 36 (2000).
- Lücke, Sonja, "Zyklusabhängige Veränderungen am Eileiter des Rindes und Charakterisierung boviner Eileiterepithelzellen in Suspensionskultur" ("Cycle-dependent changes in the bovine oviduct and characterization of bovine epithelial cells in suspension culture"). Ludwig-Maximilians-Universität München, Dr. Vet. Med. dissertation, p. 43 (2005).
- Cuentas, Edwin Roger Parra, "Pneumonias intersticiais idiopáticas: da patogênese e do remodelamento aos determinantes anátomoclínico-radiológicos de prognóstico e sobrevida com ênfase ao componente vascular" ("Idiopathic interstitial pneumonia: from pathogenesis and remodeling to anatomo-clinical-radiological determinants of prognosis and survival with emphasis on the vascular component"). Universidade de São Paulo, Ph.D. dissertation, p. 90 (2006).
- 14. Economou, Andrew David, "Phylogenetic and Developmental Studies into the Evolution of an Insect Novelty". University College London, Ph.D. dissertation, p. 70 (2008).

- Wiedemann, Christine Hyun-Zu, "Die Wirkung von Chlormadinonacetat auf humane Melanozyten" ("The effect of chlormadinone acetate on human melanocytes"). Ludwig-Maximilians-Universität zu München, Dr. med. dissertation, p. 26 (2010).
- Glasco, Derrick M., "The Role of Wnt/Planar Cell Polarity Signaling in Mouse Facial Branchiomotor Neuron Migration". University of Missouri Columbia, Ph.D. dissertation, p. 92 (2011).
- Goller, Katja Verena, "Pathogens in free-ranging African carnivores: evolution, diversity and co-infection". Humboldt-Universität zu Berlin, Dr. rer. nat. dissertation, p. 99 (2011).
- Neascu, Cristian Dan, "Untersuchungen zur Funktion der Knorpelproteine Ucma und Matrilin-1 im Zebrafisch" ("Investigations on the function of the cartilage proteins Ucma and Matrilin-1 in zebrafish"). Universität zu Köln, Dr. rer. nat. dissertation, p. 59 (2013).
- Powell, Davalyn R., "The Role of Prdm1a in Zebrafish Neural Crest Development". University of Colorado, Ph.D. dissertation, p. 36 (2014).
- Leonel, Ellen Cristina Rivas, "Efeitos da exposição neonatal ao bisfenol A e ao 17-β estradiol sobre a glândula mamária de fêmeas adultas de gerbilo da Mongólia" ("Effects of neonatal exposure to bisphenol A and 17-β estradiol on the mammary gland of adult female Mongolian gerbils"). Universidade Estadual Paulista "Júlio de Mesquita Filho", Ph.D. dissertation, p. 95 (2015).
- Incarbone, Marco, "In vivo study of the suppression of cell-autonomous and systemic RNA silencing by the Peanut clump virus protein P15". Université de Strasbourg, Ph.D. dissertation, p. 108 (2016).
- 22. AmbuAli, Aisha, "Morphological and Functional Aspects of Feeding in The Freshwater Fish Louse *Argulus foliaceus* (Linnaeus, 1758)". University of Stirling, Ph.D. dissertation, p. 74 (2017).
- O'Brien, Katie Alice, "Investigation into metabolic profile changes in environmental hypoxia and the potential for dietary nitrate to alleviate hypoxic stress". King's College London, Ph.D. dissertation, p. 203 (2017).

- 24. Alberstat, Erin Jarvis, "On Monsters, Patterns and Appendages: Targeted CRISPR/Cas9 genome editing in the amphipod crustacean *Parhyale hawaiensis* to explore Hox gene interactions in the establishment of appendage diversity". University of California Berkeley, Ph.D. dissertation, p. 120 (2018).
- van den Burg, Nicole Monica Dawn, "Reducing allergic airway inflammation with high-density microprojection array skin patches". University of Queensland, Ph.D. dissertation, p. 128 (2018).
- Bovolini, José Antonio Franchi, "The pathophysiological role of physical inactivity and fat-diet to development of metabolic syndrome in an animal model". University of Porto, Ph.D. dissertation, p. 88 (2019).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at <u>SigmaAldrich.com/techservice</u>.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at <u>SiamaAldrich.com/terms</u>.

Contact Information

For the location of the office nearest you, go to <u>SigmaAldrich.com/offices</u>.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. © 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. F4648dat Rev 03/22 PCS.RBG.GCY.MAM

