

## Product Information

### Z-Leu-Glu(OMe)-Glu(OMe)-Asp(OMe) Fluoromethyl Ketone

Product Number **C 8859**

Storage Temperature  $-20^{\circ}\text{C}$

#### Product Description

Molecular Formula:  $\text{C}_{32}\text{H}_{45}\text{FN}_4\text{O}_{12}$

Molecular weight: 696.7

Z-Leu-Glu(OMe)-Glu(OMe)-Asp(OMe) Fluoromethyl Ketone is the methylated, cell permeable derivative of the caspase inhibitor Z-Leu-Glu-Glu-Asp Fluoromethyl Ketone (Z-LEED-FMK, Z = benzyloxycarbonyl).

Z-LEED-FMK is an inhibitor of caspase 13. Caspase-13 has also been called ERICE (Evolutionarily Related Interleukin-1 $\beta$  Converting Enzyme) and is related to caspase-1, -4, and -5.<sup>1</sup> There is a question whether this enzyme occurs in human tissue.<sup>2</sup>

Apoptosis, or programmed cell death, plays an essential role in development, homeostasis, and defense of multicellular organisms. Among the many known effectors of apoptosis the interleukin-converting enzyme (ICE)-related, cysteine aspartic-specific proteases, or caspases, play a crucial role in apoptosis in almost every cell type.<sup>3,4</sup> At least 14 different caspases have been identified which differ in their substrate specificities.

Methylation of the acidic amino acids Glu and Asp enhances the cell membrane permeability of Z-LEED-FMK. Once in the cell, endogenous esterase activity hydrolyzes the methyl groups to form the biological active form. For *in vitro* studies an esterase needs to be included in the reaction mix to generate the active form of the molecule.

FMK is a trapping group responsible for irreversible inhibition and is also non-cytotoxic. Inhibition occurs when the the FMK group covalently bonds to the  $-\text{SH}$  of an adjacent cysteine residue on the target protein.

Z-LEED-FMK is supplied as a white lyophilized powder.

#### Preparation Instructions

Prepare 20 mM stock solutions in dry ( $\geq 99.9\%$ ) DMSO to maintain product stability. Also soluble in DMF.

#### Storage/Stability

Store at  $-20^{\circ}\text{C}$ . The product is reported to be stable at room temperature for one year in a desiccator. Allow container to warm to room temperature before opening to ensure stability.

Store stock solutions at  $-20^{\circ}\text{C}$  for 6-8 months.

#### References

1. Humke, E. W., et al., ERICE, a novel FLICE-activatable caspase. *J. Biol. Chem.*, **273**, 15702-15707 (1998).
2. Koenig, U., et al., Evidence that caspase-13 is not a human but a bovine gene. *Biochem. Biophys. Res. Commun.*, **285**, 1150-1154 (2001).
3. Nicholson, D.W., and Thornberry, N.A., Caspases: killer proteases. *Trends. Biochem. Sci.*, **22**, 299 (1997).
4. Cohen, G.M., Caspases: the executioners of apoptosis. *Biochem. J.*, **326**, 1 (1997).

BD 8/21/01

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.