

**Product No. I-3144**  
**Insulin-like Growth Factor-II**  
**Human, Recombinant**  
 Expressed in *E. coli*

### Product Description

Insulin-like Growth Factor-II (IGF-II) was first isolated from human serum by Froesch, et al.,<sup>1</sup> as a factor displaying insulin-like activities that were not suppressed by antibodies to insulin. It had been discovered that growth hormone-dependent factors in serum stimulate the incorporation of <sup>35</sup>S into cartilage<sup>2</sup> and that calf serum factors induced cellular division in chick fibroblasts.<sup>3</sup> In 1972, the term "somatomedin" was introduced in an unsuccessful attempt to unify the nomenclature of these hormone-dependent factors.<sup>4</sup> In 1987, a consensus among an international group of scientists endorsed the use of the terms insulin-like growth factors (IGF-I and IGF-II),<sup>5</sup> originally proposed by Rinderknecht and Humbel.<sup>6</sup> Hence, IGF-I and IGF-II have had several synonyms: nonsuppressible insulin-like activity (NSILA), sulfation factor activity (SFA), and multiplication stimulating activity (MSA). Because IGF-II was not regulated by growth hormone, only IGF-I was known as a somatomedin.

Human IGF-II contains 67 amino acids and shares similar structural features with IGF-I, including a 62% sequence homology.<sup>7</sup> In human plasma, IGF-I and IGF-II are associated with IGF-binding proteins<sup>8,9</sup> that transport the polypeptides and partially regulate their actions *in vivo*.<sup>10</sup> In addition to the insulin receptor, IGF-II binds to two forms of IGF receptors, both of which are widely distributed in different tissues and cultured cells.<sup>11</sup> IGF-II is mitogenic for a variety of cultured cells, including mouse 3T3 cells,<sup>12</sup> normal rat kidney cells,<sup>13</sup> human or chicken fibroblasts,<sup>14,15</sup> and MCF-7 human breast carcinoma cells.<sup>16</sup>

### Performance Characteristics

The mitogenic activity of IGF-II was measured by serum-free cell proliferation assay using the bovine kidney cell line MDBK. The EC<sub>50</sub> is defined as the effective concentration of growth factor that elicits a 50% increase in cell growth in a cell based bioassay.

### Product Information

Expressed in *E. coli*  
 Molecular Weight: 7.5 kD  
 Purity: ≥97% by SDS-PAGE and N-terminal analysis  
 EC<sub>50</sub>: 5 - 20 ng/ml  
 Package Size: 50 µg

Formulation: Lyophilized from a 0.2 µm-filtered solution of 10 mM acetic acid, pH 3.4.

Carrier Protein: 2.5 mg bovine serum albumin (BSA)

Sterility: 0.2 µm-filtered, aseptic fill

Endotoxin: ≤0.1 ng/µg IGF-II

### Reconstitution and Use

Reconstitute the contents of the vial using 0.2 µm-filtered 10 mM acetic acid containing 0.1% HSA or BSA to not less than 10 µg/ml. This may be diluted immediately before use in sterile-filtered PBS containing 0.1% - 1.0% BSA.

### Storage

Prior to reconstitution store vial at -20°C. After reconstitution, freeze in working aliquots at -20°C for no longer than 6 months. Prolonged storage and repeated freezing and thawing is **not** recommended.

### References

1. Froesch, E., et al., *J. Clin. Invest.*, **42**, 1816 (1963).
2. Salmon, W., et al., *J. Lab. Clin. Med.*, **49**, 825 (1957).
3. Pierson, R., et al., *J. Cell. Physiol.*, **79**, 319 (1972).
4. Daughaday, W., et al., *Nature*, **235**, 107 (1972).
5. Daughaday, W., et al., *J. Clin. Endocrinol. Metab.*, **65**, 1075 (1987).
6. Kinderknecht, E., and Humbel, R., *Proc. Natl. Acad. Sci.*, **73**, 4379 (1976).
7. Froesch, E., et al., *Ann. Rev. Physiol.*, **47**, 443 (1985).
8. Hardouin, S., et al., *J. Clin. Endocrinol. Metab.*, **69**, 1291 (1989).
9. Zapf, J., et al., *Arch. Biochem. Biophys.*, **168**, 638 (1975).
10. Baxter, R., and Martin, J., *Prog. Growth Factor Res.*, **1**, 49 (1989).
11. Baxter, R., et al., *Adv. Clin. Chem.*, **25**, 49 (1986).
12. Nissley, S., et al., *J. Cell Physiol.*, **89**, 393 (1976).
13. Marquardt, H., et al., *J. Biol. Chem.*, **255**, 9177 (1980).
14. Rechler, M., et al., *J. Clin. Endocrinol. Metab.*, **44**, 820 (1977).
15. Rechler, M., et al., *Nature*, **259**, 134 (1976).
16. Karel, K., et al., *Cancer Research*, **48**, 4083 (1988).