

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone (800) 325-5832 (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

# **ProductInformation**

b-(1® 3,4,6)-Galactosidase, Positionally specific, a mixture from *Streptococcus pneumoniae* and *Xanthomonas sp.*, recombinant, expressed in *E. coli* 

Product Number **G 1288** Storage Temperature 2–8 °C


CAS<sup>#</sup> 9031-11-2 EC 3.2.1.23

Synonyms: Lactase;  $\beta$ -D-Galactoside galactohydrolase;  $\beta$ -D-Galactopyranosidase;  $\beta$ -Lactosidase

## **Product Description**

Two major classes of oligosaccharides (glycans) may be attached to glycoproteins. N-Linked glycans are attached to the amide side-chain of some asparagine (Asn) residues, which form part of the consensus sequence AsnXaaSer/Thr, while O-linked glycans may be added to the hydroxyl side chain of serine or threonine residues. The terminal residues on the glycan chains are commonly sialic acids, which can be removed by the use of a broad-spectrum neuraminidase enzyme. After removal of sialic acids, the galactose residues are exposed. These may be linked to the core glycan in several different positions, the most common of which is via a  $\beta$ -1 $\rightarrow$ 4 bond.

Recombinant  $\beta$ -(1 $\rightarrow$ 3,4,6)-galactosidase, expressed in *Escherichia coli*, is a highly purified enzyme, which releases  $\beta$ -1 $\rightarrow$ 3,  $\beta$ -1 $\rightarrow$ 4, and  $\beta$ -1 $\rightarrow$ 6 linked galactose from the non-reducing end of glycans and glycoproteins. However, the presence of fucose linked to the penultimate N-acetylglucosamine will block cleavage of the galactose residue attached to the same N-acetylglucosamine.



Due to its broad specificity, this enzyme is an extremely useful reagent for detailed structural analysis of glycans in conjunction with other positionally specific galactosidase enzymes.

# Components

 $\beta\text{-}(1{\to}3,4,6)\text{-galactosidase}$  (Product No. G 1288) - The enzyme is supplied in 20 mM Tris HCI, pH 7.5, containing 25 mM NaCI.

Unit Definition: One unit will hydrolyze 1  $\mu$ mole of p-nitrophenyl  $\beta$ -D-galactopyranoside per minute at pH 5.0 at 37 °C.

Protease activity was not detected.

5X Reaction Buffer (Product No. E 5879) – 250 mM sodium phosphate, pH 5.0

# **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### Storage/Stability

It is recommended to store the product at 2–8  $^{\circ}\text{C}.$  Do Not Freeze

### **Procedure**

Add 2  $\mu$ l of enzyme to 100  $\mu$ g of asialoglycoprotein or 1 nmole of oligosaccharide. Add 50 mM sodium phosphate buffer, pH 5.0, and incubate for 1 hour at 37 °C.

#### References

1. Paulson, J., *et al.*, J. Biol. Chem., **253**, 5617-56242 (1978).

AE,MAM 03/05-1