

## Product Information

# Aprotinin Bovine, Recombinant

Expressed in *Nicotiana* (tobacco),  $\geq 5$  TIU/mg protein,  $\geq 98\%$  (SDS-PAGE)**A6103**

## Product Description

CAS Registry Number: 9087-70-1

Synonyms: Antikrein, Antilysin(e), Basic pancreatic trypsin inhibitor (BPTI), Kallikrein-trypsin inactivator, Kunitz protease inhibitor


Molecular formula:  $C_{284}H_{432}N_{84}O_{79}S_7$ 

Molecular weight: 6,512 Da

pI: 10.5

Extinction coefficient:  $E^{1\%} = 8.4$  (280 nm, water)

Aprotinin is a protein consisting of 58 amino acids, arranged in a single polypeptide chain that is cross-linked by three disulfide bridges (depicted below). Aprotinin is a competitive serine protease inhibitor that forms stable complexes with and blocks the active sites of various enzymes, such as those listed in Table 1. The binding is reversible. Most aprotinin-protease complexes dissociate at pH >10 or pH < 3.



This product is a recombinant form of the native bovine-sequence aprotinin. Native bovine aprotinin is traditionally isolated primarily from bovine lung by methods involving fractional precipitation, gel filtration, and ion exchange chromatography. Unlike animal-derived aprotinin, this recombinant aprotinin product is isolated and purified from plant tissue (by proprietary methods).

Several theses<sup>1-5</sup> and dissertations<sup>6-12</sup> have cited use of product A6103 in their protocols.

## Reagent

Purity:  $\geq 98\%$  (SDS-PAGE)Inhibitory Activity:  $\geq 5$  TIU/mg of protein

Unit Definition: One Trypsin Inhibitor Unit (TIU) will decrease the activity of 2 trypsin units by 50%, where 1 trypsin unit will hydrolyze 1.0  $\mu$ mole of *N*- $\alpha$ -benzoyl-DL-arginine *p*-nitroanilide (BAPNA) per minute at pH 7.8 and 25 °C.

Another commonly used unit of activity is the KIU (Kallikrein Inhibitor Unit). Different conversion factors between aprotinin units have been reported:

- 1 TIU = 1,300 KIU
- 1 TIU = 1,025 KIU<sup>14</sup>

## Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

## Preparation Instructions

Aprotinin is freely soluble in water (5 mg/mL) and in aqueous buffers of low ionic strengths.<sup>15</sup>

## Storage/Stability

Store the lyophilized powder at 2-8 °C. When stored at 2-8 °C, the product retains activity for at least 2 years.

Dilute solutions of aprotinin are generally less stable than concentrated ones. Solution stability is pH-dependent, although a pH range of 1-12 can be tolerated.<sup>16</sup> Repeated freeze-thaw cycles should be avoided.

The Cys<sup>14</sup>-Cys<sup>38</sup> disulfide bridge is readily split by reducing agents like 2-mercaptoethanol.<sup>16</sup> Due to its compact tertiary structure, aprotinin is relatively stable against denaturation due to high temperature, organic solvents, or proteolytic degradation (See Table 2). Only thermolysin has been found capable of degrading aprotinin after heating to 60-80 °C.<sup>16</sup>

The high basicity of aprotinin causes it to adhere to commonly used dialysis tubing and even gel filtration matrices. However, the use of acetylated materials and concentrated salt solutions (such as ≥ 0.1 M NaCl in buffer) minimizes this problem.<sup>16</sup>

Sterilization may be achieved by filtration through a 0.2 µm filter.

## References

1. Malysa, Agnes, "The Role of Spy1 in Mammary Morphogenesis". University of Windsor, M.Sci. thesis, p. 26 (2010).
2. Winkler, Claudia, "Ca<sup>2+</sup>- and Na<sup>+</sup>-regulating proteins in Nuclei Isolated from Human Myocardium". Karl-Franzens-University Graz, M.Sci. thesis, p. 18 (2010).
3. Flemmer, Robert T., "The anti-apoptotic role of Mcl-1 during early nervous system development". Memorial University of Newfoundland, M.Sc. thesis, p. 42 (2019).
4. Glass, Lara Joelle, "Brain invaders: IgG and IgG receptors in schizophrenia". University of New South Wales, M.Sc. thesis, p. 31 (2019).
5. Scott, Monique, "Endosomal SNX Proteins are Required for Recycling of the Epithelial Sodium Channel". University of Otago, M.Sc. thesis, p. 44 (2019).
6. Nock, Gemma, "Investigating the isoform-specific role of phosphoinositide 3-kinase p110delta". University College London, Ph.D. dissertation, p. 84 (2008).
7. Biteau, Flore, "Production de protéines recombinants par des plantes carnivores génétiquement transformées: Application à *Drosera rotundifolia* et transfert de la technologie à *Nepenthes alata*" ("Production of recombinant proteins by genetically transformed carnivorous plants: Application to *Drosera rotundifolia* and technology transfer to *Nepenthes alata*"). Université de Lorraine, Ph.D. dissertation, p. 55 (2009).
8. Schönleitner, Patrick, "Expression and function of Na<sup>+</sup> handling proteins in young β1-adrenergic receptor transgenic mice at an early stage of cardiac remodeling". Medical University of Graz, Dr. med. univ. dissertation, p. 24 (2013).
9. Rodriguez, Edwin Hadalid, "A method for annotating protein-coding sequences from ribosome profiling data". University of California San Francisco, Ph.D. dissertation, p. 19 (2014).
10. Clara, Rosmarie, "Intestinal cellular bioenergetics, glucagon-like peptide-1 and control of eating". ETH Zürich, Dr. sc. dissertation, pp. 74, 75 (2015).
11. Ugras, Scott Edward, "Proteostasis Responses to Endogenous Alpha-Synuclein Aggregation in the Brain". University of Pennsylvania, Ph.D. dissertation, p. 59 (2016).
12. Huang, Kuan-Lin, "Multi-omics Portraits of Cancer". Washington University in St. Louis, Ph.D. dissertation, p. 86 (2018).
13. Chamorro, María Cristina Casero, "Endolithic life in the Atacama desert: microbial ecology and adaptation strategies through a multidisciplinary approach". Universidad Autónoma de Madrid, Ph.D. dissertation, pp. 165, 166 (2019).
14. Hewlett, G., *Biotechnology (NY)*, 565-566 (1990).
15. *Biochemica Information*, 1st ed. (J. Keesey, ed.). Boehringer Mannheim Biochemicals (Indianapolis, IN), p. 111 (1987).
16. Fritz, H., and Wunderer, G., *Arzneimittelforschung*, **33(4)**, 479-494 (1983).
17. Zollner, H. (ed.), *Handbook of Enzyme Inhibitors*, 2nd ed., Part B. VCH Verlagsgesellschaft (Weinheim, Germany), p. 572 (1993).
18. Kunitz, M., and Northrop, J.H., *J. Gen. Physiol.*, **19(6)**, 991-1007 (1936).
19. Sherman, M.P., and Kassell, B., *Biochemistry*, **7(10)**, 3634-3641 (1968).
20. Zyznar, E.S., *Life Sci.*, **28(17)**, 1861-1866 (1981).

**Table 1. Inhibition by Aprotinin**

<b>Enzyme (Source), Condition</b>	<b>Inhibition (K<sub>i</sub>)</b>
Acrosin	Weak inhibition <sup>15</sup>
Chymotrypsin	K <sub>i</sub> = 9 nM <sup>17</sup>
CMP- <i>N</i> -Acetyl-neuraminatate lactosylceramide α-2,3-sialyltransferase	74% Inhibition at 300 nM <sup>17</sup>
Elastase (human leukocyte), pH 8.0	K <sub>i</sub> = 3.5 μM <sup>16</sup>
Kallikrein (pancreatic), pH 8.0	K <sub>i</sub> = 1.0 nM <sup>16</sup>
Kallikrein (plasma)	K <sub>i</sub> = 30 nM; 100 nM <sup>15</sup>
Kallikrein (tissue)	K <sub>i</sub> = 1 nM <sup>15</sup>
Kallikrein (urine)	K <sub>i</sub> = 1.7 nM <sup>15</sup>
Plasmin (porcine), pH 7.8	K <sub>i</sub> = 4.0 nM <sup>16</sup>
Plasminogen activator	K <sub>i</sub> = 8 μM; 27 μM <sup>15</sup>
Trypsin (bovine), pH 8.0	K <sub>i</sub> = 0.06 pM <sup>16</sup>
Trypsinogen (bovine), pH 8.0	K <sub>i</sub> = 1.8 μM <sup>16</sup>
Tryptase TL-2	16% Inhibition at 10 μM <sup>15</sup>
Urokinase (human), pH 8.8	K <sub>i</sub> = 8.0 μM <sup>16</sup>

**Table 2. Aprotinin solution stability**

<b>Solvent</b>	<b>Concentration</b>	<b>Temperature</b>	<b>% Loss/Time</b>
Sterile water with 0.9% NaCl and 0.9% benzyl alcohol, pH 5.7-6.2	10 mg/mL	0-5 °C	< 4.3%/year
2.5% Trichloroacetic acid	N/A	80 °C	No loss <sup>18</sup>
pH < 12.6	N/A	N/A	No loss observed after 24 hours <sup>19</sup>
pH < 12	N/A	N/A	Irreversibly denatured <sup>20</sup>
pH 7-8	0.065-1.95 μg/mL	4 °C	About 1 week <sup>15</sup>
pH 7-8	0.065-1.95 μg/mL	-20 °C	> 6 months <sup>15</sup>

---

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

### Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

### Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

### Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.  
© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

A6103dat Rev 06/22 RBG,MAC,GCY,MAM

The logo for MilliporeSigma, featuring the word "Millipore" in a red, sans-serif font above the word "SIGMA" in a bold, red, sans-serif font.