

3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Product Information

Shields and Sang M3 Insect Media

Insect tissue culture media have been formulated to mimic the main properties of the hemalymph of specific insects. A survey of the formulas of these media reveals significant differences in their composition.¹ Different media designed for the same species may differ more than two media designed for insects of different orders. Various media have been devised for the *in vitro* culture of *Drosophila* cells and tissues.^{2,3} The most extensively used are Echalier and Ohanessian's D-22 medium, Schneider's medium, and Shields and Sang M3 medium. *Drosophila* cells have been employed to study a variety of biological processes⁴ including genetics, endocrinology, physiology, and cell biology as well as recombinant protein expression⁵.

When supplemented with 5–20% heat-inactivated fetal bovine serum, Shields and Sang M3 medium has been found to support the rapid growth of both primary and established cultures of *Drosophila melanogaster* cells.^{6,7} The medium has been used for the growth and maintenance of the cell lines originally derived by Sang and coworkers from *Drosophila* embryos and cell lines from other laboratories (e.g. Schneider's lines 1-3, and Echalier and Ohanessian's K line and its derivatives).

	S8398	S3652
	[Powder]	[1×]
COMPONENT	g/L	g/L
Inorganic Salts		
CaCl ₂ (anhydrous)	0.760274	0.760274
MgSO₄ (anhydrous)	2.149128	2.149128
KHCO₃		0.5
NaH ₂ PO ₄	0.88	0.88
Amino acids		
β-Alanine	0.25	0.25
L-Alanine	1.5	1.5
L-Arginine	0.5	0.5
L-Asparagine (anhydrous)	0.3	0.3
L-Aspartic Acid	0.3	0.3
L-Cystine · HCI	0.02	0.02
L-Glutamic Acid · K	7.88	7.88
L-Glutamic Acid · Na	6.53	6.53
L-Glutamine	0.6	0.6
Glycine	0.5	0.5
L-Histidine	0.55	0.55
L-Isoleucine	0.25	0.25
L-Leucine	0.4	0.4
L-Lysine · HCl	0.85	0.85
L-Methionine	0.25	0.25
L-Phenylalanine	0.25	0.25
L-Proline	0.4	0.4
L-Serine	0.35	0.35
L-Threonine	0.5	0.5
L-Tryptophan	0.1	0.1
L-Tyrosine · 2Na	0.3601	0.3601
L-Valine	0.4	0.4

Vitamins and others		
Choline Chloride	0.05	0.05
Oxalacetic Acid	0.25	0.25
D(+)-Glucose	10.0	10.0
BIS-TRIS	1.05	1.05
Yeast Extract	1.0	1.0
Add		
KHCO₃	0.5	_

References

- Mitsuhashi, J., Media for insect cell cultures. in: Advances in Cell Culture vol. 2, Maramorosch, K., ed., Academic Press (New York, NY: 1982) pp 133-196.
- Echalier, G., *In vitro* culture of *Drosophila* cells and applications in physiological genetics. in: Invertebrate Tissue Culture, Kurstak, E., and Maramorosch, K., eds., Academic Press (New York, NY: 1976) pp. 131-150.
- Kuroda, Y., *Drosophila* tissue cultures: Retrospect and prospect. in: Invertebrate Cell Culture Appliations, Maramorosch, K., and Mitsuhashi, J., eds., Academics Press (New York, NY: 1982) pp. 53-104.
- 4. Schneider, I., and Blumenthal, A., in: Biology and Genetics of *Drosophila* vol. 2B, Ashburner, M., and Wright, T.R.F., eds., Academic Press (New York, NY: 1978) pp 265-315.
- Van der Straten, A. et al., Efficient expression of foreign genes in cultured *Drosophila melanogaster* cells using hygromycin B. selection. in: Invertebrate Cell System Applications Vo. 1, Mitsuhashi, J., ed., CRC Press (Boca Raton, FL: 1989) pp. 183-195.
- 6. Shields, G., and Sang, J.H., Improved medium for culture of *Drosophila* embryonic cells. *Drosophila* Information Service, **52**, 161 (1977).
- Shields, G., Dubendorfer, A., and Sang, J.H. Differentiation *in vitro* of larval cells from early embryonic cells of *Drosophila melanogaster*. J. Embryol. Exp. Morphol., 33, 159-175 (1975).

VNC, JG, PD, JF, ALF, MAM 03/15-1

©2015 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.