

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation


HYGROMYCIN B

from Steptomyces hygroscopicus

Product Numbers **H0654**, **H3274**, **H5527**, **H7772** and **H9773**

CAS #: 31282-04-09

Product Description



H7772, H3274, H9773: Lyophilized powder H0654: Aqueous solution 0.2 micron filtered H5527: Aqueous solution, gamma-irradiated

Molecular formula: C₂₀H₃₇N₃O₁₃

Formula weight: 527.54

Melting point: 160-180°C (with decomposition).

pK_a: 7.1 and 8.8.

Hygromycin B is an antibiotic substance isolated from *Streptomyces hygroscopicus*. Its mode of action is the inhibition of protein synthesis, by inducing the misreading of the mRNA template in the prokaryote (e.g. E. coli, at $100 \, \mu g/ml$), lower eukaryotes (e.g., yeast, at $200 \, \mu g/ml$) and higher eukaryotes (e.g., mammalian cells in culture). It selectively penetrates cells that have been rendered permeable by virus infection and combined with its ability to inhibit translation, it is an effective antiviral agent.

Hygromycin B has been used in a variety of cell culture applications to select drug-resistant stable transfectants after transfer of the hygromycin phosphotransferase gene. Specific references include transformation of Aspergillus species, Agaricus bisporus, Agrobacterium, Histoplasma capsulatum, and

Penicillium urticae. ¹¹ For use as a selection agent, the suggested concentration range is 100-800 μg/ml; specifically, 100 μg/ml for prokaryotes; 200 μg/ml for lower eukaryotes, and 150-400 μg/ml for higher eukaryotes.

For protocols that describe usage of Hygromycin B in units, 1 X 10⁶ units is approximately equal to 900 mg H7772, but each lot should be titered for each application.

Preparation Instructions

Hygromycin B is soluble in water at >50 mg/mL. ² It is also soluble in methanol or ethanol. Solutions should be sterilized by filtration, not by autoclaving. Hygromycin B is dissolved in water at 50 mg/mL, then sterile-filtered to produce H 0654, or γ -irradiated to produce H 5527.

Hygromycin B solutions have been reported to lose activity on freezing. Since solutions are stable refrigerated, freezing should be avoided.

Storage/Stability

Hygromycin B products should be stored as supplied at $2-8^{\circ}$ C. The dry solid is stable for at least five years if stored at $2-8^{\circ}$ C. Hygromcyin B solutions (Prod. Nos. H0654 and H5527) are stable as supplied for two years if stored at $2-8^{\circ}$ C.

References

- . Gonzalez, A. et al., *Biochim. Biophys. Acta*, 521, 459 (1978).
- 2. Rastinejad, F. et al., Cell, 75, 1107-1117 (1993).
- 3. Gritz, L. and Davies, J., Gene, 25, 179 (1983).

- 4. Malpartida, F. et al., *Biochem. Biophys. Res. Commun.*, 177, 6 (1983).
- 5. Giordano, T.J. and McAllister, W.T., *Gene*, 88, 285-288 (1990).
- 6. Wendt, S. et al., Curr. Genet., 17, 21-24 (1990).
- 7. Cullen, D. et al., Gene, 57, 21-26 (1987).

- 8. Van de Rhee, M.D. et al., *Mol. Gen. Genet.*, 250, 252-258 (1996).
- 9. Perl, A. et al., Nat. Biotechnol., 14, 624-628 (1996).
- 10. Woods, J.P. et al, *Infect. Immun.*, 66, 1697-1707 (1998).
- 11. Kiuchi, N. et al., *Agric. Biol. Chem.*, 55, 3053-3057 (1991).

CKV 01/98, revised DMG 12/02