
M Air T™ Calibration Kit User Guide

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Notice

The information in this document is subject to change without notice and should not be construed as a commitment by Millipore Corporation. Millipore Corporation assumes no responsibility for any errors that may appear in this document. This manual is believed to be complete and accurate at the time of publication. In no event shall Millipore Corporation be liable for incidental or consequential damages in connection with or arising from the use of this manual.

Copyright 2000, Millipore Corporation. All rights reserved. This book or parts thereof may not be reproduced in any form without the written permission of the publishers.

Millipore is a registered trademark of Millipore Corporation or an affiliated company.

M Air T is a trademark of Millipore Corporation or an affiliated company.

Milli-Q is a registered trademark of Millipore Corporation or an affiliated company.

Contents

Introduction	1
Calibration Kit Components	1
Check and Calibration	2
Working with Multiple Tester Sieves	4
Atmospheric Conditions	5
Flow Rate Curve	6
Calculating the Flow Rate	6
Calculating the Float Height	7
Calculation Using the Calibration Certificate	9
Interpreting the Results	11
Cleaning Procedure	11
Calibrating the M Air T Tester	12
Determining the Theoretical Height	12
Calibration Method	14
Technical Assistance	17
Standard Warranty	18

Notice

This procedure allows you to check and calibrate your M Air T Air Tester with an accuracy of +/-10% provided your measuring instruments have a resolution of:

- +/-1.1 mbar for the pressure measurement
- +/-0.5°C for the temperature measurement

At your request, we can perform a calibration accurate to +/-5%.

Introduction

The M Air T™ Calibration Kit allows you to check and calibrate the Millipore M Air T Tester with an accuracy of $\pm 5\%$. Millipore calibrates the calibration kit before delivery with an accuracy of $\pm 2.5\%$.


NOTE: Millipore recommends that you send the calibration kit at least once a year to our closest subsidiary to calibrate it against a reference flowmeter.

Calibration Kit Components

The calibration kit includes:

- A float type flowmeter
- A calibration sieve
- A calibration certificate (including conversion tables plus the table and flow rate curve equation specific to your kit)
- A calibration key


NOTE: Optional calculation software can be provided by our technical department.

Check and Calibration

Use the following procedure to check the calibration status of the M Air T Air Tester with a micro-perforated tester sieve.


1. Measure the atmospheric pressure (P in mbars) and temperature (T in °C) at the time of the test.
2. Insert a prefilled cassette into the M Air T Tester head.
3. Install the micro-perforated tester sieve on the head and lock it.
4. Press the On/Off button.
5. Make sure the unit indicates 1000 liters. (If not, refer to the User Guide supplied with the M Air T Air Tester and adjust the value to 1000 liters.)
6. Lock the M Air T Air Tester onto the Calibration kit as shown below:

Check and Calibration, continued

7. Start the Tester by pressing the Start/Delay twice in rapid succession. (Pressing the Start / Delay button once activates the count-down function.)
8. During suction, record the float height in mm displayed on the front of the flowmeter float. Sampling takes place at two different flow rates. Note the height during the first 500 liters (nominal flow rate: 139.5 L/min), then during the last 500 liters (nominal flow rate: 181.8 L/min).

NOTE: Take the reading on the top of the float. If the float oscillates during the measurement, note the average height.

9. Convert each height (mm) into its equivalent volumetric flow rate. To do so, use:
 - The equations provided in the following sections, or
 - The conversion tables supplied with the calibration certificate, or
 - The conversion software (supplied on request).

NOTE: This information is specific to the serial number of each calibration kit.

Working with Multiple Tester Sieves

If you are using more than one tester sieve, you can save time and effort by proving that all the sieves you are using are within calibration values, regardless of the sampler you are using.

First, you must check the specifications of your M Air T sampler. Do this by following the procedure outlined in the previous “Check and Calibration” section, but instead of using a tester sieve, use the calibration sieve. Complete this procedure for each sampler that you are using. If your flow rates are not within acceptable limits (see the “Interpreting the Results” section) for one or more of your samplers, you must calibrate the sampler(s) following the steps outlined in the “Calibrating the M Air T Tester” section.

Then, once all your samplers have been tested and found to meet calibration specifications, you can test each of your sieves. Using one of the samplers that has already been calibrated, follow the procedure outlined in the previous “Check and Calibration” section, this time using a tester sieve, as directed. Repeat this procedure with each of your tester sieves, using the same M Air T sampler each time.

When you’ve completed the tests for each of the tester sieves, you can be assured of staying within calibration values while using any of your samplers with any of your sieves.

Atmospheric Conditions (applied)

The standard Millipore atmospheric conditions for this unit are:

- Atmospheric pressure: 1000 mbars
- Temperature: 20°C

If your environmental conditions (atmospheric pressure and temperature) are different, use the following equation to obtain the flow rate corresponding to your conditions:

$$Q_0 = Q_{\text{Ref}} \sqrt{\frac{D_{\text{Ref}}}{D_0}}$$

Where:

Q_0 : Volumetric flow rate under your measurement conditions

Q_{Ref} : Volumetric flow rate at the standard pressure and temperature (1000 mbars , 20°C)

D_{Ref} : 1.18883 kg/m³ = Air density at the standard pressure and temperature (1000 mbars, 20°C)

D_0 : Air density under your measurement conditions

The air density is calculated either from the tables supplied with the calibration certificate (see the “Calculation Using the Calibration Certificate” section), or by applying the approximated equation below:

$$\text{Density of dry air} = \frac{\text{Pressure (mbars)}}{1013.25} \times \frac{1.293}{1 + (0.00367 \times \text{Temperature (}^\circ\text{C)})}$$

The pressure and temperature in this equation are those measured under your environmental conditions.

Flow Rate Curve (at standard temperature and pressure)

The equation of the flow rate curve (in liters per minute) versus the flowmeter flow height (in mm) established at standard Millipore temperature and pressure is:

$$Q_{\text{Ref}} = Ah^2 + Bh + C$$

Where:

Q_{Ref} : Volumetric flow rate (in liters per minute) at the standard pressure and temperature (1000 mbars, 20°C)

h: Flowmeter float height (in mm)

A, B and C: Constants specific to each calibration kit

Calculating the Flow Rate (from the float height)

Below is an example of a typical calculation:

Your measurement conditions:

■ Atmospheric pressure: 1030 mbars

■ Temperature: 15°C

Determine the air density ratio:

$$\sqrt{\frac{D_{\text{Ref}}}{D_0}}$$

Where:

$$D_{\text{Ref}}: \frac{1000}{1013.25} \times \frac{1.293}{1 + (0.00367 \times 20 \text{ [}^\circ\text{C]})} = 1.18883$$

$$D_0: \frac{\text{Pressure (mbars)}}{1013.25} \times \frac{1.293}{1 + (0.00367 \times \text{Temperature [}^\circ\text{C]})}$$

I.e. under your measurement conditions:

$$\sqrt{\frac{D_{\text{Ref}}}{D_0}} = \sqrt{\frac{1.18883}{\frac{1030}{1013.25} \times \frac{1.293}{1 + (0.00367 \times 15)}}} = 0.9769$$

Determine the flow rate:

$$Q_0 = Q_{\text{Ref}} \sqrt{\frac{D_{\text{Ref}}}{D_0}}$$

Where:

$$Q_{\text{Ref}} = A \times h^2 + B \times h + C$$

For a float height of 30 mm, the volumetric flow rate under the measurement conditions is calculated as follows:

$$Q_{\text{Ref}} = A \times 30^2 + B \times 30 + C$$

$$Q_0 = (A \times 30^2 + B \times 30 + C) \times 0.9769$$

Calculating the Float Height (from the flow rate)

Below is an example of a typical calculation.

Your measurement conditions:

- Atmospheric pressure: 1030 mbars
- Temperature: 15°C

Determine the air density ratio:

$$\sqrt{\frac{D_{\text{Ref}}}{D_0}}$$

For the method, see above.

Calculate the equivalent flow rate Q_{Ref} at the standard Millipore pressure and temperature:

$$Q_0 = Q_{\text{Ref}} \sqrt{\frac{D_{\text{Ref}}}{D_0}}$$

$$\text{I.e. } Q_{\text{Ref}} = Q_0 \times (1/0.9769)$$

$$Q_{\text{Ref}} = Q_0 \times \sqrt{\frac{D_0}{D_{\text{Ref}}}}$$

For a flow rate Q_0 of 139.5 L/min (example)

$$\begin{aligned} Q_{\text{Ref}} &= 139.5 \times (1/0.9769) \\ &= 142.8 \text{ L/min} \end{aligned}$$

Calculate the float height:

$$Q_{\text{Ref}} = A \times h^2 + B \times h + C$$

Second-degree equation whose solution is:

$$\text{Height (mm) } h = \frac{-B + \sqrt{(B^2 - (4 \times A \times (C - Q_{\text{Ref}}))}}{2 \times A}$$

I.e. in our example:

$$Q_{\text{Ref}} = 142.8$$

Use the values of A, B and C indicated on your certificate.

Calculation Using the Calibration Certificate

Below is an example of a typical calculation using the calibration certificate.

Your measurement conditions:

- Atmospheric pressure: 1030 mbars
- Temperature: 15°C

Determine the air density ratio:

$$\sqrt{\frac{D_{\text{Ref}}}{D_0}}$$

You can obtain this value from the conversion table (Attachment 1) supplied with the calibration certificate. See the example below:

Temperature °C	Atmospheric pressure mbar	Air density D _o (Kg/m ³)	Square root of D _{Ref} /D _o
15	980	1,18532	1,00
15	985	1,19137	1,00
15	990	1,19741	1,00
15	995	1,20346	0,99
15	1000	1,20951	0,99
15	1005	1,21556	0,99
15	1010	1,22160	0,99
15	1015	1,22765	0,98
15	1020	1,23370	0,98
15	1025	1,23975	0,98
15	1030	1,24579	0,98
15	1035	1,25184	0,97

The value found here is 0.98.

Interpreting the Results

The flow rate values must be between:

- 125.5 and 153.5 L/min, i.e. 139.5 L/min +/-10% for the first speed
- 163.6 and 200 L/min, i.e. 181.8 L/min +/-10% for the second speed

If this is not the case, there are two possible causes:

- The sieve is partially clogged and should be cleaned (see procedure below)
- The tester needs to be recalibrated (See the “Calibration Method” section for an outline of the procedure.)

In all cases, clean the sieve before performing a new calibration.

Cleaning Procedure

1. Immerse the sieve in an ultrasonic bath containing 50% tap water + 50% ultrapure water (Milli-Q®).
2. Run the ultrasonic bath for one hour.
3. Remove the sieve from the bath and dry it with filtered compressed air.

When the sieve has been cleaned, retest it by the above procedure.

Calibrating the M Air T Tester

Before calibrating the M Air T Tester, you must determine the correct float height at two specific flow rates under your atmospheric conditions.

Determining the Theoretical Height (at both flow rates)

Measure the atmospheric pressure and temperature before performing the calibration:

For example, your measurement conditions may be:

- Atmospheric pressure: 980 mbars
- Temperature: 19°C

Determine the air density ratio:

$$\sqrt{\frac{D_{\text{Ref}}}{D_0}}$$

Find this value on the conversion table (Attachment 1) of the calibration certificate.

Temperature °C	Atmospheric pressure mbar	Air density D ₀ (Kg/m ³)	Square root of D _{Ref} /D ₀
19	980	1,16905	1,01
19	985	1,17502	1,01
19	990	1,18098	1,00
19	995	1,18695	1,00
19	1000	1,19291	1,00
19	1005	1,19887	1,00
19	1010	1,20484	0,99
19	1015	1,21080	0,99
19	1020	1,21677	0,99
19	1025	1,22273	0,99
19	1030	1,22870	0,98
19	1035	1,23466	0,98

The value found here is 1.01

Determine the theoretical height (at both flow rates):


To determine the float height corresponding to 139.5 L/min (nominal flow rate corresponding to first speed), take the approximate value of 139.7 from the flow rate scale table (see Attachment 2).

Float height at the calibration unit, mm	Flow rate l/min						
	Square root	Square root	Square root	Square root	Square root	Square root	Square root
	$D_{Ref}/D_o=$ 0,97	$D_{Ref}/D_o=$ 0,98	$D_{Ref}/D_o=$ 0,99	$D_{Ref}/D_o=$ 1,00	$D_{Ref}/D_o=$ 1,01	$D_{Ref}/D_o=$ 1,02	$D_{Ref}/D_o=$ 1,03
##	130,5	131,8	133,2	134,5	135,9	137,2	138,6
##	131,2	132,6	133,9	135,3	136,6	138,0	139,3
##	131,9	133,3	134,7	136,0	137,4	138,7	140,1
##	132,7	134,0	135,4	136,8	138,1	139,5	140,9
##	133,4	134,8	136,1	137,5	138,9	140,3	141,6
##	134,1	135,5	136,9	138,3	139,7	141,0	142,4
##	134,9	136,3	137,6	139,0	140,4	141,8	143,2
##	135,6	137,0	138,4	139,8	141,2	142,6	144,0
##	136,4	137,8	139,2	140,6	142,0	143,4	144,8
##	137,1	138,5	139,9	141,4	142,8	144,2	145,6
##	137,9	139,3	140,7	142,1	143,6	145,0	146,4

NOTE: This is the height which should be obtained when calibrating the M Air T Tester for the first speed.


Calibration Method

1. Insert a prefilled cassette into the M Air T Tester head.
2. Replace the M Air T Tester sieve with the calibration sieve supplied with the calibration kit.
3. Lock the M Air T Tester onto the calibration kit as shown:

Calibration Method, continued

4. Unscrew the M Air T Tester rear panel using the allen key provided with the M Air T tester kit.
5. Switch the Tester into calibration mode by replacing the protective cap with the calibration key (supplied with the kit) as shown:

6. Screw on the rear panel and press the On/Off button.
7. Record the value displayed on the Tester screen (this value corresponds to the first speed currently stored). (The setting for the first speed ranges from 3.00 to 4.00 by increments of 0.01).
8. Press the Litre button to adjust the flow rate until the flowmeter float is in front of the mm graduation corresponding to a flow rate of 139.5 L/min.

NOTE: Holding down the Litre button when passing the 10s digit causes the increment to be set to 0.1.

Litre button held down continuously		
3.57 ➡ 3.58 ➡ 3.59 ➡ 3.60	Then	3.70 ➡ 3.80 ➡ 3.90 ➡ 4.00
Releasing the Litre button when passing the 10s digit causes the increment to be set to 0.01.		
Litre button held down	Released	Held down
3.57 ➡ 3.58 ➡ 3.59 ➡ 3.60	Then 3.70 ➡ 3.80	3.81 ➡ 3.82 ➡ 3.83 ➡ ➡ ➡

Calibration Method, *continued*

9. When the float is correctly aligned, press the Start/Delay button to verify and record the reading.
10. Then press the Start/Delay button again to switch the Tester to the second speed
11. Record the value displayed on the Tester screen (this value corresponds to the second speed currently stored).

NOTE: The setting for this first speed ranges from 4.00 to 5.00 by increments of 0.01.

12. Press the Litre button (see explanation above) to adjust the flow rate until the flowmeter float is in front of the mm graduation corresponding to a flow rate of 181.8 L/min.
13. When the float is correctly aligned, press the Start/Delay button to validate and record the reading.
14. Switch off the Tester by pressing the On/Off button.
15. Unscrew the M Air T Tester rear panel. Then switch the tester into Sampling mode by replacing the calibration key with the protective cap. Screw the rear panel back onto the Tester.
16. Check the M Air T tester flow rate and the micro-perforated tester sieve as described in the previous sections of this manual.

NOTE: The old calibration values can be restored until the new values have been validated by pressing the Start/Delay button. To restore the old values, switch the Tester off and back on by the On/Off button.

Technical Assistance

For more information, contact the Millipore office nearest you. In the U.S., call **1-800-MILLIPORE** (1-800-645-5476). Outside the U.S., see your Millipore laboratory catalogue for the phone number of the office nearest you. You can reach us by e-mail at tech_service@millipore.com or visit our web site (www.millipore.com).


Millipore Corporation is pleased to provide internet access to Material Safety Data Sheets (MSDS) for its products that contain hazardous materials. To obtain any MSDS documents that may be associated with this product, go to the MSDS page of our website (www.millipore.com/msds.nsf/home).

Standard Warranty

Millipore Corporation ("Millipore") warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. **MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of Millipore products appearing in Millipore's published catalogues and product literature may not be altered except by express written agreement signed by an officer of Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, Millipore's sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies Millipore promptly of any such breach. If after exercising reasonable efforts, Millipore is unable to repair or replace the product or part, then Millipore shall refund to the customer all monies paid for such applicable product or part. **MILLIPORE SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER INDIRECT DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY CUSTOMER FROM THE USE OF ITS PRODUCTS.**

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