



# Phosphor-I $\kappa$ B (Ser32) STAR ELISA Kit

Sufficient reagents for 96 assays per kit

Catalog No. 17-486

FOR RESEARCH USE ONLY  
Not for use in diagnostic procedures.

USA & Canada Phone: +1(800) 645-5476

In Europe, please contact Customer Service:

France: 0825.045.645;

Spain: 901.516.645 Option 1

Germany: 01805.045.645

Italy: 848.845.645


United Kingdom: 0870.900.46.45

For other locations across the world please visit [www.millipore.com/offices](http://www.millipore.com/offices)


## Test Principles and Background

The UPSTATE<sup>®</sup> colorimetric STAR (Signal Transduction Assay Reaction) ELISA kit is a solid phase sandwich enzyme linked immunosorbent assay that provides a fast, sensitive method to detect specific levels of signaling targets in whole cell extracts. The I $\kappa$ B plate is coated with a specific mouse monoclonal I $\kappa$ B capture antibody on the microwells of the 96-well clear plate. Sample lysate or the standard included in the kit are incubated in the microwells allowing I $\kappa$ B antigen to be captured in the plate wells. The plate is then washed to remove any un-bound nonspecific material. A specific rabbit anti-phospho-I $\kappa$ B antibody is added to detect the captured I $\kappa$ B on the plate well that is phosphorylated on Ser32. The unbound detection antibody is washed away followed by incubation with an HRP-conjugated anti-rabbit antibody. This allows for a sensitive enzymatic detection of the sample. After the addition of TMB substrate and stop solution the absorbance is measured at 450 nm using a plate reader.

The entire assay takes less than 5 hours to complete with minimal hands-on time. Many of the reagents are supplied in ready-to use formulations for ease of use. The kit also includes a standard that is run as both a positive control and to generate a standard curve for I $\kappa$ B measurement.



NFκB signaling is a critical regulator not only of immune function, but also of proliferation versus apoptosis in response to various stimuli. In a majority of unstimulated cells, the NFκB transcription factors exist in their inactive form and are retained in the cytoplasm by the bound inhibitory IκB proteins that prevent it from entering the nucleus, thus keeping NFκB in an inactive state. Upon stimulation by multiple inducers including viruses or cytokines, such as TNFα, IL-1, or PMA, IκBα is rapidly phosphorylated by IKK (IκB Kinase), which phosphorylates IκB on serines 32 and 36, allowing recognition by the ubiquitin ligase machinery, leading to its polyubiquitination and degradation, resulting in the release of the NFκB complex. Once IκB is degraded, NFκB is able to initiate transcription. The NFκB complex then translocates to the nucleus where it activates gene transcription. NFκB induces the transcription of its own inhibitor, IκBα, causing an autoregulatory mechanism of NFκB activity and generating the inactive form of NFκB. The newly formed nuclear NFκB- IκBα complexes are then exported out to the cytoplasm, thereby reestablishing the cytoplasmic pool of inactive NFκB complexes primed for another round of activation to take place. The wide variety of genes regulated by NFκB includes those encoding cytokines, chemokines, adhesion molecules, acute phase proteins, and inducible effector enzymes. Mutant IκB in which serines 32 and 36 are changed to alanines, is not phosphorylated, and therefore not degraded. Cells expressing this protein are not able to activate NFκB, providing a useful tool to study the role of NFκB in various pathways and processes.



For Research Use Only.

---

## Kit Components

- Capture Plate pre-coated with anti-I $\kappa$ B antibody: (Part No. 17-486A) One pre-coated 96-stripwell immunoplate sealed in a foil pouch.
- Anti-phospho-I $\kappa$ B (Ser32) detection antibody: (Part No. 17-486B) One bottle (6 mL) of anti-phospho-I $\kappa$ B (Ser32) detection antibody containing sodium azide, ready to use.
- ELISA Diluent: (Part No. 17-486C) One bottle (25 mL) of ELISA Diluent containing sodium azide, ready to use.
- 25X ELISA Wash Buffer: (Part No. 17-486D) One bottle (50 mL) of 25X ELISA Wash Buffer.
- Anti-Rabbit IgG HRP conjugate: (Part No. 17-486E) One vial (125  $\mu$ L) of 100X anti-rabbit HRP conjugate.
- HRP Diluent: (Part No. 17-486F) One bottle (25 mL) of HRP Diluent.
- TMB Solution: (Part No. 17-486G) One bottle (25 mL) of stabilized tetramethylbenzidine (TMB), ready to use.
- Stop Solution: (Part No. 17-486H) One bottle (25 mL) of stop solution, ready to use.
- Phosphorylated I $\kappa$ B (Ser32) Standard: (Part No. 17-486I) Two vials of phosphorylated I $\kappa$ B standard (Ser32), lyophilized.
- Plate Covers: Two plate covers.

---

## Materials Not Supplied

- Microplate reader able to measure absorbance at 450 nm
- Adjustable pipettes to measure volumes ranging from 1  $\mu$ L to 1,000  $\mu$ L
- Deionized water
- Wash bottle or automated microplate washer
- Log-Log graph paper or software for data analysis
- Tubes to prepare standard dilutions
- Timer
- Plate Shaker
- Polypropylene vials

---

## Storage and stability

Store kit and unopened components at 2°C-8°C. Do not use kit beyond its expiration date.

---

## Reagent Precautions

- **Safety Warnings and Precautions:** This kit is designed for research use only and is not recommended for internal use in humans or animals. All chemicals should be considered potentially hazardous and principles of good laboratory practice should be followed.
- This kit is not for use in diagnostic procedures.
- **Caustic Material: Stop Solution.** The stop solution contains 2N sulfuric acid which is harmful if swallowed or inhaled; avoid contact with skin and eyes; wash areas of contact immediately with water. **Caution: Eye, hand, face, and clothing protection should be worn when handling this material.**
- Substrate Solution F is harmful if inhaled or ingested. Additionally, avoid skin, eye or clothing contact with the substrate reagents
- The instructions provided have been designed to optimize the kit's performance. Deviation from the instructions may result in suboptimal performance of the kit and a failure to produce accurate data.

---

## Specimen Collection and Handling

- All kit reagents should be brought to room temperature (20°C to 25°C) prior to use.
- Do not use reagents beyond the expiration date of the kit.
- Do not mix or interchange reagent from various kit lots.
- **Manual Plate Washing:** Vigorous washing and complete removal of all liquid by aspiration at the end of each washing step is very important to obtain low background values. Gentle agitation during the wash steps or a 2-3 minute soak may reduce background values.
- The desiccant enclosed in the 96-well capture plate pouch will keep the plate stable when stored at 2° to 8°C should the pouch lose its seal during shipping.

---

## Assay Sensitivity, detection limits and species reactivity

Sensitivity: 1 Unit/mL

Range of Detection: 1.6 to 100 Units/mL

Species Reactivity: Human, mouse and rat

**NOTE:** This data is presented for reference use only and should not be used to interpret actual assay results. A standard curve must be generated for each assay.

---

## Preparation of sample

- Stimulate cultured cells to induce I $\kappa$ B activation.
- Remove culture media and wash cells twice ice-cold with 1X TBS (Tris Buffered Saline) or PBS (Phosphate Buffered Saline). Discard supernatant
- Add 1X RIPA with protease inhibitors to the cells (5-10 ml per 150 mm tissue culture plate). 10mL of 1X RIPA plus protease inhibitors can be prepared by adding 10  $\mu$ L of 1  $\mu$ g/ $\mu$ L Leupeptin, 10  $\mu$ L of 1  $\mu$ g/ $\mu$ L Aprotinin, 10  $\mu$ L of 1  $\mu$ g/ $\mu$ L Pepstatin, 100  $\mu$ L of 100mM PMSF and 1 mL of 10X RIPA (Cat. No. 20-188) to 8.97 mL of deionized water.
- Scrape cells from plate with a rubber policeman.
- Transfer cells in RIPA buffer to a microcentrifuge tube and incubate on ice for 15 minutes. Vortex tube for 10 seconds or sonicate briefly for 10 seconds.
  - Clarify lysate by centrifugation at 12,000 rpm for 10 minutes at 4°C in a microcentrifuge prior to use.
- The assay tolerates a maximum of 20  $\mu$ L per well of lysate prepared in 1X RIPA diluted in ELISA Diluent.
- Cell extract containing SDS must be diluted to 0.01% SDS using ELISA Diluent prior to use.
- Collect the supernatant and calculate protein concentration using a Bradford Assay or by densitometry.
- Keep samples cold and use immediately if not freezing away for later use. Samples can be frozen and stored at -80° C. It is suggested to use the cell lysate be immediately used or use within 6 months if storing at -80° C. Avoid repeated freeze thaws. Further information of lysate preparation protocols can be obtained at <http://www.millipore.com> Cell Lysate Extracts-General Protocols.

## Reagent Preparation

### 1. 1X Wash Buffer

Bring the 25X ELISA Wash Buffer to room temperature and mix to ensure that any precipitated salts have re-dissolved. For 500 mL of Wash Buffer, combine 20 mL of 25X ELISA Wash Buffer and 480 mL distilled or deionized water. Stir to homogeneity. Wash Buffer can be stored for up to 4 weeks at 2-8°C. Discard the Wash Buffer if it becomes turbid or if a precipitate develops.

### 2. Anti-Rabbit IgG HRP Conjugate

Dilute the anti-Rabbit IgG HRP Conjugate 100-fold with HRP Diluent immediately before use. Prepare 1 mL for each strip of wells used.

### 3. Standard


**Note:** When opening lyophilized Standard, remove rubber stopper gently as the lyophilizate may have become dislodged during shipping.

Reconstitute the standard with the volume of ELISA Diluent specified on the vial label to give a concentration of 100 Units/mL. Gently swirl the vial and allow the vial to sit for 10 minutes to ensure the material is completely reconstituted. The standard should be reconstituted immediately before the assay. This stock material (tube #1) is then used to generate a standard curve. A suggested 2-fold dilution scheme is as follows:

a) Label 7 test tubes #2-7 and "0 dose". Add 250  $\mu$ L of the ELISA Diluent to tubes #2-7 and "0 dose".

b) Add 250  $\mu$ L of the stock Standard solution [100 Units/mL] to tube #2 and vortex. This is Standard tube #2 with a concentration of 50 Units/mL.

c) Standards #2-7 are then prepared by performing a 2-fold serial dilution of the preceding standard. Refer to Fig. 1. For example, to make Standard #3, remove 250  $\mu$ L of Standard #2 and add it to tube #3 and vortex and so on. Do not add any Standard to the "0 Dose" Standard tube.



**Figure 1: Suggested 2-fold Serial Dilution of Standard**

**Note:** The Standard curve can set up with a different serial dilution scheme by making appropriate adjustments to the dilution pattern

---

## Assay Protocol

1. Prepare the reagents as described in the Reagent Preparation section.
2. Place the desired number of strips in the strip well plate holder. (Re-bag the extra strips and return unused strips to refrigerator for future use.)
3. Add 100  $\mu$ L of either the Standards 1 through 7 or the samples to wells. Add 100  $\mu$ L of the zero dose to the control wells. It is recommended that standards and samples be run in duplicate. Wells reserved for TMB blanks should be left empty.

**Note:** Do not add standard or sample lysate to wells reserved for TMB blanks.

**Note:** A standard curve must be run at each setting.

4. Seal the plate with a plate sealer. Incubate the plate for 2 hours at room temperature or at 2-8°C overnight (on a shaker if possible).

### 5. IMPORTANT WASH STEP:

**Gently remove the plate sealer and wash the plate at least 4 times. A thorough washing of the plate is extremely important to reduce background. We recommend using a multi-channel pipette to fill each well with 250  $\mu$ L of diluted Wash Buffer. Fluid removal from the wells is best accomplished by inverting the plate over a sink and flicking the fluid out of the wells and then blotting the plate on clean paper towels. Using the multichannel pipet add 250  $\mu$ L of Wash Buffer to each well; flick and blot the plate. Repeat this procedure for a total of 4 times.**

*For users of automatic plate washers:* It is important to ensure that the wash apparatus is properly maintained and operating correctly. Tubing and tips can easily become clogged, leading to incomplete washing and inadequate aspiration of wells. The result may be poor precision and an unsuitable standard curve. For best results, we recommend at least 4 wash cycles.


6. Add 50  $\mu$ L of the detection antibody to each well. Cover the plate and incubate at room temperature for 1 hour (on shaker if possible). Wash as described in Step 5.
7. Add 100  $\mu$ L of a 1:100 dilution of the anti-Rabbit IgG HRP Conjugate to each well. Cover the plate and incubate at room temperature for 45 minutes (on shaker with mild agitation if possible). Wash as described in Step 5. Remove all fluid from the wells and blot the wells dry.
8. Add 100  $\mu$ L of the TMB Solution to each well. Incubate at room temperature in the dark for 10 to 45 minutes, monitor the color development. Stop the reaction by adding 100  $\mu$ L of Stop Solution to each well. Immediately read the plate at 450 nm. Plate should be read within 1 hour of adding the stop solution.
9. The plate reader may be blanked against a TMB blank prepared by adding 100  $\mu$ L of stop solution to 100  $\mu$ L of the TMB solution.

**CAUTION:** Bubbles in the wells will cause inaccurate readings. Ensure that all bubbles are removed prior to taking the absorbance reading. **NOTE:** For very low starting protein levels, samples can be placed at 37°C during the final incubation to obtain greater sensitivity.



---

## Calculations of Results:



**Figure 2. Typical phospho- IκB (Ser32) Standard Curve**

100 μL of progressive 2 fold dilutions of the phospho- IκB (Ser32) standard included in the kit and run as described in the assay instructions.

**NOTE:** This data is presented for reference use only and should not be used to interpret actual assay results. A standard curve must be generated for each assay.

---

## Reference

1. DiDonato, J.A., Mercurio, F., and Karin, M. (1995). *Mol. Cell. Biol.* **15**:1302-1311.
2. Huang, T and Miyamoto, S (2001). *Mol. Cell. Biol.* **21**:4737-4747.
3. Ghosh, S and Karin, M. (2002). *Cell* **109**:S81-S96.
4. Baeuerle, P.A. and Baltimore, D. (1988). *Cell* **53**:211-217.
5. Miyamoto, S., Chiao, P.J., and Verma, I.M. (1994). *Mol. Cell. Biol.* **14**:3276-3282.
6. Human PBMC ( $1 \times 10^6$  cells/mL) were stimulated with 100 ng/mL LPS at 37°C for overnight. Cell culture supernatants were collected and assayed for the levels of natural human IL-8. IL-8 concentration was 45,650pg/mL in LPS-stimulated samples and 160pg/mL in unstimulated samples.

---

## Warranty

**EMD Millipore Corporation** (“EMD Millipore”) warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. **EMD MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of EMD Millipore products appearing in EMD Millipore’s published catalogues and product literature may not be altered except by express written agreement signed by an officer of EMD Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, EMD Millipore Corporation’s sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies EMD Millipore Corporation promptly of any such breach. If after exercising reasonable efforts, EMD Millipore Corporation is unable to repair or replace the product or part, then EMD Millipore shall refund to the Company all monies paid for such applicable Product. **EMD MILLIPORE CORPORATION SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY COMPANY CUSTOMER FROM THE USE OF ITS PRODUCTS.**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

(c) 2009 - 2013: Merck KGaA, Darmstadt. All rights reserved. No part of these works may be reproduced in any form without permission in writing