#### sigma-aldrich.com

3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

# **Product Information**

N1-(7-Nitro-2,1,3-benzoxadiazol-4-yl)-N1,N2,N2-tris(2-pyridinylmethyl)-1,2-ethanediamine

Catalog Number **N1040** Storage Temperature –20 °C

#### CAS RN 1111625-98-9

Synonyms: *N*-(2-(Bis((pyridin-2-yl)methyl)amino)ethyl)-7-nitro-*N*-((pyridin-2-l)methyl)benzo[c][1,2,5]oxadiazole-4-amine, NBD-TPEA

# **Product Description**



Molecular formula: C<sub>26</sub>H<sub>24</sub>N<sub>8</sub>O<sub>3</sub> Molecular weight: 496.52

 $\begin{array}{l} \lambda_{\text{max}} : \mbox{480 nm} \\ \mbox{Fluorescence properties:} \\ \lambda_{\text{Ex}} : \mbox{470 nm} \\ \lambda_{\text{Em}} : \mbox{534 nm} \end{array}$ 

Zinc is of critical importance in many biological processes, such as gene expression, enzyme regulation, and neurotransmission, and there is considerable interest in new and more effective fluorescent probes for Zn<sup>2+</sup> detection.<sup>1</sup> Fluorescence measurement using chelators and ion probes provides a high level of sensitivity *in vivo* and *in vitro* for the determination of the presence and concentration of ions.

*N*1-(7-Nitro-2,1,3-benzoxadiazol-4-yl)-*N*1,*N*2,*N*2-tris (2-pyridinylmethyl)-1,2-ethanediamine has a high sensitivity for zinc ions. Using fluorescent detection, a linear responsive is observed for  $Zn^{2+}$  concentrations in the range of 1–10  $\mu$ M (see Figures 1 and 2). It has good selectivity for  $Zn^{2+}$  in the presence of competing metal ions, such as Ca<sup>2+</sup>, Mg<sup>2+</sup>, Fe<sup>2+</sup>, Fe<sup>3+</sup>, and Mn<sup>2+</sup> (see Figures 3 and 4).<sup>2</sup>

This product has been used *in vivo* to determine  $Zn^{2+}$  concentration in pancreatic  $\beta$  cells<sup>2</sup> and zebra fish larvae.<sup>3</sup>

# **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

# **Preparation Instructions**

This product is soluble at 0.5 mg/ml in DMSO.

# Storage/Stability

The product ships on wet ice and storage at –20  $^\circ\text{C}$  is recommended.

#### References

- Nolan, E.M., and Lippard, S.J, Small-molecule fluorescent sensors for investigating zinc metalloneurochemistry. *Acc. Chem. Res.*, 42, 193-203, (2009).
- Xu, Z., et al., An NBD- based colorimetric and fluorescent chemosensor for Zn<sup>2+</sup> and its use for the detection of intracellular zinc ions. *Tetrahedron*, 65, 2307-2312, (2009).
- Qian, F., et al., Visible light excitable Zn<sup>2+</sup> fluorescent sensor derived from an intermolecular charge transfer fluorophore and its in vitro and in vivo application, *J. Amer. Chem. Soc.*, **131**, 1460-1468, (2009).

PD,MS,VNC,MAM 05/10-1

#### Figure 1.

Fluorescence spectra of the product (Catalog Number N1040) in the presence of increasing zinc ion concentration



A 2.49 mM stock solution of the product was prepared in DMSO. This was diluted to 10  $\mu$ M in 0.1 HEPES buffer, pH 7.2. Aliquots of zinc standard were added to the final concentration shown. Excitation wavelength 470 nm, Excitation slit width 2.5 nm, Emission slit width 5 nm.





All metal ions at 30 mM except zinc ion concentration, which was 3 mM. Product concentration was 10 mM in 0.1 M HEPES buffer, pH 7.2.

Excitation wavelength 470 nm, Excitation slit width 2.5 nm, Emission slit width 5 nm.

#### Figure 3.

Ultraviolet absorption spectra of the product (Catalog Number N1040) in the presence of increasing zinc ion concentration



Final concentration of the product is 0.1 mM in HEPES buffer (1:9, DMSO/water, 100 mM KNO<sub>3</sub>; 0.1 M HEPES, pH 7.2), with aliquots of zinc nitrate added (100  $\mu$ M initial concentration).

#### Figure 4.

Ultraviolet absorption spectra of the product (Catalog Number N1040) in presence of various metal ions



Final concentration of the product is 0.1 mM in HEPES buffer (1:9, DMSO/water, 0.1 M HEPES, pH 7.2). Perchlorate salt solutions (10 mM) of  $Zn^{2+}$ ,  $K^+$ ,  $Ca^{2+}$ ,  $Fe^{2+}$ ,  $Na^+$ ,  $Mn^{2+}$ , and  $Fe^{3+}$  in demineralized water were prepared. 500 µl of product solution was added to 500 µl of demineralized water and an aliquot (40 µl) of the metal ion solution added.

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.