

Product Information

GUANOSINE 5'-TRIPHOSPHATE SODIUM Sigma Prod. No. G8752 and G8877

CAS NUMBER: 36051-31-7

SYNONYMS: GTP

PHYSICAL PROPERTIES:

Appearance: White, lyophilized powder

Formula: $C_{10}H_{13}N_5O_{14}P_3Na_3$


Formula weight: 523.2 (free acid)

$E^{mM}_{253nm} = 13.7$ (0.1 M Phosphate buffer at pH 7.0)

Abs (250/260) : 1.17

(280/260) : 0.66

pK_a values: 3.3, 6.3, 9.3¹



METHOD OF PREPARATION:

Guanosine 5'-triphosphate is prepared by enzymatic phosphorylation of guanosine 5'-monophosphate and purified by ion exchange chromatography.

DESCRIPTION:

GTP functions as a carrier of phosphates and pyrophosphates involved in channeling chemical energy into specific biosynthetic pathways. It also serves as an energy-rich precursor of mononucleotide units in the enzymatic biosynthesis of DNA and RNA.²

STABILITY / STORAGE AS SUPPLIED:

This product is routinely shipped on dry ice. Storage at -20°C is recommended.

Guanosine 5'-Triphosphate is a very unstable compound. It decomposes at least 2% within 6 months when stored at -0°C and at least 4% within one year stored at -0°C. 2% decomposition occurred within four days when stored at 37°C.

These data were obtained from samples stored at -0°C. We highly recommend storing this product at -20°C to achieve greater shelf-life.

GUANOSINE 5'-TRIPHOSPHATE SODIUM
Sigma Prod. No. G8752 and G8877

SOLUBILITY / SOLUTION STABILITY:

Sigma routinely tests the solubility at 100 mg/mL in water yielding a clear, colorless solution. It is best to prepare GTP solutions immediately before use. If a stock solution is required, we recommended preparing a 10 mg/mL solution in water and storing aliquots at -20°C. The stock solution should be kept 3-6 months at most.

REFERENCES:

1. *Data for Biochemical Research*, 3rd ed., 88-89 (1987).
2. Lehninger, A.L., *Biochemistry*, 2nd ed., 316 (1977).

Sigma warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.