


Product Information

Enterobactin from *Escherichia coli*

Catalog Number **E3910**
Storage Temperature $-20\text{ }^{\circ}\text{C}$

CAS RN 28384-96-5
Synonym: Enterochelin



Product Description

Molecular formula: $\text{C}_{30}\text{H}_{27}\text{N}_3\text{O}_{15}$
Molecular weight: 669.55

Iron mobilization and uptake by microbes is mediated by low molecular weight complexing agents named siderophores.¹ Enterobactin is a catechol [a benzene-diol, $\text{C}_6\text{H}_4(\text{OH})_2$] type siderophore produced in small quantities by *Escherichia coli* and related enteric bacteria when grown on iron deficient media,² and is the most powerful ferric ion complexing agent known.^{1,3}

Since it is highly hydrophobic, in order to act as a siderophore, enterobactin undergoes modifications by the *iroA* gene cluster inside the mammalian host before it is secreted.⁴ In *Escherichia coli*, enterobactin synthesis inhibition occurs through the binding of Fe^{2+} to Fur and Diphtheria toxin repressor proteins (DtxR).²

Enterobactin is a very effective sequestering agent for iron, forming an unusual macro-bridged hexacoordinate trianion. Therefore, it can easily remove iron from proteins, insoluble iron complexes, and other siderophores.^{3,6} Studies of the chemistry, regulation, synthesis, recognition, and transport of enterobactin make it the best-understood siderophore.⁵

Complexes of enterobactin with scandium (Sc^{3+}) and indium (In^{3+}) were shown to have antibacterial effect against *Klebsiella pneumoniae*, similar to that obtained with kanamycin sulfate. The Sc^{3+} -enterobactin complex was found to be active at $0.2\text{ }\mu\text{M}$ and appears to form an equilibrium mixture with Fe^{3+} -enterobactin complex.⁷ The In^{3+} -enterobactin complex does not produce complete bacteriostasis but rather a marked increase in generation time.

Purity: $\geq 98\%$ (HPLC)

Preparation instructions

Soluble at 10 mg/ml in DMSO or acetonitrile:water (9:1).

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Store the product sealed at $-20\text{ }^{\circ}\text{C}$. Under these conditions the product is stable for at least 2 years.

References

1. Ecker, D.J., et al., Recognition and transport of ferric enterobactin in *Escherichia coli*. *J. Bacteriol.*, **167**, 666-673 (1986).
2. Braun, V., and Braun, M., Iron transport and signaling in *Escherichia coli*. *FEBS Lett.*, **529**, 78-85 (2002).
3. Williams, P.H., and Carbonetti, N.H., Iron, siderophores, and the pursuit of virulence: independence of the aerobactin and enterochelin iron uptake systems in *Escherichia coli*. *Infect. Immun.*, **51**, 942-947 (1986).
4. Fischbach, M.A., et al., How pathogenic bacteria evade mammalian sabotage in the battle for iron. *Nat. Chem. Biol.*, **2**, 132-138 (2006).
5. Raymond, K.N., et al., Enterobactin: An archetype for microbial iron transport. *Proc. Natl. Acad. Sci. USA*, **100**, 3584-3588 (2003).

6. Hollifield, W.C. Jr., and Neilands, J.B., Ferric enterobactin transport system in *Escherichia coli* K-12. Extraction, assay, and specificity of the outer membrane receptor. *Biochemistry*, **17**, 1922-1928 (1978).
7. Rogers, H.J., et al., Antibacterial effect of scandium and indium complexes of enterochelin on *Klebsiella pneumoniae*. *Antimicrob. Agents Chemother.*, **18**, 63-68 (1980).

VNC,DWF,MAM 04/10-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.