

#### SAFC

3050 Spruce Street Saint Louis, MO 63103 USA Phone: 1-800-244-1173 • 1-314-534-4900

Fax: 1-888-289-8449 Global Email: safcglobal@sial.com

# **Product Information**

# **EX-CELL® CD CHO Fusion**

# Serum-Free Medium for CHO Cells Chemically Defined, Animal-Component Free

# without L-glutamine

# CATALOG NO. 14365C

#### Description

EX-CELL CD CHO Fusion is a chemically defined, animal-component free medium developed for the long-term growth of Chinese Hamster Ovary (CHO) cells. The absence of any large macro molecules allows for isolation and purification of secreted proteins from the cells. This medium is supplied without L-glutamine to aid in media stability, to avoid L-glutamine degradation that causes ammonia build-up and to provide an appropriate medium for the culture of CHO cells using glutamine synthetase selection (example: SAFC's CHOZN GS knockout cell line). This medium does not contain hypoxanthine or thymidine to allow for its use with dihydrofolate reductase(DHFR-) gene amplification systems (example: SAFC's CHOZN DHFR knockout cell line).

### **Formulation**

The formula for EX-CELL CD CHO Fusion is proprietary to SAFC. For additional information call our Technical Services department.

#### **Precautions**

Use aseptic technique when handling or supplementing this medium. This product is for research or for further manufacturing use. THIS PRODUCT IS NOT INTENDED FOR HUMAN OR THERAPEUTIC USE.

#### Storage

Store liquid medium at 2 to 8 °C, protected from light. Do not use after the expiration date.

# **Indications of Deterioration**

Medium should be clear and free of particulate and flocculent material. Do not use if liquid medium is cloudy or contains precipitates. Other evidence of deterioration may include color change, pH shift or degradation of physical or performance characteristics.

# **Preparation Instructions**

EX-CELL CD CHO Fusion is formulated with sodium bicarbonate and without L-glutamine. For applications requiring the use of L-glutamine, supplement prior to use with 4 - 8 mM L-glutamine by adding 20 - 40 mL/L of a 200 mM solution (Catalog No. 59202C). SAFC recommends L-glutamine supplementation of the working volume only. Supplements, such as antibiotics, can be added to the sterilized medium using aseptic technique. Storage conditions and shelf life of the product may be affected by the nature of the supplement.

### Methods for Use

The following procedure is suggested for cells coming from other formulations that are not EX-CELL CD CHO Fusion.

- 1. Subculture actively growing cells by planting new cultures at 4 x 10<sup>5</sup> cells/mL in 20 30 mL of EX-CELL CD CHO Fusion.
- Subculture stocks every three to four days at a seeding density of 4 x 10<sup>s</sup> cells/mL.
- Continue stocks for four to six passages until viabilities stabilize at >90%.
- Once cells are fully adapted to EX-CELL CD CHO Fusion seeding densities can be adjusted to lower densities for initiating new cultures.

## **Culture Techniques**

Once cultures are fully adapted, the cells should be passed every 3 - 4 days at a seeding density of at least 2-4 x 10<sup>5</sup> cells/mL. An optimal seeding density should be determined by the researcher for each application and cell type. When passing the cells, medium carryover should not exceed 25% of the final volume. If carryover exceeds 25%, centrifugation is recommended. Cells propagated in serum-free or protein-free media are extremely fragile. Standard techniques of centrifugation must be modified to include low-speed centrifugation to prevent damage to cells that have been propagated in serum-free medium.



#### SAFC

3050 Spruce Street Saint Louis, MO 63103 USA Phone: 1-800-244-1173 • 1-314-534-4900

1-888-289-8449

Global Email: safcglobal@sial.com

## Cryopreservation

# Freezing:

Cells can be frozen in EX-CELL® CD CHO Fusion without the reintroduction of serum.

- 1. Choose a culture in logarithmic growth with viabilities above 90%.
- 2. Prepare a freezing medium consisting of 45% cold EX-CELL CD CHO Fusion medium, 45% spent medium and 10% dimethyl sulfoxide (DMSO).
- 3. Centrifuge the cells at 200 g for 5 minutes. Remove the supernatant.
- 4. Resuspend the cells in the freezing medium at  $5 \times 10^6$  to 1 x 107 cells/mL.
- 5. Rapidly transfer 1 2 mL of this suspension to sterile cryovials.
- 6. Place the vials at -20 C for 3 4 hours, then transfer to -70 °C for 16 - 24 hours.
- For long-term storage, transfer the vials to liquid nitrogen vapor.

#### Thawing:

- 1. Rapidly thaw a vial of frozen cells in a 37 °C water bath.
- 2. Transfer the cells aseptically to a centrifuge tube containing 10 mL of chilled EX-CELL CD CHO Fusion medium.
- 3. Using low-speed centrifugation, pellet the cell suspension at 200 g for 5 minutes and carefully decant the supernatant without disturbing the cell pellet.
- 4. Resuspend the cells in 5 mL of EX-CELL CD CHO Fusion medium.
- 5. Count the cells for viability and transfer to a sterile tissue culture flask at a seeding density of 4 x 10<sup>5</sup> cells/mL.
- 6. When cell densities reach 1-2 x 10° cells/mL, passage the cells using standard cell culture techniques.

#### Characteristics

- Appearance Clear solution
- Osmolality (as supplied) 292 322 mOsm/kg H<sub>2</sub>O
- pH (as supplied) -7.1 7.7
- Sterility No microbial growth detected

#### **Product Profile**

SAFC EX-CELL CD CHO Fusion was compared to five chemically defined competitor formulations designed for CHO cells (Table 1). Three proprietary rlgG-producing CHO cell lines were used for the purposes of these comparisons.

NOTE: All cell lines were adapted to each formulation by passing six times before evaluating growth and productivity. Cell line 2 could not be adapted to GIBCO CD CHO.

Table 1: CHO Formulations Used for Comparison Assays

Medium	Manufacturer	Cat No.
EX-CELL CD CHO Fusion	SAFC	14365C/24365C
CD CHO Medium	GIBCO (Life Technologies)	10743
CD OptiCHO™ Medium	GIBCO (Life Technologies)	12681
PowerCHO®-2	BioWhittaker (Lonza)	12-771Q
CDM4CHO <sup>TM</sup>	HyClone (Thermo Scientific)	SH30558
IS CHO-CD4™	Irvine Scientific	9110

# Warranty, Limitation of Remedies

WARTAINTY, Limitation of nemeries

SAFC warrants to the purchaser for a period of one year from date of delivery that this product conforms to its specifications. Other terms and conditions of this warranty are contained in SAFC written warranty, a copy of which is available upon request. ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCUDING THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXCLUDING TO a case will SAFC be liable for any special, incidental, or consequential damages arising out of this product or the use of this product by the customer or any third party based upon breach of warranty, breach of contract, negligence, strict tort, or any other legal theory. SAFC expressly disclaims any warranty against claims by any third party by way of infringement or the like. THIS PRODUCT IS INTENDED FOR PURPOSES DESCRIBED ONLY third party by way of infringement or the like. THIS PRODUCT IS INTE AND IS NOT INTENDED FOR ANY HUMAN OR THERAPEUTIC USE

Additional Terms and Conditions are contained in the product Catalog, a copy of which is available upon request.

©2012 Sigma-Aldrich Co., LLC. All rights reserved. SAFC and SIGMA-ALDRICH are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. EX-CELL is a registered trademarks of Sigma-Aldrich Co. LLC. PowerCHO-2 is a registered trademark of Lonza Group Ltd. CD OptiCHO Medium is a trademark of Life Technologies Corporation. DMACHO is a trademark of Thermo Fischer Scientific, Inc. IS CHO-CO4 is a trademark of Irefundose Scientific Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the SAFC website at www.safcglobal.com and/or on the reverse side of the invoice or packing slip.

Issued August 2012 P14365 1108 1208

79469-509909



#### SAFC

3050 Spruce Street Saint Louis, MO 63103 USA Phone: 1-800-244-1173 • 1-314-534-4900

Fax: 1-888-289-8449 Global Email: safcglobal@sial.com





