

Product Information

Lysozyme from chicken egg white for Molecular Biology

Catalog Number **L7651**
Storage Temperature $-20\text{ }^{\circ}\text{C}$

CAS RN 12650-88-3
EC 3.2.1.17
Synonyms: Muramidase; Lysozyme c;
Mucopolysaccharide N-acetylmuramoylhydrolase

Product Description

Lysozyme is a single chain polypeptide of 129 amino acids cross-linked with four disulfide bridges.¹ It hydrolyzes $\beta(1\rightarrow4)$ linkages between N-acetylmuramic acid and N-acetyl-D-glucosamine residues in peptidoglycan and between N-acetyl-D-glucosamine residues in chitodextrin.^{2,3} The enzyme is often used for lysing bacterial cells by hydrolyzing the peptidoglycan present in the cell walls. Gram-positive cells are quite susceptible to this hydrolysis as their cell walls have a high proportion of peptidoglycan. Gram-negative bacteria are less susceptible due to the presence of an outer membrane and a lower proportion of peptidoglycan. However, these cells may be hydrolyzed more easily in the presence of EDTA that chelates metal ions in the outer bacterial membrane.^{4,5}

This lysozyme preparation is purified from chicken egg white, crystallized three times, dialyzed, and supplied as a lyophilized powder. Protein content by UV absorbance is $\geq 90\%$ with the remainder ($\sim 10\%$) being buffer salts such as sodium acetate and sodium chloride.

This highly purified enzyme preparation has been used in mass spectrometry as a protein mass calibration standard and in structural studies of proteins.⁶⁻⁸ It is suitable for use as a lysing agent in the purification of plasmid DNA using a boiling lysing technique.⁹

Molecular mass:¹⁰ 14,307 Da (amino acid sequence)

Isoelectric point (pI):¹¹ 11.35

Extinction co-efficients:
 $E^{1\%}$ (281.5 nm):¹² 26.4 in 0.1 M potassium chloride
 E^{mM} (280 nm):¹³ 36

Optimal pH:

The activity of lysozyme is a function of both pH and ionic strength. The enzyme is active over a broad pH range (6.0–9.0). At pH 6.2, maximal activity is observed over a wider range of ionic strengths (0.02–0.100 M) than at pH 9.2 (0.01–0.06 M).¹³

Inhibitors:

Lysozyme is inhibited by indole derivatives, which bind to and distort the active site, and imidazole, which induces the formation of a charge-transfer complex.¹⁴ It is also inhibited by surface-active agents such as sodium dodecyl sulfate, sodium dodecanate, and dodecyl alcohol. Other compounds of these types with carbon chains of 12 or more carbons in length will also inhibit lysozyme.¹⁵

Substrates:

The natural substrate for lysozyme is the peptidoglycan layer of bacterial cell walls. However, a variety of low molecular mass substrates including murein degradation products as well as synthetic compounds have been used for various photometric, isotopic, and immunological lysozyme assays.¹⁶

The following low molecular mass lysozyme substrates are available:

- 4-Methylumbelliferyl β -D-N,N',N''-triacetyl-chitotrioside (Catalog Number M5639, a fluorogenic substrate)
- 4-Nitrophenyl β -D-N,N',N''-triacetylchitotriose (Catalog Number N8638)

Lysozyme activity: $\geq 40,000$ units/mg protein

Unit definition: One unit will produce a change in A_{450} of 0.001 per minute at pH 6.24 at $25\text{ }^{\circ}\text{C}$, using a suspension of *Micrococcus lysodeikticus* as substrate, in a 2.6 ml reaction mixture (1 cm light path).

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

For *E. coli* cell lysis, use a freshly prepared lysozyme solution (10 mg/ml) in 10 mM Tris-HCl, pH 8.0.⁹

The product is also soluble in water (10 mg/ml) yielding a clear to slightly hazy colorless solution. Aqueous solutions should retain activity for at least one month when stored between 2–8 °C.

Storage/Stability

The product, as supplied, should be stored at –20 °C. When stored at –20 °C, the enzyme retains activity for at least 4 years.

Solutions (pH 4–5) remain active for several weeks if refrigerated.

Procedure

The following procedure is for the lysis of *E. coli*. It may be used as a guideline for other species. The optimal pH for *E. coli* cell lysis is 8.0±0.1.⁹

1. Incubate *E. coli* (strain ATCC 37017) bearing the pBR322 plasmid overnight in Terrific Broth (Catalog Number T0918) with 25 µg/ml tetracycline (Catalog Number T3383) and 25 µg/ml ampicillin (Catalog Number A9518).
2. Centrifuge 1–2 ml samples of the overnight culture.
3. Resuspend the pellets in 350 µl of STET buffer (10 mM Tris-HCl, pH 8.0, with 0.1 M NaCl, 1 mM EDTA, and 5% [w/v] TRITON® X-100).
4. Add 25 µl of a freshly prepared lysozyme solution (10 mg/ml in 10 mM Tris-HCl, pH 8.0).
5. Mix by vortexing for 3 seconds.
6. Incubate the lysis mixture for 30 minutes at 37 °C
7. After incubation, place the tube containing the lysis mixture in a boiling water bath for exactly 40 seconds.
8. Centrifuge the lysis mixture at 14,000 × *g*.
9. Remove the pellet (cell debris) from the tube using a sterile toothpick.
10. Plasmid DNA from the supernatant may then be purified and analyzed.

References

1. Jolles, P., *Angewandte Chemie, International Edition*, **8**, 227-239 (1969).
2. Rupley, J.A., *Biochim. Biophys. Acta*, **83**, 245-255 (1964).
3. Holler, H., et al., *Biochem.*, **14**, 2377-2385 (1975).
4. Schutte, H., et al., *Biotech. Applied Biochem.*, **12**, 599-620 (1990).
5. Vazquez, -Laslop, N., et al., *J. Bact.*, **183**, 2399-2404 (2001).
6. Galvani, M., et al., *Electrophoresis*, **22**, 2058-2065 (2001).
7. Abgar, S., et al., *Eur. J. Biochem.*, **267**, 5916-5925 (2000).
8. Sethuraman, A., et al., *Proteins: Structure, Function and Bioinformatics*, **56**, 669-678 (2004).
9. Sambrook, J., et al., in *Molecular Cloning, a Laboratory Manual*, Cold Spring Harbor Laboratory Press, (Cold Spring Harbor, NY: 1989) p 1.29.
10. Canfield, R.E., *J. Biol. Chem.*, **238**, 2698-2707 (1963).
11. Wetter, L.R., et al., *J. Biol. Chem.*, **192**, 237-242 (1951).
12. Aune, K.C., et al., *Biochem.*, **8**, 4579-4585 (1965).
13. Davies, R.C., et al., *Biochim. Biophys. Acta*, **178**, 294-305 (1969).
14. Swan, I., *J. Mol. Bio.*, **65**, 59-62 (1972).
15. Smith, G., and Stoker, C., *Arch. Biochem. Biophys.*, **21**, 383-394 (1949).
16. Holtje, J.V., *EXS*, **75**, 105-110 (1996).

TRITON is a registered trademark of the Union Carbide Corp.

CS, RBG, ALC, MAM 03/08-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.