



EX-CELL® 302 Serum-Free Medium for CHO Cells

without L-glutamine

Catalog No. 14324C

Product Description

EX-CELL® 302 is a serum-free liquid medium, which has been specifically developed for the long-term growth of transformed Chinese Hamster Ovary (CHO) cells in suspension for the expression of antibodies or protein products. EX-CELL® 302 is an appropriate medium for use with the DHFR- or Glutamine Synthetase, or the GS System™, because it does not contain hypoxanthine, thymidine or L-glutamine. If EX-CELL® 302 is to be used with CHO cells transformed using the GS selection system the media needs to be supplemented with GS Supplement 50X (Catalog No. 58672C) for the additional amino acids and nucleosides. For use with DHFR- selection system we recommend the addition of 4 mM L-glutamine for optimal growth.

Due to current regulatory concerns about the sources of raw materials, EX-CELL® 302 was developed using only recombinant human proteins that have molecular weights less than 10 kD. The total protein concentration found in EX-CELL® 302 is less than 1 mg/L. Pluronic® F68 has been added at a final concentration of 0.1% to EX-CELL® 302 to protect against shear damage in sparged bioreactor systems.

Catalog No. 14324C replaces Catalog No. 14312 and includes an alternate source of soy hydrolysate to that found in the original EX-CELL® 302 formulation. With more consistent performance and improved filtration characteristics, the alternate hydrolysate will improve the overall performance and consistency of EX-CELL® 302. Comparability testing utilizing the previous soy hydrolysate and the replacement hydrolysate demonstrated comparable growth-promoting characteristics.

Formulation

The formulation for EX-CELL® 302 is proprietary to SAFC Biosciences. For additional information please call our Technical Services department.

Precautions

Use aseptic technique when handling or supplementing this medium. This product is for research or for further manufacturing use. THIS PRODUCT IS NOT INTENDED FOR HUMAN OR THERAPEUTIC USE.

Storage

Store liquid medium at 2° to 8° C, protected from light. Do not use after the expiration date.

Indications of Deterioration

Medium should be clear and free of particulates and flocculent material. Do not use if medium is cloudy or contains precipitate. Other evidence of deterioration may include color change, pH shift or degradation of physical or performance characteristics.

Preparation Instructions

EX-CELL® 302 is formulated with sodium bicarbonate and without L-glutamine. For applications requiring the use of L-glutamine, supplement with 4 mM L-glutamine by adding 20 mL/L of a 200 mM solution (Catalog No. 59202C) prior to use. We recommend L-glutamine supplementation of the working volume only. Supplements, such as antibiotics, can be added to the medium using aseptic technique. Storage conditions and shelf life of the product may be affected by the nature of the supplements.

Methods for Use

Adaptation

Mammalian cells can be adapted to serum-free conditions by direct adaptation from serum-containing media or by gradual weaning. Both procedures require healthy, viable cultures in mid-logarithmic growth phase. During the adaptation phase growth rates will usually be somewhat slower than growth in media supplemented with serum.

Gradual Weaning to Serum-Free Media (Recommended):

1. Passage cells from the serum-supplemented medium directly into EX-CELL® 302 + 5% gamma irradiated Fetal Bovine Serum (FBS) (Catalog No. 12107C) using normal seeding densities.
2. Once confluent (80 - 90%), rinse with 20 mL of Dulbecco's Phosphate Buffered Saline (DPBS Modified) (Catalog No. 59321C) and trypsinize with 5 mL of trypsin. Add 15 mL of EX-CELL® 302 + 2% FBS. Obtain a cell count and pellet the cells by centrifugation.
3. Passage cells into a 75 cm² flask containing 30 mL EX-CELL® 302 + 2% FBS at a density of 1×10^5 cells/mL.
4. Allow cells to reach 80 - 90% confluency. Retain any detached cells. Observe viability of detached cells. If detached cells are not viable, discard. If viable, reserve and combine with trypsinized cells.
5. Trypsinize cells as in Step 2, reducing the amount of trypsin to 2 mL. Inactivate trypsin with EX-CELL® 302 + 1% FBS. Centrifuge.
6. Resuspend cells in EX-CELL® 302 + 1% FBS at a density of 2×10^5 cells/mL.
7. Observe cells daily. Gently rap the flask against the palm of the hand to detach cells, return to the incubator until cells reach 80 - 90% confluency (~ 4 - 5 days).
8. Once confluent, retain detached cells and rap flask 2 - 3 times on the palm of the hand or a padded surface. Avoid causing bubbles. Triturate cells to further remove cells from the flask and to break up clumping cells. Obtain a cell count and pellet the cells by centrifugation.
9. Resuspend cells at a density of 8×10^5 cells/mL in 75 cm² flasks containing EX-CELL® 302 without serum until cells reach a density of $1.5\text{--}2 \times 10^6$ cells/mL (~ 3 - 4 days).
10. Triturate cells to break up clumps. Obtain a cell count and pellet the cells by centrifugation.
11. Resuspend in EX-CELL® 302 at a density of 8×10^5 cells/mL. Incubate until cells reach a density of $1.5\text{--}2 \times 10^6$ cells/mL. Take note of cell clump size (number of cells/clump).
12. Passage cells as outlined in Steps 9 - 11.
13. Repeat Step 12 except seed 250 mL shaker flasks at 6×10^5 cells/mL in 60 mL at 60 rpm.
14. After the cultures have been established in EX-CELL® 302, passage the cells for an additional 2 passages in 250 mL shaker flasks.
15. Cells can now be scaled up as necessary using standard procedures and densities.

NOTES: Adjust protocol as necessary to ensure proper adaptation to shakers. Start with shaker speed of 60 rpm and adjust over the next few passes.

If cell clumping is extensive, select for single cells by letting the large clumps settle to the bottom of the flask and pipetting the single cells for passage. Be sure to count total cells and then count the retained single cell suspension to ensure proper seed density.

Always centrifuge cells when passing during the adaptation.

Culture density should reach $2\text{--}3 \times 10^6$ cells/mL over 3 - 4 days.

Direct Adaptation to Serum-Free Media:

1. Passage the cells into pre-warmed (37 C) EX-CELL® 302 at 1.5-2X the recommended seeding density ($2\text{--}5 \times 10^5$ cells/mL).
2. Refeed the culture after 48 hours with a 100% exchange of fresh EX-CELL® 302.
3. Allow the cultures to achieve a minimum density of 1×10^6 cells/mL before subculturing.
4. Subculture into fresh EX-CELL® 302 as in Step 1 using normal seed densities.

Culture Techniques

Once cultures are fully adapted, the cells should be passed every 3 - 4 days at a seeding density of at least 2×10^5 cells/mL. An optimal seeding density should be determined by the researcher for each application and cell type.

When passing the cells, carryover should not exceed 25% of the final volume. If carryover exceeds 25%, centrifugation is recommended. Cells propagated in serum-free medium are extremely fragile. Standard techniques for centrifugation must be modified to include low-speed centrifugation to prevent damage to cells that have been propagated in serumfree medium.

During adaptation, normal trypsin concentrations may be used, but incubations should be carried out at 4° C, and exposure time should be minimal. We recommend the use of a soybean trypsin inhibitor (0.1%), or sedimentation by centrifugation to remove the trypsin. Soybean trypsin inhibitor should be used with caution, as it is toxic to some cell types. Cells may also be dislodged with NO-ZYME™ (Catalog No. 59226C), a non-enzymatic dissociating agent.

Cryopreservation

Freezing:

Cells can be frozen in EX-CELL® 302 without the reintroduction of serum.

1. Choose cultures in logarithmic growth with viabilities above 90%.
2. Prepare a freezing medium consisting of 45% cold EX-CELL® 302 medium, 45% spent medium and 10% dimethyl sulfoxide (DMSO).
3. Centrifuge the cells at 200 g for 5 minutes. Remove the supernatant.
4. Resuspend the cells in the freezing medium at 5×10^6 to 1×10^7 cells/mL.
5. Rapidly transfer 1 - 2 mL of this suspension to sterile cryovials.
6. Place the vials at -20°C for 3 - 4 hours, then transfer to -70°C for 16 - 24 hours.
7. For long-term storage, transfer the vials to liquid nitrogen vapor.

Thawing:

1. Rapidly thaw a vial of frozen cells in a 37°C water bath.
2. Transfer the cells aseptically to a centrifuge tube containing 10 mL of cold EX-CELL® 302 medium.
3. Using low-speed centrifugation, pellet the cell suspension at 200 g for 5 minutes and carefully decant the supernatant without disturbing the cell pellet.
4. Resuspend the cells in 5 mL of EX-CELL® 302 medium.
5. Count the cells for viability and transfer to a sterile shaker flask at a seeding density of 6×10^5 cells/mL.
6. When the culture has reached a density of 1×10^6 cells/mL, passage the cells using standard cell culture techniques.

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

To place an order or receive technical assistance

In Europe, please call Customer Service:

France: 0825 045 645

Germany: 069 86798021

Italy: 848 845 645

Spain: 901 516 645 Option 1

Switzerland: 0848 645 645

United Kingdom: 0870 900 4645

For other countries across Europe,
please call: +44 (0) 115 943 0840

Or visit: **MerckMillipore.com/offices**

For Technical Service visit:

MerckMillipore.com/techservice

MerckMillipore.com

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

