



## Product Information

### Cyclosporin A from *Tolypocladium inflatum*

Product Number **C 3662**  
Storage Temperature 2-8 °C

CAS RN: 59865-13-3

Synonyms: Cyclosporine; Antibiotic S 7481F1;  
Ciclosporin; CsA

Molecular formula: C<sub>62</sub>H<sub>111</sub>N<sub>11</sub>O<sub>12</sub>  
Molecular weight: 1202.61

Melting point: 148-151°C<sup>1</sup>  
[α]<sub>D</sub><sup>20</sup>: -244° (c = 0.6 in chloroform)<sup>1</sup>

#### Product Description

Cyclosporin A (CsA) is a non-polar cyclic oligopeptide metabolite from the fungus *Tolypocladium inflatum*. It has too narrow a spectrum of antifungal activity to be very useful as an antibiotic, but it possesses potent immunosuppressive properties, affecting primarily T-lymphocytes<sup>2</sup>. It has been shown to inhibit the functioning of several nuclear proteins involved in T-cell activation at the level of mRNA transcription.<sup>3</sup>

Cyclosporin is the primary tool used to prevent rejection following solid organ and bone marrow transplantation. It forms a complex with its intracellular receptor cyclophilin, which can then bind to calcineurin, inhibiting its enzymatic activity.<sup>4,5,6</sup> CsA was found to suppress the replication of hepatitis C virus genome in cultured hepatocytes<sup>7</sup>

In a study of its specific disruption of renal function (noting its hepatotoxicity) and of gene transcription, CsA was administered to rats intramuscularly at a dose of 7.5 mg/kg using a mixture of 100 mg CsA per mL in 90% olive oil, 10% ethanol.<sup>8</sup> In another study, CsA at concentrations >10 nM protected isolated hepatocytes against the action of phalloidin.<sup>9</sup> Measuring the concentration of CsA in solution by HPLC was shown to be significantly temperature-dependent, due to interconversion of CsA between two forms.<sup>10</sup> An extensive list of references has been reported, including a comprehensive review of analytical properties.<sup>11</sup>

Sigma also offers C 1832, Cyclosporin A, which has received additional testing for molecular biology applications. In Jurkat cells (a leukemic T-cell line), the production of interleukin-2 was inhibited by 90% in the presence of 1 µg/ml C 1832.<sup>12</sup>

#### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### Preparation Instructions

Sigma assays involve only organic solvents: 10 mg/mL in methylene chloride, 6 mg/mL in chloroform, 10 mg/mL in ethanol, 50 mg/mL in DMSO. Solutions are clear, colorless to faint yellow.<sup>12</sup> Cyclosporin is reported to be "slightly soluble in water and saturated hydrocarbons."<sup>1</sup>

Stock solutions in ethanol or DMSO should be stored at -20 °C. Cyclosporin is stable in solution if protected from light, but its concentration may drop due to adsorption to the container walls.

[The concentration of] "Cyclosporin was stable over 72 hours following dilution in glucose 5% or glucose/amino-acid solutions and storage at room temperature in the dark; similar stability was seen following dilution in lipid emulsion, but dilutions in sodium chloride 0.9% were considered to be stable only for 8 hours." In this study, the solution was "stable if the initial cyclosporine concentration remained at 90% or above."<sup>13,14</sup>

#### Storage/Stability

Store the product desiccated and protected from light at 2-8 °C. Under these conditions the product is stable for 2 years. It should be re-evaluated for suitability in user's application every two years.

## References

1. Merck Index, 12th ed., No. 2821 (1996).
2. Dreyfuss, M., et al., Eur. J. Appl. Microbiol. Biotechnol., **3**, 125-33 (1976).
3. Emmel, E.A., et al., Science, **246**, 1617-20 (1989)
4. Clipstone, N.A., et al., J. Biol. Chem., **269**, 26431-437 (1994).
5. Nichols, R.A., et al., J. Biol. Chem., **269**, 23817-823 (1994).
6. Nelson, P.A., et al., J. Immunol., **150**, 2139-47 (1993).
7. Watashi, K, et al., Hepatology, **38**, 1282-8 (2003)
8. Morris, S.M., et al., J. Biol. Chem., **267**, 13768-771 (1992).
9. Ziegler, K. and Frimmer, M., Biochim. Biophys. Acta, **805**, 174-180 (1984).
10. Nishikawa, T., et al., Chromatographia, **38**, 359-363 (1994).
11. Hassan, M.M. and Al-Yahya, M.A., Analytical Profiles of Drug Substances, **16**, 145-206 (1987).
12. Sigma quality control data
13. Martindale: The Extra Pharmacopoeia, 30th Ed., (Pharmaceutical Press, 1993) p. 467.
14. McLeod, H.L., et al., Ann. Pharmacother., **26**, 172-175 (1992).

NDH/PHC 12/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.