

Magnesium Assay Kit

Catalogue number MAK507

Product Description

Magnesium (Mg) is one of the most abundant and essential minerals in mammals. Magnesium is involved in more than 300 biochemical reactions in the body and plays important roles in muscle and nerve functions, heart rhythm, immune system and bone formation. Simple, direct and automation-ready procedures for measuring magnesium concentration in biological samples are becoming popular in metabolomics research and drug discovery. The Magnesium Assay Kit is designed to measure magnesium directly in biological samples without any pretreatment. A calmagite dye in the kit forms a colored complex with magnesium. The intensity of the color, measured at 500 nm, is directly proportional to the magnesium concentration in the sample. The optimized formulation minimizes interference by potential substances.

The linear detection range of the kit is 0.1 mg/dL (41 µM) to 3 mg/dL (1.2 mM) Mg²⁺. The kit is suitable for determination of magnesium in biological samples such as plasma, serum, and urine. The kit is also suitable for determination of magnesium in deproteinated samples such as milk and environmental samples such as water and soil.

Components

The kit is sufficient for 250 colorimetric assays in 96-well plates.

- | | |
|---|--------|
| • Reagent A | 25 mL |
| Catalogue Number MAK507A | |
| • Reagent B | 25 mL |
| Catalogue Number MAK507B | |
| • EDTA Solution | 1.5 mL |
| Catalogue Number MAK507C | |
| • Standard (10 mg/dL Mg ²⁺) | 1 mL |
| Catalogue Number MAK507D | |

Reagents and Equipment Required but Not Provided

- Pipetting devices and accessories (example, multichannel pipettor)
- Spectrophotometric multiwell plate reader.
- Clear bottom 96-well plates (example., Corning Costar). Cell culture or tissue culture treated plates are not recommended.
- 1.5 mL microcentrifuge tubes.
- Trichloroacetic Acid Solution (Catalogue number T0699)

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

The kit is shipped at room temperature. Store all components at 4 °C.

Preparation Instructions

Briefly centrifuge small vials prior to opening. Equilibrate all components to room temperature prior to use.

Procedure

All Samples and Standards should be run in duplicate. This procedure is written for 96-well plates.

Note: The procedure may be used with cuvettes by scaling up the volumes 5x.

Sample Preparation

Serum and urine samples can be measured directly.

Note: EDTA and other Mg^{2+} chelators interfere with this assay. This assay cannot be applied to plasma samples obtained with EDTA.

Milk and Lipid/Protein Rich Samples:

1. Mix equal volumes of sample and 10% TCA.
2. Incubate at room temperature for 5 minutes.
3. Centrifuge for 2 minutes at 14,000 rpm.
4. Proceed using the supernatant for the assay.

Add 5 μ L of each sample to separate wells in a clear, flat-bottom 96-well plate.

Standard Preparation

1. Dilute Standard to 2 mg/dL by mixing 40 μ L of 10 mg/dL Standard with 160 μ L purified water.
2. Transfer 5 μ L of diluted standard into the 96-well plate.

Working Reagent

1. Prepare enough working reagent for each sample and standard.
2. Combine equal volumes of Reagent A and Reagent B.

Assay Reaction

1. Add 200 μ L of Working Reagent to each Sample and Standard well. Tap plate to mix briefly and thoroughly.
2. Incubate for 2 minutes at room temperature.
3. Read optical density at 500 nm. This is the OD value for the samples and standard.
4. Add 10 μ L EDTA solution to all wells and tap plate to mix thoroughly.
5. Incubate 2 minutes at room temperature.
6. Read optical density at 500 nm. This is the OD value for the blanks.

Results

Calculate the concentration of Mg^{2+} in the sample using the following equation:

$$Mg^{2+} \text{ Conc.} = \frac{OD_{\text{Sample}} - OD_{\text{Blank}}}{OD_{\text{Mg}} - OD_{\text{MgBlank}}} \times 2 \text{ (mg/dL)}$$

Where:

OD_{Sample} and OD_{Blank} are the optical density readings before and after the addition of EDTA.

OD_{Mg} and OD_{MgBlank} are the optical density readings of the standard (2 mg/dL) before and after the addition of EDTA.

Conversions: 1 mg/dL Mg^{2+} equals 411 μ M, 0.001% or 10 ppm.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.
All other trademarks are the property of their respective owners. Detailed information on
trademarks is available via publicly accessible resources.

© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
mak507pis Rev 04/24

