

# GenElute™-E Single Spin Checklist for Tissue DNA 96 Kit

## EC396 Preparation Before Starting

- Heat the thermal shaker or heating block/chamber to 60 °C.
- Set the centrifuge to 1,000 x g.

## Lysis

- Add 1–20 mg of tissue sample per well of Lysis Plate.
- Prepare Lysis Master Mix, add 135 µL per well of Lysis Plate.

Number of samples	1	96 (+20%)
Tissue Lysis Buffer <b>LB</b>	130 µL	14,976 µL
SmartLyse™ T Protease <b>P</b>	5 µL	576 µL
Final Volume	135 µL	15,552 µL

## Plate Preparation during lysis

Seal Lysis Plate tightly with Adhesive Foil. Incubate



60 °C



80 °C



Thermal Shaker maximum agitation

- Remove Adhesive Foil. Prepare RNase Digest Master Mix, add 16 µL per well of Lysis Plate.

Number of samples	1	96 (+20%)
Clearing Solution T <b>CS</b>	15 µL	1,728 µL
RNase A Tissue <b>R</b>	1 µL	115.2 µL
Final Volume	16 µL	1,843.2 µL

- Mix by pipetting. Centrifuge Lysis Plate for 3 minutes at maximum speed.

## Preparation of Purification Plate (during 60 °C and 80 °C incubation)

- Detach lower and upper sealing foils from purification plate.
- Place Purification Plate on top of a Conditioning Plate.
- Centrifuge 1 minute at 1,000 x g to collect void buffer.
- Place conditioned Purification Plate on top of Storage Plate.

## Purification of DNA

- Transfer lysis supernatant from Lysis Plate to Purification Plate.
- Centrifuge 1 minute at 1,000 x g to collect DNA into the Storage Plate.
- Collected DNA is ready to use.