

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

2,4-Dichlorophenoxyacetic acid Plant Cell Culture Tested

Product Number **D 7299**Store at Room Temperature

Product Description

Molecular Formula: C₈H₆Cl₂O₃ Molecular Weight: 221.0 CAS Number: 94-75-7 Melting point: 140.5 °C¹ Synonym: 2,4-D

This product is plant cell culture tested at a concentration of 2 mg/L.

- 2,4-Dichlorophenoxyacetic acid is a synthetic auxin. It is the auxin most commonly used for inducing callus in culture; it is also used for maintaining cells in the dedifferentiated state. A typical working concentration is in the 1–50 μ M range.²
- 2,4-D is the metabolite of 2,4-Dichlorophenoxybutyric acid in higher plants.^{3,4}
- 2,4-D has been used at a concentration of 3 mg/L in the culture of *Lilium speciosum* Thunb. Var. *gloriosoides* Baker, a perennial that is a lily variety.⁵

The IR spectrum (Nujol mull) has been published.⁶

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in 95% ethanol (100 mg/ml) yielding a clear to slightly hazy, faint yellow solution. Heat may be required to get the product to dissolve at this concentration.

Storage/Stability

Dilute aqueous solutions of this compound may be sterilized by autoclaving.

References

- 1. The Pesticide Manual, 9th ed., Worthing, C.R. and Hance, R.J. (The British Crop Protection Council, Surrey, UK, 1991) pp. 218-219.
- Plant Cell Culture: A Practical Approach, 2nd ed., Dixon R.A and Gonzales, R.A, Eds. (IRL Press, Oxford, 1994) p. 10.
- Hayashi, M. et al., 2,4-Dichlorodphenoxybutyric Acid–Resistant Mutants of Arabidopsis Have Defects in Glyoxysomal Fatty Acid β–Oxidation, The Plant Cell, 10, 183-195 (1998).
- 4. Wain, R.L., and Wightman, F., The growth-regulating activity of certain ω –substituted alkyl carboxylic acids in relation to their β –oxidation within the plant. Proc. R. Soc. Lond. Biol. Sci. 142, 525-536 (1954).
- 5. Chang, C. et al., A tissue culture protocol for propagation of a rare plant, *Lilium speciosum* Thunb. var. *gloriosoides* Baker, Bot. Bull. Acad. Sin. 41, 139-142 (2000).
- 6. Sigma Library of FT-IR Spectra, **Ed. 1**, **Vol. 2**, p. 450A.

IRB/NSB 4/03