Instructions for Use

Elastic Stain

Procedure No. HT25

Intended Use

Sigma-Aldrich Elastic Stain kit is intended for histologic staining of elastin in tissue sections. Elastic Stain reagents are for "In Vitro Diagnostic Use". For professional use only. The data obtained from this manual, qualitative procedure serves only as an aid to diagnosis of human specimens and should be reviewed in conjunction with other clinical diagnostic tests or information.

Elastin, a type of connective tissue, is protein occurring either as branching fibers or as sheets.¹ It is highly refractive, elastic and usually thinner than collagen fibers.

Sections are overstained in hematoxylin-iodine-ferric chloride solution in a dye-lake reaction. Ferric chloride and iodine function primarily as mordants, additionally acting as oxidizers assisting the conversion of hematoxylin dye to hematein.² The cationic dye-lake is trapped by iodine allowing the staining of elastic fibers by dye-lake and retarding extraction of stain from elastic fibers during the differentiating step.³ Since sections are overstained, differentiation is necessary and is accomplished by use of excess mordant (a dilute ferric chloride solution) to break tissue-mordant-dye complex. Van Gieson solution is used as a counterstain to stain collagen fibers red and other components yellow.

Demonstration of elastic tissue is useful in cases of emphysema (i.e., atrophy of elastic tissue), arteriosclerosis (i.e., thinning and loss of elastic fibers) and various other vascular diseases.

Reagents

Hematoxylin Solution, Alcoholic (Cat. No. HT251-200ML)

Hematoxylin, 5% (w/v), C.I. 75290 in SD 3A alcohol. Danger. Highly flammable liquid and vapour. Harmful if swallowed. Causes mild skin irritation. Causes damage to organs.

Ferric Chloride Solution (Cat. No. HT252-125ML)

Aqueous solution of ferric chloride, 29% (w/v). Danger. May be corrosive to metals. May be harmful if swallowed or in contact with skin. Causes eye irritation. Causes serious eye damage. Toxic to aquatic life with long lasting effects. Avoid release to the environment.

Weigert's Iodine Solution (Cat. No. HT253-180ML)

Aqueous solution of iodine, 1% (w/v), and potassium iodide, 2% (w/v). Warning. Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. Harmful to aquatic life. Wear protective gloves. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.

Van Gieson Solution (Cat. No. HT254-250ML)

Acid fuchsin, 0.05% (w/v), C.I. 42685 in 95% saturated picric acid solution.

Special Materials Required but Not Provided

- Positive control slides, such as Sigma Elastic TISSUE-TROL™ (Cat. No. TTR003) should be included in each run
- Alcohol, 95% and 100%
- Xylene or xylene substitute

Storage and Stability

Store Elastic Stain reagents at room temperature (18–26°C). Reagent labels bear expiration dates.

Working Elastic Stain Solution and Working Ferric Chloride Solution should be prepared fresh daily.

Preparation

Working Elastic Stain Solution is prepared by adding the following reagents, in the order indicated, into a Coplin jar. Mix well.

Hematoxylin Solution, Alcoholic	20 mL
Ferric Chloride Solution	3 mL
Weigert's Iodine Solution	8 mL
Deionized Water	5 mL

Working Ferric Chloride Solution (for differentiation) is prepared by adding 3 mL of Ferric Chloride Solution and 37 mL deionized water to a Coplin jar. Mix well.

Van Gieson Solution is provided ready for use.

Precautions

The IVDs included in this kit are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Elastic TISSUE-TROL™ control slides are paraffin embedded human tissue containing elastic fibers and should be considered potentially infectious.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Any well fixed tissue (10% neutral buffered formalin or Zenker's solution is preferred).

Notes

It is not necessary to remove mercury deposits if tissues have been fixed in a mercury containing fixative, since it will be removed by the staining solution. Embed in paraffin and cut at 5 microns.

Procedure

1. Deparaffinize slides to deionized water.
2. Place in Working Elastic Stain Solution for 10 minutes.
3. Rinse in deionized water.
4. Differentiate in Working Ferric Chloride Solution.
5. Rinse in tap water.
6. Check microscopically. If over differentiated, return to Working Elastic Stain Solution.
7. Rinse in 95% alcohol to remove iodine.⁴
8. Rinse in deionized water.
9. Stain in Van Gieson Solution 1–3 minutes.
10. Rinse in 95% alcohol.
11. Dehydrate to Xylene or xylene substitute and mount.

Performance Characteristics

Target Structure	Staining Result
Elastic Fibers	Black to Black
Nuclei	Blue to Black
Collagen	Red
Muscle and Other	Yellow

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT251	Hematoxylin Solution, Alcoholic	Elastic	3 of 3	3 of 3	3 of 3	3 of 3
		Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Collagen	3 of 3	3 of 3	3 of 3	3 of 3
		Muscles	3 of 3	3 of 3	3 of 3	3 of 3
HT252	Ferric Chloride Solution	Elastic	3 of 3	3 of 3	3 of 3	3 of 3
		Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Collagen	3 of 3	3 of 3	3 of 3	3 of 3
		Muscles	3 of 3	3 of 3	3 of 3	3 of 3
HT253	Weigert's Iodine Solution	Elastic	3 of 3	3 of 3	3 of 3	3 of 3
		Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Collagen	3 of 3	3 of 3	3 of 3	3 of 3
		Muscles	3 of 3	3 of 3	3 of 3	3 of 3
HT254	Van Gieson Solution	Elastic	3 of 3	3 of 3	3 of 3	3 of 3
		Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Collagen	3 of 3	3 of 3	3 of 3	3 of 3
		Muscles	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT25A:

H225: Highly flammable liquid and vapor.

H290: May be corrosive to metals.

H315: Causes skin irritation.

H318: Causes serious eye damage.

H371: May cause damage to organs (Eyes).

H373: May cause damage to organs (Thyroid) through prolonged or repeated exposure.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P260: Do not breathe dust/ fume/ gas/ mist/ vapors/ spray.

P280: Wear protective gloves/ protective clothing/ eye protection/ face protection.

P305 + P351 + P338 + P310: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/ doctor.

P370 + P378: In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer

References

- 1. Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, p 175
- 2. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, (MO), 1980, pp 196-197
- 3. Selected Histochemical and Histopathological Methods, SW Thompson, Editor, CC Thomas, Springfield (IL), 1966, p 794
- 4. Histopathologic Technic and Practical Histochemistry. 3rd ed. RD Lillie, Editor, McGraw-Hill, New York, 1965, p 551

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 5.0	2016
Rev. 6.0	2022
Rev. 7.0	2022
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards.	

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271 Darmstadt,
Germany

The Initial M, TISSUE-TROL and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

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Referenzen

- Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, p 175
- Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, (MO), 1980, pp 196-197
- Selected Histochemical and Histopathological Methods, SW Thompson, Editor, CC Thomas, Springfield (IL), 1966, p 794
- Histopathologic Technic and Practical Histochemistry. 3rd ed. RD Lillie, Editor, McGraw-Hill, New York, 1965, p 551

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Riferimenti

- Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, p 175.
- Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, (MO), 1980, pp 196-197.
- Selected Histochemical and Histopathological Methods, SW Thompson, Editor, CC Thomas, Springfield (IL), 1966, p 794.
- Histopathologic Technic and Practical Histochemistry. 3rd ed. RD Lillie, Editor, McGraw-Hill, New York, 1965, p 551.

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Referenser

- 1. Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, p 175
- 2. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, (MO), 1980, pp 196–197
- 3. Selected Histochemical and Histopathological Methods, SW Thompson, Editor, CC Thomas, Springfield (IL), 1966, p 794
- 4. Histopathologic Technic and Practical Histochemistry. 3rd ed. RD Lillie, Editor, McGraw-Hill, New York, 1965, p 551

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Referências

- Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, p 175
- Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, (MO), 1980, pp 196–197
- Selected Histochemical and Histopathological Methods, SW Thompson, Editor, CC Thomas, Springfield (IL), 1966, p 794
- Histopathologic Technic and Practical Histochemistry. 3rd ed. RD Lillie, Editor, McGraw-Hill, New York, 1965, p 551

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Reference

1. Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, Londýn 1985, str. 175
2. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, (MO), 1980, pp 196–197
3. Selected Histochemical and Histopathological Methods, SW Thompson, Editor, CC Thomas, Springfield (IL), 1966, p 794
4. Histopathologic Technic and Practical Histochemistry. 3rd ed. RD Lillie, Editor, McGraw-Hill, New York, 1965, p 551

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Referanslar

- Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, s. 175
- Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, (MO), 1980, s. 196-197
- Selected Histochemical and Histopathological Methods, SW Thompson, Editor, CC Thomas, Springfield (IL), 1966, s. 794
- Histopathologic Technic and Practical Histochemistry. 3rd ed. RD Lillie, Editor, McGraw-Hill, New York, 1965, s. 551

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MDSS GmbH