For life science research only. Not for use in diagnostic procedures.



# **Nick Translation Mix**

Version: 10
Content Version: July 2020

For generation of highly sensitive probes for fluorescence in situ hybridization (FISH).

**Cat. No. 11 745 808 910** 200 μl 50 labeling reactions

Store the product at -15 to -25°C.

1.	General Information	3
1.1.	Contents	3
1.2.	Storage and Stability	
	Storage Conditions (Product)	3
1.3.	Additional Equipment and Reagent required	3
1.4.	Application	3
2.	How to Use this Product	4
2.1.	Before you Begin	4
	Sample Materials	
	Templates for labeling reaction	
	General ConsiderationsLabeling of in situ probes	
	Labeling of PCR products or varying amounts of template	
	Removal of unincorporated free fluorophore-dNTP	
	Nucleotide sets	
	Working Solution	
	Preparation of fluorophore labeling mix	
2.2.	Protocols	
	Standard labeling reaction  Determination of fragment length of labeled probe	
3.	Additional Information on this Product	
3.1.	Test Principle	
	In situ hybridization	
3.2.	Quality Control	7
4.	Supplementary Information	8
4.1.	Conventions	8
4.2.	Changes to previous version	8
4.3.	Ordering Information	8
4.4.	Trademarks	9
4.5.	License Disclaimer	9
4.6.	Regulatory Disclaimer	9
4.7.	Safety Data Sheet	9
4.8.	Contact and Support	9

## 1. General Information

#### 1.1. Contents

Vial / Bottle	Label	Function / Description	Content
1	Nick Translation Mix, 5x conc.	<ul> <li>Premixed solution for <i>in situ</i> probes.</li> <li>Stabilized reaction buffer in 50% glycerol (v/v) and DNA Polymerase I, DNase I.</li> </ul>	1 vial, 200 μl

# 1.2. Storage and Stability

## **Storage Conditions (Product)**

When stored at -15 to -25°C, the product is stable through the expiry date printed on the label.

Vial / Bottle	Label	Storage
1	Nick Translation Mix, 5x conc.	Store at −15 to −25°C.  Avoid repeated freezing and thawing. To avoid contamination, aliquot the Nick Translation Mix solution into 2 to 3 vials.

# 1.3. Additional Equipment and Reagent required

#### For preparation of fluorophore labeling mix

- Deoxynucleoside Triphosphate Set\*
- Sterile, double-distilled water
- Fluorescein-12-dUTP\* or Tetramethyl-Rhodamine-5-dUTP\*

#### For standard labeling reaction

- Heating block
- 0.5 M EDTA, pH 8.0
- 5x fluorophore labeling mix
  - See section, Working Solution.
- Sterile, double-distilled water

#### For determination of fragment length of labeled probe

- Agarose minigel
- Gel loading buffer
- Heating block
- DNA Molecular Weight Marker\*
- 0.5 M EDTA, pH 8.0

# 1.4. Application

The Nick Translation Mix is designed for direct fluorophore labeling of *in situ* probes.

- Fluorescein-12-dUTP\* and Tetramethyl-Rhodamine-5-dUTP\*, or other commercially available fluorophore-labeled nucleotides can be combined with the Nick Translation Mix for *in situ* probes.
- Direct flurophore-labeled *in situ* probes are used for the detection of multi-copy or very large hybridization targets on metaphase chromosomes or interphase nuclei.

## 2. How to Use this Product

# 2.1. Before you Begin

## **Sample Materials**

## **Templates for labeling reaction**

- Supercoiled and linearized plasmid DNA.
- Supercoiled and linearized cosmid DNA.
- Purified PCR products.
- *i* Denaturing of the template before nick translation is not required.

#### **General Considerations**

#### Labeling of in situ probes

For a standard labeling reaction using 1  $\mu$ g template in 20  $\mu$ l total reaction volume, 4  $\mu$ l of 5x-concentrated fluorophore labeling mix are required.

#### Labeling of PCR products or varying amounts of template

- When using higher amounts of template DNA (>1 μg) or PCR products as templates, analyze the fragment length distribution of the labeled probe before stopping the reaction with EDTA and heat.
- The total volume of the labeling reaction should be 20 μl/μg template containing 4 μl Nick Translation Mix/μg template.
- Purify PCR products with the High Pure PCR Product Purification Kit\* before use as nick translation templates to remove excess unincorporated dNTPs.

#### Removal of unincorporated free fluorophore-dNTP

For some *in situ* applications, it might be advantageous to remove unincorporated dNTPs. This can be performed by either ethanol precipitation, gel filtration using Quick Spin Columns for radiolabeled DNA purification, Sephadex G-50\*, or using the High Pure PCR Product Purification Kit\*.

#### **Nucleotide sets**

Use the following nucleotide set and fluorophore-dUTP conjugates for the labeling:

- Deoxynucleoside Triphosphate Set\*, containing individual vials of dATP, dCTP, dGTP and dTTP at a concentration of 100 mM.
- Fluorescein-12-dUTP\*, 25 µl, 1 mM
- Tetramethyl-Rhodamine-5-dUTP\*, 25 μl, 1 mM
- (i) Any other commercially available labeled nucleotides can be used in combination with the Nick Translation Mix for in situ probes.

# **Working Solution**

Solution	Composition	Use	Storage and Stability
EDTA	0.5 M ethylenediaminetetraacetic acid, pH 8.0	Stops the reaction.	+15 to +25°C

## Preparation of fluorophore labeling mix

For 50  $\mu$ l of a 5x-concentrated fluorophore labeling mix sufficient for 12 standard labeling reactions, combine in a sterile reaction tube:

Reagent	Volume [µl]	Concentration
dATP	5	1:40 dilution, 2.5 mM
dCTP	5	1:40 dilution, 2.5 mM
dGTP	5	1:40 dilution, 2.5 mM
dTTP	3.4	1:40 dilution, 2.5 mM
Fluorescein-12-dUTP or Tetramethyl-Rhodamine-5-dUTP	4	Undiluted, 1 mM
Add sterile, double-distilled water to a final volume of 50 $\mu$ l; mix thoroughly.	X	_

i If using fluorophore-dNTPs from other suppliers, it might be advantageous not to add the corresponding unlabeled dNTP.

## 2.2. Protocols

### **Standard labeling reaction**

- 1 The input DNA can be either linear or supercoiled. Purified PCR products can also be labeled.
- 1 In a reaction vial, carefully add the following reagents:

Reagent	Volume [μl]
1 μg template DNA in sterile, double-distilled water.	12
5x fluorophore labeling mix	4
Nick Translation Mix	4

- Mix and centrifuge briefly.
- 2 Incubate for 90 minutes at +15°C.
- 3 Stop the reaction by adding 1 μl 0.5 M EDTA, pH 8.0, and heat to +65°C for 10 minutes.

#### 2. How to Use this Product

## **Determination of fragment length of labeled probe**

- 1 After incubation at +15°C, place the reaction on ice.
- Remove a 3 μl aliquot per 20 μl reaction volume from the reaction, add gel loading buffer, and denature at +95°C for 3 minutes.
  - Place on ice for 3 minutes.
  - Run the sample on an agarose minigel along with a DNA Molecular Weight Marker.
  - The probe should range between 200 and 500 nucleotides in length.
- 3 If necessary, reincubate the reaction at +15°C and check the fragment size again.
- When correct probe length is achieved, stop the reaction by adding 1 μl 0.5 M EDTA, pH 8.0 per 20 μl reaction volume and heat to +65°C for 10 minutes.
  - *(i)* The labeled fragments obtained in the standard labeling reaction show a length distribution maximum in the range of 200 to 500 nucleotides.

Version: 10

## 3. Additional Information on this Product

# 3.1. Test Principle

The nick translation method is based on the ability of DNase I to introduce randomly distributed nicks into DNA at low enzyme concentrations in the presence of MgCl<sub>2</sub>.

*E. coli* DNA Polymerase I synthesizes DNA complementary to the intact strand in a  $5' \rightarrow 3'$  direction using the 3'-OH termini of the nick as a primer. The  $5' \rightarrow 3'$  exonucleolytic activity of DNA polymerase I simultaneously removes nucleotides in the direction of synthesis. The polymerase activity sequentially replaces the removed nucleotides with isotope-labeled or hapten-labeled deoxyribonucleoside triphosphates. At low temperature (+15°C), the unlabeled DNA in the reaction is thus replaced by newly synthesized labeled DNA.

#### In situ hybridization

In *in situ* hybridization experiments, the fragment length distribution severely influences the efficiency of hybridization. The use of probes showing fragment lengths above the optimal range of 200 to 500 nucleotides usually results in enhanced spotty background signals due to nonspecific sticking of the probe to the glass surface. It also can result in reduced accessibility to the target nucleic acid, such as metaphase chromosomes or cellular and tissue targets. On the other hand, use of excessively short probes will result in poor hybridization efficiency and sensitivity. This is due to fast rehybridization kinetics of short fragments yielding a high proportion of "snapback" probe DNA that reduces the amount of probe being available for hybridization to the target. Thus, the level of DNase I is of high relevance in probe labeling for *in situ* applications.

# 3.2. Quality Control

For lot-specific certificates of analysis, see section, Contact and Support.

# 4. Supplementary Information

## 4.1. Conventions

To make information consistent and easier to read, the following text conventions and symbols are used in this document to highlight important information:

Text convention and symbols			
1 Information Note: Additional information about the current topic or procedure.			
⚠ Important Note: Information critical to the success of the current procedure or use of the product.			
1 2 3 etc.	Stages in a process that usually occur in the order listed.		
1 2 3 etc.	Steps in a procedure that must be performed in the order listed.		
* (Asterisk)	The Asterisk denotes a product available from Roche Diagnostics.		

# 4.2. Changes to previous version

Layout changes. Editorial changes.

# 4.3. Ordering Information

Product	Pack Size	Cat. No.
Reagents, kits		
Fluorescein-12-dUTP	25 nmol, 25 μl, 1 mM	11 373 242 910
Tetramethylrhodamine-5-dUTP	custom fill	11 542 907 103
DNA Molecular Weight Marker VI	50 μg, 1 A <sub>260</sub> , 200 μl, 50 gel lanes	11 062 590 001
DNA Molecular Weight Marker VIII	50 μg, 1 A <sub>260</sub> , 200 μl, 50 gel lanes	11 336 045 001
High Pure PCR Product Purification Kit	1 kit, up to 50 purifications	11 732 668 001
	1 kit, up to 250 purifications	11 732 676 001
Quick Spin Columns for radiolabeled	20 columns	11 273 965 001
DNA purification	50 columns	11 273 973 001
Deoxynucleoside Triphosphate Set	4 x 250 μl, 4 x 25 μmol, 100 mM, 6,250 reactions at 20 μl final volume.	11 969 064 001
	4 x 1,250 μl, 4 x 125 μmo, 100 mM, 37,500 reactions at 20 μl final volume.	03 622 614 001

#### 4.4. Trademarks

HIGH PURE is a trademark of Roche.

All other product names and trademarks are the property of their respective owners.

#### 4.5. License Disclaimer

For patent license limitations for individual products please refer to: **List of biochemical reagent products**.

# 4.6. Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

# 4.7. Safety Data Sheet

Please follow the instructions in the Safety Data Sheet (SDS).

# 4.8. Contact and Support

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support Site**.

To call, write, fax, or email us, visit **sigma-aldrich.com**, and select your home country. Country-specific contact information will be displayed.