

ProductInformation

p-AMINOPHENYLMERCURIC ACETATE Sigma Prod. No. A9563

 NH_2 -(C_6H_4)-Hg-O-CO-C H_3

CAS NUMBER: 6283-24-5

SYNONYMS: APMA; (Acetato)(p-aminophenyl)mercury; p-(acetoxymercuri)aniline;

p-aminophenylmercuric acetate

PHYSICAL DESCRIPTION:

Appearance: white powder with a faint yellow cast

Molecular formula: C₈H₉HgNO₂

Formula weight: 351.8 Melting range 163-165°C¹

STORAGE / STABILITY AS SUPPLIED:

APMA is stable at room temperature for at least 3 years.²

SOLUBILITY / SOLUTION STABILITY:

APMA is soluble in DMSO at least at 10 mg/mL. In 100% acetic acid, a 50 mg/mL solution is clear light yellow.² Stock solutions at 10 mM in 0.05 N NaOH have been reported,³ but in water only about 5 mM is possible. Solutions of APMA should be prepared just before use.⁶

GENERAL USAGE:

The product is used as an activator of several of metalloproteases at 0.5 mM³ to 1.5 mM. Incubations were conducted at 37°C for 4 to 16 hours. 3,4,5,6,7,8,9

REFERENCES:

- 1. Sigma MSDS.
- 2. Sigma quality control.
- 3. Collier, I.E. et al., J. Biol. Chem., 263, 6579-6587 (1988).
- 4. Goldberg, G.I. et al, Proc. Natl. Acad. Sci., USA, 86, 8207-8211 (1989).
- 5. Nagase, H. et al., Biochem., 29, 5783-5789 (1990).
- 6. DeClerck, Y.A. et al., J. Biol. Chem., 266, 3893-3899 (1991).
- 7. Blaeser, J. et al., Eur. J. Biochem., 202, 1223-1230 (1991).
- 8. Okada, Y. et al., J. Biol. Chem., 267, 21712-21719 (1992).
- 9. Housley, T.J. et al., J. Biol. Chem., 268, 4481-4487 (1993).

Sigma warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.