Instructions for Use

Fetal Hemoglobin Kit

Procedure No. 285

Intended Use

Fetal Hemoglobin reagents are intended for use as a general laboratory histochemical stain for the acid elution determination of fetal hemoglobin in blood smears. Fetal Hemoglobin stain reagents are for "In Vitro Diagnostic Use." For professional use only. This manual, semi-quantitative procedure detects fetal hemoglobin in erythrocytes of human blood smear specimens.

Advances in techniques for protein isolation and characterization led to the discovery that there are several distinguishing properties that make it possible to differentiate fetal from adult hemoglobin. Among these is the resistance of fetal hemoglobin (hemoglobin F) to acid elution. When blood smears are immersed in acid buffer, for example, adult hemoglobin is eluted from the erythrocytes, whereas fetal hemoglobin is not. If blood smears are treated in this manner and subsequently stained, erythrocytes having hemoglobin F will take up the stain, while those containing only adult hemoglobin appear as "ghosts".

The slide technique for demonstrating fetal hemoglobin in terms of its resistance to acid elution was originally proposed by Kleihauer et al.,¹ and later modified by Shepard et al.² The Sigma produced represents a further improvement in this approach as described by Oski and Naiman.³

It is becoming increasingly common in cases of Rh incompatibility to suppress immune reactions to red blood cells entering maternal circulation from the fetus. The amount of specific gamma globulin, containing anti Rh(D) to be administered, is calculated by assessing the magnitude of fetal-maternal hemorrhage.⁴

According to described technique, blood smears, which have been properly dried and fixed, are immersed in a citrate buffer pH 3.3 at 37 °C. Adult hemoglobin A (HbA) dissolves out of the cells, whereas fetal hemoglobin (HbF) which is acid resistant, remains intracellular and can be stained for microscopic examination.

Reagents

Citrate Phosphate Buffer Concentrate (Cat. No. 2851)

Sodium citrate, 0.7 mol/L, and sodium phosphate, 0.6 mol/L, chloroform added as preservative

Acid Hematoxylin Solution (Cat. No. 2852)

Hematoxylin, certified, 1 g/L, C.I. 75290, aluminum ammonium sulfate, sodium iodates and stabilizers, pH 3.3

Eosin B Solution (Cat. No. 2853)

Eosin B, 0.1%, aqueous solution, C.I. 45400. Sodium azide, 0.1%, added as preservative

Special Materials Required but Not Provided

- Microscope
- Microscopes slides, cover slips
- Staining rack/Coplin jars
- Water bath, 37 °C
- Ethanol Fixative (Cat. No. 2858) 80% v/v ethyl alcohol

Storage and Stability

Store Citrate Phosphate Buffer Concentrate in refrigerator (2-8 °C). Discard if there is evidence of microbial growth.

Store Citrate Phosphate Buffer Solution in refrigerator (2-8 °C). Stable for 2 weeks. Use a fresh aliquot each day. Discard if there is evidence of microbial growth.

Store Acid Hematoxylin Solution and Eosin B Solution at room temperature (18-26 °C). Solutions may be reused if they are stored in tightly sealed staining jars in subdued lights.

Ethanol Fixative should be stored at room temperature. Store tightly sealed and as a flammable liquid. Solution may be reused but should be discarded if fixation is not adequate.

Deterioration

Discard Acid Hematoxylin Solution when the time required for suitable staining exceeds 8 minutes.

Preparation

Citrate Phosphate Buffer Solution is prepared by diluting 1 volume of Citrate Phosphate Buffer Concentrate with 9 volumes of water.

Acid Hematoxylin Solution, Eosin B Solution and Ethanol Fixative are ready to use.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

All samples must be clearly labeled. Suitable instruments must be used for obtaining and preparing samples. Follow manufacturer's instructions for use.

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Either capillary or venous blood may be used. Capillary blood may be transferred directly to a clean microscope slide. Venous blood should be added to a tube containing EDTA or oxalate. For convenience, use 1-2 drops of 2% EDTA Solution (Cat. No. 2854) per mL of blood (1 drop = 1 mg). Although blood-EDTA mixtures have been reported to be satisfactory for use up to 2 weeks when refrigerated.² Smears should be prepared within 24 hours from blood collected in oxalate.⁵ Using samples from the newborn, it is recommended that the blood be diluted with 0.85% saline, since such specimens have a high content of HbF. Blood smears are not stable and must be tested immediately after preparation.

Notes

For quality control purposes, it is recommended that blood from a normal adult (HbA) and from a newborn or infant (HbF) be included in each series of tests. Barr and Shafer⁶ report that fixed positive control slides from cord blood in EDTA can be preserved up to 1 year at -20 °C, in a sealed cardboard box. However, EDTA negative controls do not elute completely after being stored for more than 2 months at -20 °C. These investigators suggest preparing positive and negative smears on the same slide; thus, providing clear and rapid contrasts as reference in reading test slides.

Normal Ranges⁷

Age	Fetal Hemoglobin (%)
At Birth	50-90
< 2 years	0-4
> 2 years	0-2

Procedure

1. Citrate Phosphate Buffer Solution should be warmed to 37 °C in a Coplin jar or staining dish.
2. Using clean, labeled microscope slides, make thin blood smears. Prepare CONTROL slides using positive HbF blood (cord-blood) and normal adult blood. Air dry approximately 10 minutes.
3. Fix slides by immersing in Ethanol Fixative for 5 minutes, rinse thoroughly with tap water and air dry.
4. Immerse TEST and CONTROL slides in pre-warmed Citrate Phosphate Buffer Solution at 37 °C for 5 minutes. Agitate after 1 and 3 minutes of immersion. Degree of agitation may be varied to achieve most desirable results. Rinse thoroughly with distilled water and air dry **completely** to avoid staining artifacts.
5. Stain the slides for 3 minutes in Acid Hematoxylin Solution. Rinse slides thoroughly with distilled water and shake off excess water.
6. Counterstain slides for 4 minutes in 0.1% Eosin B Solution. Rinse thoroughly with distilled water and air dry.
7. Place **dry** coverslip on slide and examine using oil immersion (1000X). The absence of HbF is evident by the presence of ghost cells while retained HbF causes cells to appear bright red. Do **not** apply oil directly to slide. (**Note:** The 400X magnification may be used, but the resulting larger field may be more difficult to count.)

Performance Characteristics

The proportion of erythrocytes containing fetal hemoglobin may be estimated several ways. When studying maternal blood evidence of HbF-containing cells, Oski and Naiman³ recommended the following:

1. Count total number of erythrocytes in 5 fields and determine the average number per field.
2. Count the number of deeply stained HbF-containing erythrocytes in about 30 fields and determine the average number per field.
3. Calculate percentage of HbF erythrocytes on the bases of the total number of erythrocytes per field.

Results are reported as the percent HbF present.

According to Oski and Naiman³ this method is capable of detecting as little as 0.1 mL of fetal blood in maternal circulation.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
2851	Citrate Phosphate Buffer Concentrate	Fetal Hemaglobin	3 of 3	3 of 3	3 of 3	3 of 3
2852	Acid Hematoxylin Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
2853	Eosin B Solution	Fetal Hemaglobin	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

285C, 285D:

H319 Causes serious eye irritation.

H350 May cause cancer.

P201 Obtain special instructions before use.

P264 Wash skin thoroughly after handling.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

- P280 Wear eye protection/ face protection.
- P305 + P351 + P338 IF IN EYES Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P308 + P313 IF exposed or concerned Get medical advice/attention.
- P337 + P313 If eye irritation persists Get medical advice/attention.
- P405 Store locked up.
- P501 Dispose of contents/container to an approved waste disposal plant.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

- Kleihauer E, Braun H, Betke K. Demonstration von fetalem Hämoglobin in den Erythrocyten eines Blutaussstrichs [Demonstration of fetal hemoglobin in erythrocytes of a blood smear]. Klin Wochenschr. 1957 Jun 15; **35(12)**:637-8. German. doi: 10.1007/BF01481043. PMID: 13450287.
- Shepard MK, Weatherall DJ, Conley CL: Semi-quantitative estimation of the distribution of fetal hemoglobin in red cell populations. Bull Johns Hopkins Hosp. **110**:293, 1962.
- Oski FA, Naiman JL: Hematologic Problems in the Newborn, 2nd ed. Saunders, Philadelphia, 1972, pp. 62-63.
- Technical Manual. American Association of Blood Banks, 9th ed. Arlington (VA), 1985, pp. 320-322.
- Clayton EM, Feldhaus WD, Phythyon JM: The demonstration of fetal erythrocytes in the presence of adult blood cells. Am J Clin Pathol. **40**:487, 1963.
- Barr JK, Shafer JA: Preparation and storage of control slides for fetal hemoglobin determination by acid elution. Am J Med Technol. **42**:54, 1976.
- Seiverd CE: Hematology for Medical Technologist, 4th ed. Lea & Febiger, Philadelphia, 1973, pp 511-512

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 4.0	2014
Rev. 5.0	2016
Rev. 6.0	2023
	Transferred to new template with current branding. Specified for professional use in intended use and precautions. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Updated literature references. Added a revision history table. Added Warnings and Hazards. Added CH-REP information.
Rev. 7.0	2025
	Updatedd SDS information. Added CH importer information.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

- Kleihauer E, Braun H, Betke K. Demonstration von fetalem Hämoglobin in den Erythrocyten eines Blutausstrichs [Demonstration of fetal hemoglobin in erythrocytes of a blood smear]. Klin Wochenschr. 1957 Jun 15; **35(12)**:637-8. German. doi: 10.1007/BF01481043. PMID: 13450287.
- Shepard MK, Weatherall DJ, Conley CL: Semi-quantitative estimation of the distribution of fetal hemoglobin in red cell populations. Bull Johns Hopkins Hosp. **110**:293, 1962.
- Oski FA, Naiman JL: Hematologic Problems in the Newborn, 2nd ed. Saunders, Philadelphia, 1972, pp. 62-63.
- Technical Manual. American Association of Blood Banks, 9th ed. Arlington (VA), 1985, pp. 320-322.
- Clayton EM, Feldhaus WD, Phythyon JM: The demonstration of fetal erythrocytes in the presence of adult blood cells. Am J Clin Pathol. **40**:487, 1963.
- Barr JK, Shafer JA: Preparation and storage of control slides for fetal hemoglobin determination by acid elution. Am J Med Technol. **42**:54, 1976.
- Seiverd CE: Hematology for Medical Technologist, 4th ed. Lea & Febiger, Philadelphia, 1973, pp 511-512

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- Kleihauer E, Braun H, Betke K. Demonstration von fetalem Hämoglobin in den Erythrocyten eines Blutaussstrichs [Demonstration of fetal hemoglobin in erythrocytes of a blood smear]. Klin Wochenschr. 1957 Jun 15; **35(12)**:637-8. German. doi: 10.1007/BF01481043. PMID: 13450287.
- Shepard MK, Weatherall DJ, Conley CL: Semi-quantitative estimation of the distribution of fetal hemoglobin in red cell populations. Bull Johns Hopkins Hosp. **110**:293, 1962.
- Oski FA, Naiman JL: Hematologic Problems in the Newborn, 2nd ed. Saunders, Philadelphia, 1972, pp. 62-63.
- Technical Manual. American Association of Blood Banks, 9th ed. Arlington (VA), 1985, pp. 320-322.
- Clayton EM, Feldhaus WD, Pythyon JM: The demonstration of fetal erythrocytes in the presence of adult blood cells. Am J Clin Pathol. **40**:487, 1963.
- Barr JK, Shafer JA: Preparation and storage of control slides for fetal hemoglobin determination by acid elution. Am J Med Technol. **42**:54, 1976.
- Seiverd CE: Hematology for Medical Technologist, 4th ed. Lea & Febiger, Philadelphia, 1973, pp 511-512

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MDSS GmbH
MDSS CH GmbH
Chemie GmbH