



Product Information

INSULIN-LIKE GROWTH FACTOR-II (IGF-II) Human, Recombinant Expressed in *E. coli*

Product No. I 2139

Product Description

Insulin-like Growth Factor-II (IGF-II) was first isolated from human serum by Froesch et al.¹ as a factor displaying insulin-like activities that were not suppressed by antibodies to insulin. It had also been discovered that growth hormone-dependent factors in serum could stimulate the incorporation of ³⁵S into cartilage² and that calf serum factors induced cellular division in chick fibroblasts.³ In 1972, the term "somatomedin" was introduced in an unsuccessful attempt to unify the nomenclature of these factors.⁴ In 1987, a consensus among an international group of scientists endorsed the use of the terms insulin-like growth factors (IGF-I and IGF-II),⁵ originally proposed by Rinderknecht and Humbel.⁶ Hence, IGF-I and IGF-II have had several synonyms: nonsuppressible insulin-like activity (NSILA), sulfation factor activity (SFA), and multiplication stimulating activity (MSA). However, since IGF-II was not regulated by growth hormone, only IGF-I was known as a somatomedin.

Human IGF-II contains 67 amino acids and shares similar structural features with IGF-I, including a 62% sequence homology.⁷ In human plasma IGF-I and IGF-II are associated with IGF-binding proteins,^{8,9} which transport the polypeptides and partially regulate their actions *in vivo*.¹⁰ In addition to the insulin receptor, IGF-II binds to two forms of IGF receptors, both of which are widely distributed in different tissues and cultured cells.¹¹ IGF-II is mitogenic for a variety of cultured cells, including mouse 3T3 cells,¹² normal rat kidney cells,¹³ human or chicken fibroblasts,^{14,15} and MCF-7 human breast carcinoma cells.¹⁶

Performance Characteristics

The mitogenic activity of IGF-II on a breast carcinoma cell line, MCF-7, has been tested in a serum-free assay system using a modification of the procedure of Karey and Sirbasu.¹⁶ The EC₅₀ is defined as the effective concentration of growth factor that elicits a 50% increase in cell growth in a cell based bioassay.

Product Information

Expressed in *E. coli*
Molecular Weight: 7.5 kD
Purity: ≥97% by HPLC
EC₅₀: 0.5 - 25 ng/ml
Package Size: 10 µg
Formulation: Lyophilized from deionized water.
Carrier Protein: None
Sterility: sterile filtered
Endotoxin: <1 EU/vial

Reconstitution and Use

Reconstitute the contents of the vial in 10 - 100 mM acetic acid to prepare a stock solution. This may be diluted immediately before use in sterile medium of balanced salt solution containing 0.1% - 1.0% BSA or 1% - 10% serum.

Storage

Prior to reconstitution, store vial below 0 °C. After reconstitution, the product may be stored for two weeks at 2-8 °C or may be stored as aliquots at -20 °C. Prolonged storage and repeated freezing and thawing is **not** recommended.

References

1. Froesch, E., et al., J. Clin. Invest., **42**, 1816 (1963).
2. Salmon, W., et al., J. Lab. Clin. Med., **49**, 825(1957).
3. Pierson, R., et al., J. Cell. Physiol., **79**, 319 (1972).
4. Daughaday, W., et al., Nature, **235**, 107 (1972).
5. Daughaday, W., et al., J. Clin. Endocrinol. Metab., **65**, 1075 (1987).
6. Rinderknecht, E. and Humbel, R., Proc. Natl. Acad. Sci., **73**, 4379 (1976).
7. Froesch, E., et al., Ann. Rev. Physiol., **47**, 443 (1985).
8. Hardouin, S., et al., J. Clin. Endocrinol. Metab., **69**, 1291 (1989).
9. Zapf, J., et al., Arch. Biochem. Biophys., **168**, 638 (1975).
10. Baxter, R., and Martin, J., Prog. Growth Factor Res., **1**, 49 (1989).

11. Baxter, R., et al., Adv. Clin. Chem., **25**, 49 (1986).
12. Nissley, S., et al., J. Cell Physiol., **89**, 393 (1976).
13. Marquardt, H., et al., J. Biol. Chem., **255**, 9177 (1980).
14. Rechler, M., et al., J. Clin. Endocrinol. Metab., **44**, 820 (1977).
15. Rechler, M., et al., Nature, **259**, 134 (1976).
16. Kary, K. and Sirbasco, D., Cancer Res., **48**, 4803 (1988).

Pcs9/00

Sigma brand products are sold through Sigma-Aldrich, Inc.
Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.
Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.
Please see reverse side of the invoice or packing slip.