



## Product Information

### 3'-Amino-3'-deoxy-N<sup>6</sup>,N<sup>6</sup>-dimethyladenosine

Product Number **P 7130**

Storage Temperature 2-8 °C

#### Product Description

Molecular Formula: C<sub>12</sub>H<sub>18</sub>N<sub>6</sub>O<sub>3</sub>

Molecular Weight: 294.3

CAS Number: 58-60-6

Melting Point: 215-216 °C<sup>1</sup>

$\lambda_{\max}$ : 269 nm (0.1 M HCl)<sup>1</sup>

Extinction Coefficient:  $E^{1\text{cm}} = 18.6$  (0.1 mM HCl)<sup>1</sup>

Synonyms: puromycin aminonucleoside;

3'-amino-3'-deoxy-N,N-dimethyladenosine

3'-Amino-3'-deoxy-N<sup>6</sup>,N<sup>6</sup>-dimethyladenosine is the aminonucleoside portion of the antibiotic puromycin.<sup>1</sup> It has been used in nephrology research, such as studies of focal and segmental glomerulosclerosis and in the induction of nephrosis in rats.<sup>2-4</sup> The excretion of sodium and NO<sub>x</sub> metabolites in rats with puromycin aminonucleoside-induced nephrotic syndrome has been studied.<sup>5</sup> The generation of reactive oxygen species in rats in the acute phase of puromycin aminonucleoside induced nephrosis has been investigated.<sup>6</sup>

Puromycin aminonucleoside has been used to probe endothelial glycosaminoglycan synthesis in cultured glomerular endothelial cells and their relation to cell permeability.<sup>7</sup> The role of the neuron-specific ubiquitin C-terminal hydrolase protein gene product 9.5 (PGP 9.5) in rat kidney nephrogenesis has been studied in the puromycin aminonucleoside nephrosis model of rat glomerular disease.<sup>8</sup>

#### Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

#### Preparation Instructions

This product is soluble in water (50 mg/ml), with heat as needed, yielding a clear, colorless to faint yellow solution.

#### References

1. Baker, B. R., et al., Puromycin. Synthetic studies. VII. Partial synthesis of amino acid analogs. *J. Am. Chem. Soc.*, **77**, 1-7 (1955).
2. McDonough, A. A., and Biemesderfer, D., Does membrane trafficking play a role in regulating the sodium/hydrogen exchanger isoform 3 in the proximal tubule? *Curr. Opin. Nephrol. Hypertens.*, **12(5)**, 533-541 (2003).
3. Fogo, A. B., Animal models of FSGS: lessons for pathogenesis and treatment. *Semin. Nephrol.*, **23(2)**, 161-171 (2003).
4. Mattoo, T. K., and Kovacevic, L., Effect of grape seed extract on puromycin-aminonucleoside-induced nephrosis in rats. *Pediatr. Nephrol.*, **18(9)**, 872-877 (2003).
5. Ni, Z., and Vaziri, N. D., Downregulation of nitric oxide synthase in nephrotic syndrome: role of proteinuria. *Biochim. Biophys. Acta*, **1638(2)**, 129-137 (2003).
6. Gwinner, W., et al., Reactive oxygen species and antioxidant defense in puromycin aminonucleoside glomerulopathy. *J. Am. Soc. Nephrol.*, **8(11)**, 1722-1731 (1997).
7. Sorensson, J., et al., Synthesis of sulfated proteoglycans by bovine glomerular endothelial cells in culture. *Am. J. Physiol. Renal Physiol.*, **284(2)**, F373-380 (2003).
8. Shirato, I., et al., Protein gene product 9.5 is selectively localized in parietal epithelial cells of Bowman's capsule in the rat kidney. *J. Am. Soc. Nephrol.*, **11(12)**, 2381-2386 (2000).

GCY/RXR 12/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.