# Enhanced MALDI-TOF MS Spectra of Serum Peptides Using Ultrafiltration Aldo M. Pitt and Sara Gutierrez, Millipore Corporation, Life Sciences Division, Danvers, MA. USA

### Abstract

Proteomic analysis of complex samples, such as serum or plasma is frequently influenced by the presence of high protein concentrations that hinder peptide detection. These proteins suppress the ionization of native peptides during MALDI-TOF MS analysis. As a result, sample complexity reduction to lower the level of abundant proteins is rapidly becoming an essential first step of many proteomics analysis schemes. Several prefractionation strategies using chromatographic absorbents have been employed to remove abundant proteins such as albumin. As an alternative to adsorption chromatography, we have treated mammalian (human, murine, & bovine) serum or plasma samples with ultrafiltration (UF) membranes to produce relatively protein free filtrates. Experiments were performed to select the best molecular weight cutoff (MWCO) centrifugal device, establish the recoveries of spiked known peptides, and optimize the protocol. The UF filtrate was then acidified with 1% TFA and treated with reverse phase C18 and/or ion exchange ZipTip for de-salting and concentration. The results demonstrated that the 10K MWCO membrane gave the optimal results based on significantly improved detection of serum peptides in the 800-4000m/z range. The sample complexity reduction technology described provides a convenient and rapid method for the enhancement of native low molecular weight peptides in biological fluids such as serum or plasma.

## Introduction

Peptides and other low molecular weight molecules have been associated with many pathological states such as cancers, AIDS, diabetes, cardiovascular and neurological diseases (1-3). Many of the potential research & diagnostic applications of these "biomarkers" have not been fully realized because of their difficulty in analysis & detection. Many of these significant factors are often measured by radio-immunoassay protocols for R&D or diagnostic applications. Analyzing low molecular weight biomarkers (peptides) in serum and plasma has been notoriously difficult due to the vast number of contaminating salts, proteins, and lipids present. These are extremely problematic for both multidimensional liquid chromatography (MDLC) and mass spectroscopy techniques. High concentrations of proteins, lipids and salts suppress the ionization of native peptides during MALDI-TOF MS analysis. As a result, sample complexity reduction to lower the level of abundant proteins is rapidly becoming an essential first step of many proteomics analysis schemes. Ultrafiltration has been investigated and found to be unsuitable for the removal of abundant proteins like albumin for subsequent protein analysis (4). Several other pre-fractionation strategies using chromatographic absorbents have been more successfully employed to remove abundant proteins such as albumin prior to MDLC or electrophoresis. However, the low molecular weight filtrate fractions are typically difficult to analyze and are frequently lost in these other approaches. As an alternative to adsorption chromatography, we have investigated ultrafiltration (UF) techniques to produce relatively protein free filtrates and enhance the MALDI-TOF detection of serum peptides and other low molecular weight molecules.

# References

1. Basso, D., Valerio, A., Seraglia, R., Mazza, S., Piva, M.G., Greco, E., Fogar, P., Gall, N., Pedrazzoli, S., Tiengo, A., and Plebani, M. Putative Pancreatic Cancer-Associated Diabetogenic Factor: 2030 MW Peptide, Pancreas 2002 24:8-14.

2. Jungblu, P.R., Zimny-Arndt, U., Zeindl-Eberhart, E., Stulik, J., Koupilova, K., Pleißner, K-P., Otto, A., Müller, E-C., Sokolowska-Köhler, W., Grabher, G., Stöffler, G. Proteomics in human disease: Cancer, Heart and infectious diseases, Electrophoresis 1999 20:2100-2110.

3. Petricoin III, E., Ardekani, A.M., Hitt, B.A., Levine, P.J., Fusaro, V.A., Steinberg, S.M. Mills G.B. Simone C. Fishman D.A. Kohn F.C. and Liotta I.A. Use of proteomic patterns in serum to identify ovarian cancer. The Lancet 2002 359:572-577.


4. Georgiou HM, Rice GE, Baker MS. Proteomic analysis of human plasma: Failure of centrifugal ultrafiltration to remove albumin and other high molecular weight proteins. Proteomics 2001 Nov:1(12):1503-6.

#### Materials and Methods


Amicon® Ultra-4 centrifugal filter devices ( 10 K MWCO, cat #UFC801008), µC18 ZipTip (cat # ZTC18M096) and SCX ZipTip (Cat # STSCXS008) were purchased from Millipore Corporation (Bedford MA), Bovine, Mouse, and Human serum, as well as Brilliant blue G-colloidal Coomasie™, were purchased from Sigma Co. (St. Louis MO). The rest of the human plasma samples were taken from healthy male and female volunteer donors. Ficoll was added to blood samples and centrifuged at 500xg. Ficoll diluted plasma was separated and keep frozen at -20°C.

Trifluoroacetic acid (TFA), methanol, ammonium hydroxide and acetonitrile (ACN) were purchased from Fisher Co. (Pittsburg PA.). Alpha-cyano-4hydroxy cinnamic acid (CHCA) matrix was from Applied Biosystems (Foster City CA).10-20% tris-glycine SDS gels were purchased from Invitrogen Co. (Carlsbad CA)

Samples were analyzed in a linear Voyager-DE™ BioSpectrometry Workstation, Applied Biosystems (Framingham, MA), or in a reflectron Autoflex® Bruker Daltonics® (Billerica, MA.) MALDI-TOF Instruments




#### SDS-PAGE Gels of Original And Ultrafiltrate Serum Samples




Human, Bovine and Mouse sera were passed through Amicon Ultra devices using three different MWCO (10, 30 and 100K) membranes. It is clear that the 10K MWCO membrane produced a filtrate with the lowest protein concentration for all three sera (~99% rejection of the starting total protein concentration). All filtrate samples were further processed as above for MALDI evaluations.

#### Human Serum Flow Curve



Samples were spun at room temperature, in a swinging bucket rotor at 3000 x g for the times indicated for the 3 different molecular weight cut off devices. Both the Filtrate (F) and retentate (R) volumes are shown.



Percent frequency of the five most common peaks observed in

correspondent peptides reported in the literature on 10% fetal

Peptide (MW) Frequency(%

1668

1793

1949

2062

2191

76.5

94.7

68.4

82.3

0/ 7

the adult bovine serum reported in this poster and the

10% Fetal bovine serum (1) Adult Bovine serum

13.3

11.8

9.6

49.6

100

bovine serum (Basso, et. al., 2002, ref# 1).

Peptide (MW) Frequency(%

Original serum

Serum filtrate

A comparison of peptide

the filtrate.

patterns observed in Bovine

serum before (original) and in

1667

1795

1952

2065

2193

•A significantly increased number of low molecular weight molecules (800-4000 m/z) were detected in the filtrate compared to the starting serum or plasma sample. The actual identity of these observed low molecular weight molecules is still under investigation in addition to the significance, if any, of the more complex peptide patterns observed.

· A higher signal intensity and lower background were observed in the MALDI spectra obtained on the filtrate samples, compared with the starting material.

•The 10K dalton MWCO ultrafiltration devices were determined to be optimal to prepare low molecular weight fractions from serum or plasma for MALDI analysis

•All samples required ZipTip processing to concentrate and de-salt the samples prior to analysis

•The 5 most abundant peptides observed in bovine serum agreed with other literature reports. The results clearly demonstrates the increased observed frequencies after ultrafiltration for a number of peptide masses (1269, 1468, 1834 and 2211) while others were observed equally in both samples (such as 1793, 2191 and 1665).

•The sample complexity reduction technology described provides a convenient and rapid method for the enhancement and detection of native low molecular weight peptides in biological fluids such as serum or plasma.