

## Product Information

### Cholesterol Oxidase from microorganisms

Catalog Number **C8868**  
Storage Temperature  $-20\text{ }^{\circ}\text{C}$

CAS RN 9028-76-6


EC 1.1.3.6

Synonyms: Cholesterol:oxygen oxidoreductase;  
 $\beta$ -hydroxy steroid oxidoreductase; CHOD;  
 $3\beta$ -hydroxysteroid:oxygen oxidoreductase

#### Product Description

Cholesterol oxidase (CHOD) catalyzes the first step in cholesterol catabolism. Some non-pathogenic bacteria, such as *Streptomyces* are able to utilize cholesterol as a carbon source. Pathogenic bacteria, such as *Rhodococcus equi*, require CHOD to infect a host's macrophage.<sup>1</sup>

CHOD is bifunctional. Cholesterol is initially oxidized to cholest-5-en-3-one in an FAD-requiring step. The cholest-5-en-3-one is isomerized to cholest-4-en-3-one.<sup>1</sup> The isomerization reaction may be partially reversible.<sup>2</sup> The activity of CHOD depends on the physical properties of membrane to which the substrate is bound.<sup>3</sup> The net reaction is:



CHOD is used to determine serum cholesterol.<sup>4,5</sup> It is the second most widely used enzyme in diagnostic applications after glucose oxidase.<sup>6</sup> CHOD also finds application in the microanalysis of steroids in food samples and in distinguishing 3-ketosteroids from  $3\beta$ -hydroxysteroids.<sup>7</sup>

Transgenic plants expressing cholesterol oxidase are being investigated in the fight against the cotton boll weevil.<sup>8</sup> Cholesterol oxidase has also been used as a molecular probe to elucidate cellular membrane structures.<sup>3,9</sup>

Cholesterol oxidase is a monomeric flavoprotein containing FAD.<sup>1</sup>

Molecular mass:<sup>10</sup> 64 kDa


Cofactor:<sup>11</sup> FAD

Inhibitors:<sup>12</sup>  $\text{AgNO}_3$ ,  $\text{HgCl}_2$ , 0.1% SDS

Isoelectric point (pI):<sup>10</sup> 4.7


pH optimum:<sup>10</sup> 6.5 (see Figure 1)

**Figure 1.**  
pH Profile of Enzyme Activity




pH range:<sup>10</sup> Maximal activity retained from pH 5.3–7.5 (see Figure 2)

**Figure 2.**  
pH Profile of Enzyme Activity



Temperature optimum:<sup>10</sup>  $50\text{ }^{\circ}\text{C}$  (see Figure 3)

**Figure 3.**  
Temperature Profile of Enzyme Activity



This product is supplied as a lyophilized powder containing  $\geq 15\%$  protein with sucrose.

Specific activity:  $\geq 50$  units/mg protein

Unit definition: one unit will convert 1.0  $\mu$ mole of cholesterol to 4-cholesten-3-one per minute at pH 7.5 at 25 °C.

**Note:** 4-cholesten-3-one may undergo isomerization.

CHOD is assayed spectrophotometrically in a 3.0 ml reaction mixture containing 38 mM potassium phosphate, 0.009% (w/v) *o*-dianisidine, 0.017% (w/v) cholesterol, 0.33% (v/v) Triton™ X-100, 10 units of peroxidase, and 0.01–0.02 unit of cholesterol oxidase.

#### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### Preparation Instructions

CHOD is soluble in cold 50 mM potassium phosphate buffer, pH 7.5. Prepare solutions immediately before use.


#### Storage/Stability

Store product at  $-20$  °C with desiccation. When stored at  $-20$  °C, the enzyme retains activity for at least 9 months.

At 50 °C, no loss of activity is observed after 5 minutes; 15% loss of activity is observed after 60 minutes.<sup>10</sup> (see Figure 4)

**Figure 4.**

Temperature-dependent activity versus time



#### References

1. Caldinelli, L. *et al.*, Dissecting the structural determinants of the stability of cholesterol oxidase containing covalently bound flavin. *J. Biol. Chem.*, **280**, 22572-81 (2005).
2. Smith, A.G., and Brooks, C.J.W., The mechanism of the isomerization of cholest-5-en-3-one to cholest-4-en-3-one by cholesterol oxidase. *Biochem. Soc. Trans.*, **5**, 1088-90 (1977).
3. Ahn, K.-W., and Sampson, N.S., Cholesterol oxidase senses subtle changes in lipid bilayer structure. *Biochemistry*, **43**, 827-36 (2004).
4. Allain, C.C. *et al.*, Enzymatic determination of total serum cholesterol. *Clin. Chem.*, **20**, 470-75 (1974).
5. Lolekha, P.H. *et al.*, Performance of four sources of cholesterol oxidase for serum cholesterol determination by the enzymatic endpoint method. *Clin. Chim. Acta*, **339**, 135-45 (2004).
6. MacLachlan, J. *et al.*, Cholesterol oxidase: Sources, physical properties and analytical applications. *J. Steroid Biochem., Mol. Biol.*, **72**, 169-95 (2000).
7. Toyama, M. *et al.*, Alteration of substrate specificity of cholesterol oxidase from *Streptomyces* sp. by site-directed mutagenesis. *Protein Eng.*, **15**, 177-84 (2002).
8. Corbin, D.R. *et al.*, Expression and chloroplast targeting of cholesterol oxidase in transgenic tobacco plants. *Plant Physiology*, **126**, 1116-28 (2001).
9. Pal, R., *et al.*, Effect of cholesterol concentration on organization of viral and vesicle membranes. *J. Biol. Chem.*, **255**, 5802-06 (1980).
10. Supplier Data
11. Lee, S., *et al.* Purification and characterization of cholesterol oxidase from *Pseudomonas* sp. and taxonomic study of the strain. *Appl. Microbiol. Biotechnol.*, **31**, 542-546 (1989).
12. [http://brenda-enzymes.info/php/result\\_flat.php4?ecno=1.1.3.6](http://brenda-enzymes.info/php/result_flat.php4?ecno=1.1.3.6)

Triton is a trademark of The Dow Chemical Company or an affiliated company of Dow.

TD,CS,RGB,MAM 03/14-1