



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

SAPONIN FROM QUILLAJA BARK PURIFIED

Sigma Prod. No. S4521

Storage: Room Temperature

CAS NUMBER: 8047-15-2

SYNONYMS: Sapogenin Glycosides¹

STRUCTURE:

Quillaja saponaria saponin (Quillaja saponins) is a heterogenous mixture of molecules varying both in their aglycone and sugar moieties.² The main aglycone (sapogenin) moiety is quillaic acid, a triterpene of predominantly 30-carbon atoms (hydrophobic) of the Δ^{12} -oleanane type. The aglycone is bound to various sugars (hydrophilic) including glucose, glucuronic acid, galactose, xylose, apiose, rhamnose, fucose and arabinose. Sapogenin devoid of any sugars can be isolated by acid hydrolysis of saponins.^{3,4,5,6,7} The structure of a component isolated from the acylated triterpenoid saponin mixture of Quillaja saponaria was reported.⁸

PHYSICAL DESCRIPTION:

Appearance: powder³

Methods for the identification and quantitative determination of the aglycone and carbohydrate moieties of saponins have been reported.^{2,4,9,10,11}

METHOD OF PREPARATION:

Quillaja saponin is obtained from the bark of the South American soaptree, Quillaja saponaria Molina (Rosaceae family). Product S4521 is purified by ultrafiltration to reduce low molecular weight contaminants.³ A general method of preparation of the extract has been reported.²

STABILITY / STORAGE AS SUPPLIED:

Quillaja saponin is hygroscopic and is expected to be stable for at least one year at room temperature when stored dry.

SOLUBILITY / SOLUTION STABILITY:

Quillaja saponin is soluble in water yielding micelles with an average MW of 56,000. The solubility is tested at 50 mg/ml deionized water.³ The solubility in water may be increased by additions of small amounts of alkali.¹²

SAPONIN FROM QUILLAJA BARK PURIFIED
Sigma Prod. No. S4521

SOLUBILITY / SOLUTION STABILITY: (continued)

Quillaja saponin is soluble in hot alcohol and is insoluble in most organic solvents. Aqueous solutions will froth when shaken; the froth can be dispersed by alcohol or ether.¹² Quillaja saponin solutions are not autoclavable. Solutions have been found to be stable for about one month when stored at 2-8°C.³

USAGE / APPLICATIONS:

Saponin (plant source not identified) solutions have been used in the following ways: concentrations of 0.005%-0.01% were used to permeabilize cultured human intestinal epithelial cells after incubation for 30 minutes prior to immunofluorescence;¹³ a 0.05% solution of saponin was used to permeabilize paraformaldehyde-saponin-fixed human fibroblasts for staining of the intracellular fibronectin by peroxidase conjugates for light and electron microscopy;¹⁴ saponins (surface-active agents) have also been used to lyse the outer membranes of Rous sarcoma viruses, cell membranes of chicken liver and erythrocytes of human and guinea pig. Washed erythrocytes were suspended in isotonic saline and lysed by pouring into aqueous or saline solutions of saponin. Hemolysis was complete at 37°C for 30 minutes. At a saponin concentration of 0.05%, a large number of pits in the erythrocyte cell membrane was observed by electron microscopy whereby a 0.09% saponin concentration effected complete dissolution of the erythrocyte ghosts.¹⁵

Partially purified Quillaja saponins were reported to associate with hydrophobic or amphipathic proteins and lipids to form detergent/lipid/saponin complexes termed ISCOM (immunostimulating complexes). The ISCOMs, depending on the preparation, may induce serum antibody titers which are about 10-fold higher than titers produced by immunization with protein micelles alone.¹⁶ The isolation and quantification of purified Quillaja saponins and lipids in ISCOMS have been described.¹⁷ Orally fed Quillaja saponin was reported to enhance the immunopotentiating ability of an intraperitoneally administered inactivated rabies vaccine in mice.¹⁸ Purified saponin from crude Quillaja saponaria amplified antigen-specific immune responses to an experimental HIV-1 vaccine¹⁹ and potentiated an immune response elicited by albumin and venom both in mice.²⁰ Quillaja saponin has been shown to influence cholesterol metabolism and blood lipid profiles in animals. For casein-fed animals but not soy-fed, saponin significantly lowered the LDL cholesterol and LDL/HDL ratios but had no influence on the serum lipids of isolated soy protein (ISP)-fed animals. These results suggest that Quillaja saponins may interact differently with different proteins.²¹

GENERAL NOTES:

Saponins are plant glycosides which are widely distributed in plants. Saponins may also be powerful emulsifiers having hemolytic ability¹ (and therefore toxic) which may be related to their binding to cholesterol in cell membranes. Some saponins, like certain purified fractions of Quillaja saponaria, have been shown to have adjuvant activity.^{6,9,16-20,22} Other properties which may vary among the saponins include foaming when shaken in aqueous solutions, formation of molecular compounds with cholesterol and other hydroxy steroids, toxicity to fish and amphibians, formation of emulsions with oils, formation of mixed micelles with bile acids²³ and properties of a protective colloid and a bitter taste.^{4,5,24} Saponin has been used for the direct determination of metals in milk by flame atomic-absorption spectrophotometry.²⁵

SAPONIN FROM QUILLAJA BARK PURIFIED
Sigma Prod. No. S4521

REFERENCES:

1. *The Merck Index*, 12th Ed. #8513, 1996.
2. Dalsgaard, K., *Acta Vet. Scand.* 19, Suppl. 69, 7, (1978).
3. Sigma data
4. Birk, Y. "Saponins", Chapter 7, *Toxic Constituents of Plant Foodstuffs*, edited by Irvin E. Liener, p. 169, Academic Press, NY, 1969.
5. McGilroy, R.J., *The Plant Glycosides*, p. 64, Edward Arnold & Company, London.
6. Campbell, J.B. and Peerbaye, Y.A., *Research in Immunology*, 143, 526, (1992) (review).
7. Bondi, A. et al., "Forage Saponins", Chapter 11, *Chemistry and Biochemistry of Herbage*, edited by G.W. Butler and R. W. Bailey, Vol. 1, p. 511, Academic Press, NY, 1973.
8. Higuchi, R. et al., *Phytochem.*, 27, 1165, (1988).
9. Kensil, C.R. et al., *J. Immun.*, 146, 431, (1991).
10. Basu, N. and Rastogi, R.P., *Phytochem.*, 6, 1249, (1967) (review).
11. Chandel, R.S., and Rastogi, R.P., *Phytochem.*, 19, 1889, (1980) (review).
12. *Martindale, The Extra Pharmacopoeia*, 28th ed., p. 376, (1982).
13. Jalal, F. et al., *Biochem. J.*, 288, 945, (1992).
14. Hedman, K., *J. Histochem. Cytochem.* 28, 1233, (1980).
15. Dourmashkin, R.R. et al., *Nature*, 194, 1116, 1962.
16. Helenius, A. and von Bonsdorff, C.H., *Biochim. Biophys. Acta*, 436, 895, (1976).
17. Behboudi, S. et al., *Vaccine*, 13, 1690, (1995).
18. Chavali, S.R. et al., *Clin. Exp. Immunol.* 74, 339, (1988).
19. Wu, J.-Y. et al., *J. Immuno.* 148, 1519, (1992).
20. Cainelli Gebara, V.C.B. et al., *Biotechnol. Appl. Biochem.* 21, 31, (1995).
21. Potter, S.M. et al., *J. Agric. Food Chem.* 41, 1287, (1993).
22. Bomford, R. et al., *Vaccine*, 10, 572, (1992).
23. Oakenfull, D. *Aust. J. Chem.* 39, 1671, (1986).
24. Birk, T. and Peri, I. *Toxic Const. Plant Foodstuff*, I.E. Liener, ed. 161, Academic Press, New York, (1980).
25. Arpadjan, S. Stojanova, D., Fresenius Z., *Anal. Chem.* 302, 206, (1980).

Sigma warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.