



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

OVALBUMIN

Individual 2D Electrophoresis Marker

Product Number **O 4757**

Storage Temperature -70°C

Product Description

Ovalbumin is a protein with an molecular weight of approximately 45 kDa and an isoelectric point (pI) of approximately 5.1. Each vial contains approximately 200 μl of protein solution with 8 M urea and 2% 2-mercaptoethanol.

Preparation Instructions

Gels to be stained with Coomassie Blue: Add 4 μl of the desired Pharmalyte or Ampholine to the marker vial for a final ampholyte concentration of 2% (v/v). Vortex until the contents of the vial are completely mixed. Aliquot the unused portion and freeze at -70°C .

Gels to be silver stained: Dilute the marker 10-fold with sample buffer (sample buffer is prepared by mixing 800 μl of a 10 M urea solution, 20 μl of Triton X-100,

50 μl of Ampholine or Pharmalyte, and 50 μl of β -mercaptoethanol). Vortex until the contents of the vial are completely mixed. Aliquot the unused portion and freeze at -70°C .

Storage/Stability

It is recommended to store the product at -70°C . Repeated freezing and thawing is **not** recommended.

Pharmalyte is a registered trademark of Pharmacia, Inc. Ampholine and Immobiline are registered trademarks of Pharmacia LKB Biotechnology AB. Coomassie Brilliant Blue G is a trademark of Imperial Chemical Industries PLC.

MAM/KMR 02/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.