

Product Information

FLAG® Peptide

F3290

Storage Temperature 2-8 °C

Product Description

Molecular Weight: 1013.0

Sequence: N-Asp-Tyr-Lys-Asp-Asp-Asp-Asp-Lys-C (synthetic octapeptide)

The FLAG® Peptide is useful for the competitive elution of amino-terminal, Met-amino-terminal or carboxy-terminal FLAG® fusion proteins from the ANTI-FLAG® M1 or M2 monoclonal antibody in solution or bound to agarose gel. A working concentration of 100 µg/mL is commonly used to elute FLAG® fusion proteins from the ANTI-FLAG® M1 and M2 affinity resins.²⁻⁴ Five one-column volumes of the working solution are sufficient to elute most FLAG® fusion proteins.

Various publications have used the FLAG® Peptide (Cat. No. F3290) at other concentrations, which we have not necessarily tested, such as:

- 150 µg/mL⁵
- 200 µg/mL⁶
- 250 µg/mL⁷
- 300 µg/mL^{8,9}
- 0.5 mg/mL (500 µg/mL)¹⁰

The FLAG® Peptide will not elute 3X FLAG® fusion proteins. For this application, use the 3X FLAG® Peptide (Cat. No. F4799).

Precautions and Disclaimer

For research use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

To prepare a stock solution, dissolve in TBS (10 mM Tris-HCl, 150 mM NaCl, pH 7.4) to a final concentration of 5 mg/mL. Aliquot and store at -20 °C. Repeated freezing and thawing is **not** recommended.

Storage/Stability

Store the lyophilized powder at 2-8 °C. After reconstitution, aliquot and store at -20 °C. Repeated freezing and thawing is **not** recommended.

Procedure

Elution of FLAG® Fusion Proteins by Competition with FLAG® Peptide: Elute the bound FLAG® fusion protein from the ANTI-FLAG® M1 or M2 affinity resin by competitive elution with five one-column volumes of a solution containing 100 µg/mL FLAG® peptide in TBS.

References

1. Brizzard, B.L. *et al.*, *BioTechniques*, **16(4)**: 730-735 (1994).
2. Antonysamy, S. *et al.*, *Proc. Nat. Acad. Sci. USA*, **109(44)**: 17960-17965 (2012).
3. He, J. *et al.*, *Nucleic Acids Res.*, **40(13)**: 6097-6108 (2012).
4. Singh, O. *et al.*, *Protein Sci.*, **22(8)**: 1071-1077 (2013).
5. Li, W. *et al.*, *Mol. Cell. Biol.*, **36(22)**: 2855-2866 (2016).
6. Cheng, F. *et al.*, *Mol. Immunol.*, **60(1)**: 44-53 (2014).
7. Singh, G. *et al.*, *Methods*, **65(3)**: 320-332 (2014).
8. Meller, N. *et al.*, *J. Biol. Chem.*, **279(36)**: 37470-37476 (2004).
9. Aguado, L.C. *et al.*, *Nature*, **547(7661)**: 114-117 (2017).
10. Simsek, D. *et al.*, *Cell*, **169(6)**: 1051-1065.18 (2017).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.
© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.