

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

# **Product Information**

# 27672/27670 Collagenase Substrate Kit (for quantitative Collagenase-Determination)

#### Components:

Component A (Sigma 27673):

Collagenase-Substrate: Carbobenzoxy-Gly-Pro-Gly-Gly-Pro-Ala-OH

Molecular Weight 588.61 g/mol CAS Number 13075-38-2

Component B (Sigma 27674):

Standard: Gly-Pro-Ala

Molecular Weight 243.26 g/mol

CAS Number 837-83-2

#### **Enzymatic Assay of Collagenase**

The component A (syntethic substrate) is hydrolysed by the collagenase to Carbobenzoxy-Gly-Pro-Gly and Gly-Pro-Ala (component B) which reacts with ninhydrin and can be determined with spectroscopic method.

#### Principle:

Carbobenzoxy-Gly-Pro-Gly-Gly-Pro-Ala (= Component A)

Carbobenzoxy-Gly-Pro-Gly
+ Gly-Pro-Ala (= Component B)

Abbreviations used: Z = Carbobenzoxy

Conditions: Incubation T = 37°C, Incubation volume 2.5 ml, pH = 6.3, Measure T = 25°C, Measure volume 5.0 ml,

 $A_{565nm}$ , Light path = 1 cm

Method: Spectrophotometric Stop Rate Determination

#### Reagents:

		Sigma No.
1.	Citric Acid anhydrous, BioChemika Ultra	27487
2.	Calcium acetate, BioChemika Ultra	21056
3.	Sodium hydroxide, MicroSeiect	S8045
4.	Trichloric Acid, BioChemika Ultra	91228
5.	Collagenase	
6.	Ninhydrin	151173
7.	2-Methoxyethanol	33457
8.	Tin(II) chloride dihydrate, puriss p.a.	31669
9.	Sodium acetate, anhydrous, puriss.p.a.	32319
10.	Acetic acid puriss.p.a.	33209
11.	1-Propanol	402893
12.	Argon	
13.	Z-Gly-Pro-Gly-Gly-Pro-Ala Chromophore-Substrate (Component A)	27673
14.	Gly-Pro-Ala Standard (Component B)	27674

#### Solutions:

- A) Citrate Buffer; 0,1 mol/l, pH 6,3 contains Calcium acetate 0,01 mol/l
  - Prepare by dissolving 9.62 g Citric Acid and 0.79 g Calcium Acetate in ca. 300 ml H<sub>2</sub>O distilled. Adjust the pH with NaOH >2 mol/l to 6.3. Bring to a total volume of 500 ml with H<sub>2</sub>O distilled.
- B) Calcium acetate Solution; 0,01 mol/l:
  - Dissolve 158 mg Calcium acetate in 100 ml H<sub>2</sub>O distilled
- C) Collagenase-Enzyme-Solution; about 30 U/ml
  - 0.68 mg Collagenase in 25 ml Calcium acetate Solution (B)
- D) Trichloric Acid Solution; 1 mol/l:
  - Dissolve 1.63 g in 25 ml with H<sub>2</sub>O distilled.
- E) Sodium Acetate Buffer; 4 mol/l, pH 6.3:
  - Prepare by dissolving 16.4 g Sodium acetate in ca. 30 ml H<sub>2</sub>O distilled and adjust pH to 6.3 with acetic acid.
     Bring to a total volume of 50 ml with H<sub>2</sub>O distilled.
- F) Ninhydrin Reagent-Solution:
  - Dissolve 2.0 g Ninhydrin in 75 ml 2-Methoxyethanol, add 40 mg Tin(II) chloride dihydrate and bing to 100 ml with Sodium Acetate Buffer (E). It is a orange-yellow solution. Blanket solution with Argon and store about one day in the fridge. Colour of solution change to yellow.
- G) Propanol/H<sub>2</sub>O distilled; 50%
  - Bring 25 ml Propanol to an end volume of 50 ml with H<sub>2</sub>O distilled
- H) Collagenase-Substrate (Componente A) Solution; Z-Gly-Pro-Gly-Pro-Ala-OH, 0.233 mg/ml
  - Prepare by dissolving 11.5 mg of Z-Gly-Pro-Gly-Gly-Pro-Ala-OH, Sigma 27673, in 50 ml of Citrate Buffer (A).
- I) Collagenase-Substrate-Standard (Componente B) Solution; Gly-Pro-Ala, 125 μg/ml
  - Dissolve the different concentration in Citrate Buffer (A)
    - 2.500 mg in 20 ml
- $\rightarrow$  30 µg/ 5 ml in Test

# Stability and Storage of Solutions:

Store all solutions at 4°C. Solution A is about 1 week, solutions B, C, D, E and G for some weeks stable. Solution F should be prepared at least one day before use. Solution H and I should be prepared immediately before use.

#### Procedure:

## Step 1:

Pipette (in milliliters) the following reagents into suitable containers:

•	Test	Blank	Standard
Substrate Solution	2.0	2.0	-
(Component A)			
Standard Solution	-	-	2.0
(Component B)			
Calcium acetate	-	0.5	0.5
Collagenase Solution	0.5	-	-

Immediately mix by swirling and incubate for exactly 5 minutes at 37°C. After incubation give 0.3 ml into a 5 ml graduated flask.

#### Step 2:

Pipette (in milliliters) the following reagents into suitable containers:

	Test	Blank	Standard
Mixture from Step 1	0.3	0.3	0.3
Citrate Buffer	0.1	0.1	0.1

Mix immediately after addition of the buffer.

## Step 3:

Pipette (in milliliters) the following reagents into suitable containers:

	Test	Blank	Standard
Mixture from Step 2	0.4	0.4	0.4
Ninhydrin Reagent-	1.0	1.0	1.0
Solution			

Mix well and heat 20 min. at  $100^{\circ}$ C. Cool in water bath to room temperature and bring to 5 ml total volume with Propanol/H<sub>2</sub>O solution. Transfer the solutions to suitable cuvettes and record the A<sub>565nm</sub> using a suitable spectrophotometer.

# Calculations:

Units / 
$$ml = \frac{(A_{565nm}Test - A_{565nm}Blank)2.5 \cdot df}{5 \cdot 3.784 \cdot 0.3 \cdot 0.5}$$

2.5 [ml] = Total volume of assay in Step 1

df = Dilution factor

5 [min.] = Incubation time

3.784 [1/µmol] = extinction coefficient of 1 micromol of Gly-Pro-Ala in 5ml and 1 cm light path

 $(\epsilon_{565nm} = 18.92 \text{ cm}^2/\mu\text{mol})$ 

0.3 [ml] = Volume of reaction mixture in Step 1 used in Step 2+3

0.5 [ml] = Volume of enzyme used in Step 1 units/ml enzyme

Units/mg solid = 
$$\frac{\text{mg solid}}{\text{ml enzy me}}$$

#### **Units Definition:**

One unit will liberate 1  $\mu$ mol of Gly-Pro-Ala from Z-Gly-Pro-Gly-Gly-Pro-Ala (Sigma 27673) in 1 minute at pH 6.3 at 37°C.

*Note:* The standard is to check the extinction coefficient ( $\varepsilon_{565nm} = 18.92 \text{ cm}^2/\mu\text{mol}$ ).

#### References:

1. W. Grassmann, A. Nording, Z. Physiol. Chemie 322, 267 (1960)

#### **Precautions and Disclaimer:**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser mus determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invo or packing slip.