


Novagen®

iFOLD™

## Protein Refolding Systems

Take the guesswork out of finding  
the optimal refolding condition  
for your target protein!



# iFOLD™ Protein Refolding Systems at a Glance

Protein functional and structural studies often require a large amount of pure, correctly folded protein, which is commonly produced in *Escherichia coli* (*E. coli*) expression systems. However, production of foreign proteins in *E. coli* can result in the formation of inclusion bodies (IB) – dense, insoluble, aggregates of misfolded protein. IB also have useful attributes – they are easily purified, resistant to proteolysis, and can be solubilized with chaotropic agents. Defining conditions that promote refolding of a chemically denatured protein into its native conformation is empirical and often time consuming. Simultaneous and systematic evaluation of a large number of refolding conditions increases chances of identifying an optimal refolding condition for a given protein.

We have taken the tedium and guesswork out of finding the optimal refolding condition for your target protein! Our new iFOLD™ Protein Refolding Systems provide comprehensive refolding screens in a 96-well plate format and are based on an extensive literature review of successful refolding experiments and information contained in the REFOLD database (<http://refold.med.monash.edu.au>). The systems differ in the chemistry used to denature the inclusion bodies and in the refolding additives included in the 96-well plate. In addition to a 96-well plate containing 92 (System 1) or 95 (System 2) unique buffer and refolding additive combinations, both systems include the reagents needed to purify inclusion bodies and solubilize the component proteins. Significantly, all steps of the refolding screens are equally compatible with manual and high-throughput automated liquid handling.

## Comparison of components in iFOLD Protein Refolding Systems 1 and 2

	iFOLD™ System 1	iFOLD™ System 2
IB wash agent	1.0% Triton X-100	0.125 M NDSB-201
IB denaturant	4.4% N-lauroylsarcosine	7.0 M GuHCl
IB reducing agent	5 mM TCEP	10 mM TCEP
Denatured target protein: required concentration required volume	1 mg/ml 5 ml	5 mg/ml 1 ml
Total Protein per well	50 µg	50 µg
Buffer system (50 mM)	Tris-HCl, pH 7.0, 7.5, 8.0 & 8.5	MOPS, pH 7.0 HEPES, pH 7.5 EPPS, pH 8.0 TAPS, pH 8.5 CHES, pH 9.0
Refolding additives	100 mM NaCl	24 mM NaCl + 1.0 mM KCl
	250 mM NaCl	240 mM NaCl + 10 mM KCl
	1.0 mM TCEP	1.0 mM TCEP
	3.8 mM GSH + 1.2 mM GSSG	9.0 mM GSH + 1.0 mM GSSG
		6.0 mM GSH + 4.0 mM GSSG
	1.0 mM EDTA	1.0 mM EDTA
	1.0 mM each of CaCl <sub>2</sub> and MgCl <sub>2</sub>	0.25 mM each of CaCl <sub>2</sub> , MgCl <sub>2</sub> , MnCl <sub>2</sub> and ZnCl <sub>2</sub>
	0.5 M GuHCl	0.5 or 1.0 M NDSB-201
		0.5 or 1.0 M NDSB-256
	20% glycerol	1.5 M sorbitol
		0.58 M trehalose
	0.1% PEG6000	0.06% PEG3350
	0.5 M L-arginine	0.5 M L-arginine
12.5 mM methyl-β-D-cyclodextrin	10 mM methyl-β-D-cyclodextrin	

Buffer and additive concentrations are final, after addition of the target protein

Abbreviations: TCEP: tris(2-carboxyethyl)phosphine, GSH: reduced glutathione, GSSG: oxidized glutathione, EDTA: ethylenediaminetetraacetic acid, NDSB: non-detergent sulfobetaine, Sorbitol: D-sorbitol or D-glucitol, Trehalose: 1-O-α-D-Glucopyranosyl-α-D-glucopyranoside, PEG6000: polyethyleneglycol, average molecular weight of 6000 Da, PEG3350: polyethyleneglycol, average molecular weight of 3350 Da



Prices and availability are subject to change. ©Copyright 2007 EMD Chemicals, an affiliate of Merck KGaA, Darmstadt, Germany. All rights reserved. Each product is sold with a limited warranty which is provided with each purchase. Each product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. EMD Chemicals products may not be resold, modified for resale, or used to manufacture commercial products without written approval of EMD Chemicals. NOVAGEN® is a registered trademark of EMD Biosciences in Australia, the European Union, Japan, Switzerland, and the United States. D-Tube™, D-Tube96™, IB-Prep™, iFOLD™, Lysonase™ and Perfect Protein™ are trademarks of EMD Biosciences. BODIPY® is a registered trademark of Molecular Probes, Inc., Brij® is a registered trademark of ICI Americas Inc., Triton® is a registered trademark of Dow Chemical Company.

# iFOLD™ Protein Refolding Systems

## Brief Protocol

- ↓ Express target protein
- ↓ Harvest cells
- ↓ Lyse cells by sonication plus Lysonase™ Bioprocessing Agent
- ↓ Centrifuge to pellet inclusion bodies
- ↓ Wash inclusion bodies with Triton® X-100 (System 1) or NDSB-201 (System 2)
- ↓ Denature inclusion bodies with N-lauroylsarcosine (System 1) or guanidine-HCl or urea (System 2)
- ↓ Refold target protein by rapid dilution into iFOLD™ Protein Refolding System 1 or System 2 buffer matrix
- ↓ Assay\* for correctly folded and active target protein

\*Assay is determined by end user

## iFOLD™ Protein Refolding System 1

- All reagents for inclusion body purification and pre-dispensed 96-well plate-based protein refolding matrix
- Uses N-lauroylsarcosine, a chaotropic anionic detergent, to denature the purified inclusion bodies
- 92 unique buffer and refolding additive combinations for simultaneous and systematic evaluation of protein refolding conditions
- pH range 7.0 – 8.5
- Refolding additives include salts, cyclodextrin, redox agents, chaotropes, glycols
- Refolding conditions based on extensive literature review and REFOLD database (<http://refold.med.monash.edu.au>)
- High-throughput compatible

## iFOLD™ Protein Refolding System 2

- Entirely detergent-free system 2 contains inclusion body purification reagents and pre-dispensed 96-well plate-based protein refolding matrix
- Uses guanidine hydrochloride or urea (not included) to denature the purified inclusion bodies
- 95 unique buffer and refolding additive combinations for simultaneous and systematic evaluation of protein refolding conditions
- pH range 7.0 – 9.0
- Refolding additives include salts, redox agents, cyclodextrin, chaotropes, glycols, nondetergent sulfobetains (NDSBs)
- Refolding conditions based on extensive literature review and REFOLD database (<http://refold.med.monash.edu.au>)
- High-throughput compatible



Cat. No. 71719-3 [www.novagen.com/ifoldsystem2](http://www.novagen.com/ifoldsystem2)

	1	2	3	4	5	6	7	8	9	10	11	12
A	MOPS cyclodextrin EDTA GSH + GSSG (1) NaCl + KCl (1)	MOPS L-Arg EDTA GSH + GSSG (2) NaCl + KCl (2)	MOPS NDSB 256 (1) metals TCIP NaCl + KCl (2)	HEPES trehalose EDTA GSH + GSSG (1) NaCl + KCl (1)	HEPES sorbitol EDTA GSH + GSSG (2) NaCl + KCl (2)	EPPS NDSB 256 (1) GSH + GSSG (1)	NDSB 201 (1) EDTA GSH + GSSG (2)	EPPS cyclodextrin EDTA GSH + GSSG (2) NaCl + KCl (2)	TAPS sorbitol EDTA GSH + GSSG (1) NaCl + KCl (1)	TAPS sorbitol EDTA GSH + GSSG (2) NaCl + KCl (2)	CHES sorbitol EDTA GSH + GSSG (1) NaCl + KCl (1)	CHES sorbitol EDTA GSH + GSSG (2) NaCl + KCl (2)
B	MOPS L-Arg GSH + GSSG (1) NaCl + KCl (1)	MOPS NDSB 256 (1) EDTA GSH + GSSG (2) NaCl + KCl (1)	MOPS NDSB 201 (2) metals TCIP NaCl + KCl (1)	HEPES NDSB 201 (1) EDTA GSH + GSSG (1) NaCl + KCl (2)	HEPES EDTA GSH + GSSG (2)	EPPS sorbitol EDTA GSH + GSSG (1)	EPPS NDSB 256 (2) EDTA GSH + GSSG (2) NaCl + KCl (1)	EPPS metals TCIP NaCl + KCl (1)	TAPS NDSB 201 (2) EDTA GSH + GSSG (1) NaCl + KCl (1)	TAPS sorbitol EDTA GSH + GSSG (2) NaCl + KCl (2)	CHES NDSB 256 (2) EDTA GSH + GSSG (1) NaCl + KCl (2)	CHES sorbitol EDTA GSH + GSSG (2) NaCl + KCl (2)
C	MOPS PEG 3350 EDTA GSH + GSSG (1)	MOPS cyclodextrin EDTA GSH + GSSG (2) NaCl + KCl (2)	MOPS NDSB 256 (2) EDTA GSH + GSSG (1) NaCl + KCl (2)	HEPES NDSB 256 (2) EDTA GSH + GSSG (1) NaCl + KCl (2)	HEPES NDSB 256 (1) EDTA GSH + GSSG (2)	EPPS sorbitol EDTA GSH + GSSG (1)	EPPS NDSB 201 (2) EDTA GSH + GSSG (2) NaCl + KCl (1)	EPPS trehalose EDTA GSH + GSSG (2) NaCl + KCl (1)	TAPS L-Arg EDTA GSH + GSSG (1) NaCl + KCl (1)	TAPS NDSB 256 (1) EDTA GSH + GSSG (2) NaCl + KCl (2)	CHES NDSB 201 (1) EDTA GSH + GSSG (1) NaCl + KCl (1)	CHES cyclodextrin EDTA GSH + GSSG (2) NaCl + KCl (2)
D	MOPS EDTA GSH + GSSG (1) NaCl + KCl (1)	MOPS NDSB 201 (1) EDTA GSH + GSSG (2)	HEPES sorbitol EDTA GSH + GSSG (1) NaCl + KCl (1)	HEPES NDSB 256 (2) EDTA GSH + GSSG (1) NaCl + KCl (2)	HEPES PEG 3350 EDTA GSH + GSSG (2)	EPPS NDSB 256 (1) EDTA GSH + GSSG (1)	EPPS NDSB 201 (2) EDTA GSH + GSSG (2) NaCl + KCl (1)	EPPS trehalose EDTA GSH + GSSG (2) NaCl + KCl (1)	TAPS L-Arg EDTA GSH + GSSG (1) NaCl + KCl (1)	TAPS NDSB 201 (1) EDTA GSH + GSSG (2) NaCl + KCl (2)	CHES NDSB 256 (1) EDTA GSH + GSSG (1) NaCl + KCl (1)	CHES sorbitol EDTA GSH + GSSG (2) NaCl + KCl (2)
E	MOPS trehalose EDTA GSH + GSSG (1)	MOPS NDSB 201 (2) EDTA GSH + GSSG (2)	HEPES NDSB 201 (1) EDTA GSH + GSSG (1) NaCl + KCl (1)	HEPES L-Arg EDTA GSH + GSSG (1) NaCl + KCl (1)	HEPES L-Arg EDTA GSH + GSSG (2)	EPPS PEG 3350 EDTA GSH + GSSG (1)	EPPS trehalose EDTA GSH + GSSG (2) NaCl + KCl (1)	EPPS sorbitol EDTA GSH + GSSG (2) NaCl + KCl (1)	TAPS sorbitol EDTA GSH + GSSG (1) NaCl + KCl (1)	TAPS cyclodextrin EDTA GSH + GSSG (2) NaCl + KCl (2)	CHES NDSB 201 (1) EDTA GSH + GSSG (1) NaCl + KCl (1)	CHES NDSB 256 (2) EDTA GSH + GSSG (2) NaCl + KCl (2)

Components of each well are shown; the pH is indicated below. Concentrations of other components can be found in the table on the reverse side of this card.

MOPS, pH 7.0 ●  
 HEPES, pH 7.5 ●  
 EPPS, pH 8.0 ●

# Preparation of Inclusion Bodies


Typically, inclusion bodies (IB) are highly enriched for the target protein but contain some contaminating biomolecules (nucleic acids, phospholipids, lipopolysaccharides, other proteins). Usually these contaminants are reduced by washing the IB pellet with a detergent, such as Triton® X-100 or sodium deoxycholate. The iFOLD™ System 1 uses 1.0% Triton X-100 for IB washing step.

When refolded proteins are used in downstream applications that can not tolerate even trace amounts of detergent, we recommend the iFOLD System 2, which uses a non-detergent sulfobetaine, NDSB-201, to wash the IB pellet. NDSB-201 is a 201 Da zwitterionic compound that does not form micelles and is easily removed by dialysis.

## Comparison of a thioredoxin-GFP fusion protein from inclusion bodies washed with NDSB-201 or Triton® X-100


The trx-GFP fusion protein was expressed in *E. coli* strain BL21(DE3) using Overnight Express™ Instant TB medium. Cells were lysed by sonication and IB were processed according to the iFOLD 2 protocol, using 1.0% Triton X-100

(lanes 2-7) or 0.125 M NDSB-201 (lanes 9-14) for the IB washes. Samples from each step were assayed by SDS-PAGE (10-20% gradient gel) and detected with Coomassie blue.



→ NDSB-201 comparable to Triton X-100

## Refolded HRV 3C protease activity from IB washed with NDSB-201 or Triton X-100



HRV 3C IB were washed with NDSB-201 or Triton X-100, using the iFOLD System 2 protocol. Washed inclusion bodies were denatured with 7.0 M GuHCl and diluted into the iFOLD Protein Refolding Plate 2. After refolding for 20 h at 22°C, reactions were dialyzed in D-Tube96™ Dialyzers against 2 x 4.0 L buffer (20 mM Tris-HCl, pH 8.0, 150 mM NaCl, 5.0 mM DTT, and 0.03% Brij®-35) overnight at 10°C. Enzymatic activity of the refolded and dialyzed proteins was quantified by measuring cleavage of a peptide substrate (Glu-Ala-Leu-Phe-Gln-pNa: Bachem L-2050). A. Graph of HRV 3C activity for individual wells from the refolding screen. B. Data grouped by refolding buffer pH.


↓ NDSB-201 comparable to Triton X-100

## Denaturing Inclusion Bodies, System 1

With the iFOLD™ Protein Refolding System 1, the purified inclusion bodies are denatured by addition of TCEP and N-lauroylsarcosine and

refolded by rapid dilution into the iFOLD Protein Refolding Plate 1.

### Refolded HRV 3C protease activity from IB denatured with N-lauroylsarcosine



Enzymatic activity of HRV 3C protease (22 kDa) refolded at 22°C was determined by measuring  $\Delta A_{405}/\text{min}$ , indicating cleavage of a *p*-nitroanilide-labeled peptide substrate (Wang 1997).

Reference

Wang, Q. et al. 1997. *Anal. Biochem.* 252, 238.




In this refolding experiment, the iFOLD System 1 revealed that presence of residual detergent is essential for refolding HRV 3C protease. None of the 14 wells with the highest enzymatic activity contain cyclodextrin, which acts as a detergent trap.

## Denaturing Inclusion Bodies, System 2

With the iFOLD Protein Refolding System 2, the purified inclusion bodies are denatured with guanidine-HCl or urea (not included)

and refolded by rapid dilution into the iFOLD System 2 96-well buffer matrix.



### Comparison of GFP-Trx fluorescence after refolding from IB denatured with GuHCl or urea


Thioredoxin-GFP inclusion bodies were isolated and washed with 0.125 M NDSB-201 according to the iFOLD System 2 protocol, denatured with 7.0 M guanidine-HCl or 8.0 M urea, and diluted into the iFOLD System 2 refolding matrix. Following a 24 h incubation at 22°C, refolding reactions were diluted 1:4 with 50 mM Tris-Cl, pH 8.0, and the relative fluorescent intensity recorded (390 nm excitation, 510 nm emission).

## Protein Refolding Temperature

Following addition of the target protein, the refolding plate is incubated at 4–22°C for approximately 20 h. iFOLD™ System 2 buffers maintain a relatively constant pH between 4 and 22°C, facilitating refolding experiments at different temperatures to accommodate the requirements of differ-

ent target proteins. For many proteins, refolding at 22°C results in a greater fraction of active protein than refolding at lower temperatures. If the target protein is thermally unstable or particularly difficult to refold, refolding can be done at lower temperatures (4–10°C)

### Comparison of recombinant enterokinase (rEK) activity after refolding at 10°C and 22°C




Previously prepared rEK IB were denatured with GuHCl according to the iFOLD System 2 protocol. Denatured IB were diluted into an iFOLD Protein Refolding Plate 2 equilibrated at 10°C or 22°C. Following a 20-h incubation, refolding reactions were dialyzed in D-Tube96™ Dialyzers against 2 × 4.0 L buffer (20 mM Tris-Cl, pH 7.5, 150 mM NaCl) overnight at 10°C. Enzymatic activity of refolded and dialyzed reactions was quantified by measuring cleavage of a peptide substrate (Gly-Asp-Asp-Asp-Lys-β-naphthylamide; Sigma G5261).



rEK, a protease with 9 cysteines (8 participating in disulfide bonds), was refolded using iFOLD System 2 at 10°C or 22°C. Successful refolding was determined by a protease activity assay, which indicated that refolding at 10°C was significantly better than at 22°C. Each of the 12 refolding wells with detectable proteolytic activity included glutathione redox buffers.

### Comparison of MMP-12 refolding at 10°C and 22°C



Previously prepared MMP-12 IB were denatured with GuHCl according to the iFOLD System 2 protocol. Denatured IB were diluted into an iFOLD Protein Refolding Plate 2 equilibrated at 10°C or 22°C. Following a 20-h incubation, refolding reactions were dialyzed in D-Tube96™ Dialyzers against 2 × 4.0 L buffer (50 mM Tris-HCl, pH 7.5, 150 mM NaCl, 2.0 mM CaCl<sub>2</sub>, 1 μM ZnCl<sub>2</sub>, and 0.03% Brij®-35) overnight at 10°C. Enzymatic activity of refolded and dialyzed reactions was quantified by measuring cleavage of a BODIPY®-labeled elastin (Invitrogen, E 12056).



In contrast to rEK, MMP-12 refolding was significantly better at 22°C. The refolding screen also revealed that divalent metals (wells C3, C8, and F2) are required for successful MMP-12 refolding.

## Measuring Refolding

Although no universal, accessible, high-throughput method is available to monitor protein refolding, sample absorbance at 340 nm ( $A_{340}$ ) provides an initial screen for refolding efficiency (Vincentelli 2004). Low  $A_{340}$  values indicate soluble proteins. Light scattering instruments can provide a more accurate measure of protein precipitation. However, as a soluble protein does not always correlate with a correctly folded and active protein, we recommend to use an activity assay specific for the target protein.

### Reference

Vincentelli, R. et al. 2004. *Protein Science*. 13, 2782-2792.

## D-Tube 96™ and D-Tube Dialyzers

D-Tube96™ Dialyzers allow convenient, high throughput dialysis of 96 samples simultaneously. The device features the advantages of the D-Tube™ Dialyzer Mini in a 96-tube format. D-Tube Dialyzers are easy to handle tubes with dual membranes providing a large surface area for fast, efficient dialysis. The membrane is ultra clean, EDTA-treated, regenerated cellulose that is sulfur- and heavy metal-free. After screening for optimal refolding conditions with the iFOLD™ Protein Refolding Systems, D-Tube96 Dialyzers provide convenient buffer exchange for the 96 protein samples into a physiologically relevant buffer.

Individual D-Tube Dialyzers are available with molecular weight cut-offs (MWCO) from 3.5 to 14 kDa and are designed with four volume capacities: Mini (10-250 µl), Midi (50-800 µl), Maxi (100-3000 µl), and Mega (3-20 ml). Each kit contains 10 D-Tubes and one floating rack.



Product	Cat. No.	Size	Price
iFOLD™ Protein Refolding System 2	71719-3	1 system	
Components:			
• 30 mL	10X IB•Prep™ Buffer		
• 100 µL	Lysenase™ Bioprocessing Reagent		
• 0.5 mL	1.0 M TCEP		
• 10 mL	1.5 M NDSB-201		
• 10 mL	iFOLD System 2 Denaturation Buffer (50 mM Tris-HCl, 0.2 M NaCl, 2.0 mM EDTA, 7.0 M GuHCl, pH 8.0)		
• 1	iFOLD Protein Refolding Plate 2		
• 2	Aluminum plate sealers		
iFOLD™ Protein Refolding System 1	71552-3	1 system	
Components:			
• 30 ml	10X IB-Prep™ Buffer		
• 0.5 ml	1M TCEP		
• 1.5 ml	Triton® X-100		
• 0.1 ml	Lysenase™ Bioprocessing Reagent		
• 10 ml	30% N-Lauroylsarcosine		
• 50 ml	50X iFOLD Dialysis Buffer		
• 1	iFOLD Protein Refolding Plate 1		
• 2	Aluminum Plate Sealers		

### Features

- Efficient dialysis – dual membrane surface
- D-Tube96 Dialyzers are ideal for sample volume from 10–250 µl per well/tube
- D-Tube Dialyzers are ideal for sample volume from 10 µl to 20 ml per well/tube
- Convenient – use single or multichannel pipet or robotic dispensing
- High sample recovery, >97%
- Protease-, RNase-, and DNase-free
- Versatile – dialyze proteins, oligonucleotides, RNA, or DNA

Product	Cat. No.	Size	Price
D-Tube™ Dialyzer Mini, MWCO 6-8 kDa	71504-3	1 kit	
D-Tube™ Dialyzer Mini, MWCO 12-14 kDa	71505-3	1 kit	
D-Tube™ Dialyzer Midi, MWCO 3.5 kDa	71506-3	1 kit	
D-Tube™ Dialyzer Midi, MWCO 6-8 kDa	71507-3	1 kit	
D-Tube™ Dialyzer Maxi, MWCO 3.5 kDa	71508-3	1 kit	
D-Tube™ Dialyzer Maxi, MWCO 6-8 kDa	71509-3	1 kit	
D-Tube™ Dialyzer Maxi, MWCO 12-14 kDa	71510-3	1 kit	
D-Tube™ Dialyzer Mega, 10 ml, MWCO 3.5 kDa	71739-3	1 kit	
D-Tube™ Dialyzer Mega, 10 ml, MWCO 6-8 kDa	71740-3	1 kit	
D-Tube™ Dialyzer Mega, 10 ml, MWCO 12-14 kDa	71741-3	1 kit	
D-Tube™ Dialyzer Mega, 15 ml, MWCO 3.5 kDa	71742-3	1 kit	
D-Tube™ Dialyzer Mega, 15 ml, MWCO 6-8 kDa	71743-3	1 kit	
D-Tube™ Dialyzer Mega, 15 ml, MWCO 12-14 kDa	71744-3	1 kit	
D-Tube™ Dialyzer Mega, 20 ml, MWCO 3.5 kDa	71745-3	1 kit	
D-Tube™ Dialyzer Mega, 20 ml, MWCO 6-8 kDa	71746-3	1 kit	
D-Tube™ Dialyzer Mega, 20 ml, MWCO 12-14 kDa	71747-3	1 kit	
D-Tube96™ Dialyzer, MWCO 6-8 kDa	71712-3	1 kit	
D-Tube96™ Dialyzer, MWCO 12-14 kDa	71713-3	1 kit	

**Argentina**

Merck Quimica Argentina S.A.I.C.  
Tel: +54 11 4546 8100  
+54 11 4546 8156  
Fax: +54 11 4546 8175  
E-mail: gbarros@merck.com.ar  
www.merck.com.ar

**Australia**

Merck Pty. Limited  
Tel: +61 3 9728 7600  
Fax: +61 3 9728 1351  
E-mail: merck@merck.com.au  
www.merck.com.au

**Brazil**

Merck S.A.  
Toll Free: 0800 21 9292  
Tel: +55 11 3346 8500  
Fax: +55 11 3207 5040  
E-mail: quimica@merck.com.br  
www.merck.com.br

**Imprint do Brasil Ltda.**

Tel: +55 19 3772 2900  
Fax: +55 19 3273 5389  
E-mail: imprint@imprint.com.br  
www.imprint-corp.com

**Caribbean**

Merck Centroamericana S.A.  
See information under Central America

**Central America**

Merck, S.A. Guatemala  
Tel: +50 2 2410-2300  
Fax: +50 2 2434-2954  
E-mail: quimicos@merck.com.gt  
www.merck.com.gt

**Chile**

Merck S.A.  
Tel: +56 2 3400 000  
Fax: +56 2 3400 199  
E-mail: mqch@merck.cl  
www.merck.cl

**China**

Merck China  
Toll Free: 800 820 8872  
Tel: +86 21 3222 4788  
Fax: +86 21 6247 9680  
E-mail: Bioteam@merck-china.com  
www.merckbio.cn

**Colombia**

Merck Colombia S.A.  
Tel: +57 1 425 4770  
Fax: +57 1 425 5407  
E-mail: mcsa@merck.com.co  
www.merck.com.co

**Ecuador**

Merck C.A.  
Tel: +593 2 2981677  
Fax: +593 2 2981644  
E-mail: sicmerck@merck.com.ec  
www.merck.com.ec

**Guatemala**

Merck Centroamericana S.A.  
See information under Central America

**Hong Kong**

Onwon Trading Limited  
Tel: +852 2757 7569  
Fax: +852 2757 7211  
E-mail: info@onwon.com.hk  
www.onwon.com.hk

**India**

Merck Specialities Private Limited  
Tel: +91 22 6660 9184  
Tel: +91 22 6660 9000  
Fax: +91 22 2495 4590  
Fax: +91 22 2495 0307  
E-mail: life.science@merck.co.in  
www.merck.co.in

**Indonesia**

Pt. Merck Tbk.  
Toll Free: 0800 140 1253  
Tel: +62 21 841 3889  
Fax: +62 21 841 5537  
E-mail: chemicals@merck.co.id

**Israel**

Mercury Scientific & Industrial  
Products Ltd.  
Tel: +972 3 9387164  
Fax: +972 3 9387174  
E-mail: mercury@mercury-ltd.co.il  
www.mercury-ltd.co.il

**Japan**

Merck Ltd.  
Tel: +81 0120 189 390  
Fax: +81 0120 189 350  
E-mail: service@merck.co.jp  
www.merck.co.jp

**Korea**

Merck Limited  
Tel: +82 2 2185 3836  
Fax: +82 2 2185 3870  
E-mail: service@merck.co.kr  
www.merck.co.kr

**Malaysia**

Merck Sdn Bhd  
Tel: +6 03 7882 4888  
Tel: +6 03 7880 0811  
Fax: +6 03 7880 0749  
Fax: +6 03 7880 0792  
E-mail: chemlab@merck-de.com.my  
www.merck.com.my

**New Zealand**

Merck Limited  
Toll Free: 0800 46 37 25  
Tel: +64 06 356 7328  
Fax: +64 06 356 7311  
E-mail: info@merck.co.nz  
www.merck.co.nz

**Mexico**

Control Técnico Y Representaciones  
Monterrey, Nuevo, León  
Tel: +52 81 8158 0600  
Fax: +52 81 8373 2891  
E-mail: ctrscientific@infosel.net.mx

**Mexico City**

Tel: +52 55 5208 5197  
+52 55 5208 5198  
+52 55 5208 8116  
Fax: +52 55 5203 6229

**Pakistan**

Merck Marker (Pvt) Ltd.  
Tel: +92 21 455 9210  
Fax: +92 21 453 5294  
E-mail: lab@merck.com.pk  
www.merck.com.pk

**Peru**

Merck Peruana S.A.  
Tel: +51 1 6187 500  
Fax: +51 1 4372 955  
E-mail: merck.peruana@merck.com.pe  
www.merck.com.pe

**Philippines**

Merck Inc.  
Tel: +63 2 815 4067  
Tel: +63 2 814 5221  
Fax: +63 2 815 4883  
E-mail: tish.aligada@merck.ph

**Singapore**

Merck Pte. Ltd.  
Customer Hotline: +65 6890 6660  
Tel: +65 6890 6638  
Fax: +65 6890 6636  
E-mail: chem@merck-de.com.sg  
www.merck.com.sg

**Taiwan**

Merck Ltd.  
Tel: +886 2 2742 2788  
Fax: +886 2 2742 2766  
E-mail: bioservice@merck.com.tw  
www.merckbio.com.tw

**Thailand**

Merck Ltd.  
Tel: +66 2 667 8333  
Fax: +66 2 667 8338  
E-mail: customercare@merck.co.th  
www.merck.co.th

**Venezuela**

Merck S.A.  
Tel: +58 21 2235 1379  
Fax: +58 21 2237 7632  
E-mail: mven@merck.com.ve  
www.merck.com.ve

**Vietnam**

Merck Representative Office  
Tel: +84 8 932 0187  
Fax: +84 8 526 0201  
E-mail: merckvnadmin@hcm.vnn.vn

220023-2007APLA  
iFOLD Protein Refolding Systems