

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

pFLAG-CMV™-3 Expression Vector

Catalog Number **E6783** Storage Temperature –20 °C

Product Description

The pFLAG-CMV-3 Expression Vector is a 6.3 kb vector used for transient or stable expression in mammalian cells. The vector is a derivative of pCMV5¹ used for transient or stable expression and secretion of a properly inserted open reading frame as an N-terminal FLAG[®] fusion protein.

The promoter-regulatory region of the human cytomegalovirus² drives transcription of FLAG-fusion constructs. The preprotrypsin leader sequence³ precedes the FLAG sequence and directs secretion of the fusion protein into the culture medium. The aminoglycoside phosphotransferase II gene⁴ (neo^f) confers resistance to aminoglycosides such as G 418,⁵ allowing for selection of stable transfectants.

The pFLAG-CMV-3 Expression Vector is a shuttle vector containing both bacterial and SV40 origins of replication for propagation in both *Escherichia coli* and mammalian cells. Efficiency of replication and genomic integration is optimal when using host cells that express the SV40 large T antigen (e.g. COS-7). pFLAG-CMV-3 has been used for stable transfection of HEK 293 cells⁶.

The FLAG epitope is a small, hydrophilic 8 amino acid tag (DYKDDDDK)⁷ that provides for sensitive detection and high quality purification using ANTI-FLAG[®] products (visit www.sigma-aldrich.com for a complete listing). Removal of the N-terminal FLAG tag is possible using enterokinase, which cleaves following the Asp-Asp-Asp-Lys recognition site at the C-terminal end of the FLAG peptide.

The pFLAG-CMV-3-BAP Control Plasmid is a 7.7 kb derivative of pCMV5¹ used for transient expression and secretion of N-terminal FLAG bacterial alkaline phosphatase fusion protein in mammalian cells.

The promoter-regulatory region of the human cytomegalovirus² drives transcription of bacterial alkaline phosphatase. The preprotrypsin leader sequence³ precedes the FLAG sequence. The aminoglycoside phosphotransferase II gene⁴ (Neo) confers resistance to aminoglycosides such as G 418.⁵

pFLAG-CMV-3-BAP Control Plasmid is a shuttle vector for *E. coli* and mammalian cells. Efficiency of replication and genomic integration is optimal when using an SV40 T antigen-expressing host, such as COS cells.

Map positions of key features in the pFLAG-CMV-3 Expression Vector and the pFLAG-CMV-3-BAP Control Plasmid can be found at www.sigma.com/vectormaps

Components

- pFLAG-CMV-3 Expression Vector 20 μg Catalog Number E8770
 Supplied as 0.5 mg/ml in 10 mM Tris-HCl, pH 8.0, 1 mM EDTA.
- pFLAG-CMV-3-BAP Control Plasmid 20 μg Catalog Number C3972 Supplied as 0.5 mg/ml in 10 mM Tris-HCl, pH 8.0, 1 mM EDTA.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Store at -20 °C

References

- Andersson, S., et al., Cloning, structure, and expression of the mitochondrial cytochrome P-450 sterol 26-hydroxylase, a bile acid biosynthetic enzyme. J. Biol. Chem., 264, 8222-8229 (1989).
- Thomsen, D. R., et al., Cloning, structure, and expression of the mitochondrial cytochrome P-450 sterol 26-hydroxylase, a bile acid biosynthetic enzyme. Proc. Natl. Acad. Sci. USA, 81, 659-663 (1984).
- 3. Stevenson, B. J., et al., Sequence organization and transcriptional regulation of the mouse elastase II and trypsin genes. *Nucleic Acids Res.*, **21**, 8307-8330 (1986).

- 4. Brewer, C. B., Cytomegalovirus plasmid vectors for permanent lines of polarized epithelial cells. *Methods in Cell Biology*, **43**, 233-245 (1994).
- Jimenez, A. and Davies, J., Expression of a transposable antibiotic resistance element in Saccharomyces. *Nature*, 287, 869-871 (1980).
- 6. Soto, H., et al., The CC chemokine 6Ckine binds the CXC chemokine receptor CXCR3. *Proc. Natl. Acad. Sci. USA*, **95**, 8205-8210 (1998).
- 7. Hopp, T. P., et al., A short polypeptide marker sequence useful for recombinant protein identification and purification., *Bio/Technology*, **6**, 1204-1210 (1988).

FLAG is a registered trademark and pFLAG-CMV is a trademark of Sigma-Aldrich Biotechnology LP and Sigma-Aldrich Co..

AH,RS,PHC 08/10-1