

30020 D-Cycloserine (D-4-Amino-3-isoxazolidone)

CAS number: 68-41-7

Product Description:

prepared from microbial source

Appearance: White powder

Molecular Formula: C₃H₆N₂O₂

Molecular Weight: 102.1 g/mol


Mp: ~ 150°C (dec.)

E_{226 nm}: 402 (1% in water)

[α]_D²³: +115° (c = 1.0%, water)¹

[α]_D²⁰: +114 ± 3° (c = 2% in water)¹⁹

[α]₅₄₆²⁰: +138 ± 3° (c = 2% in water)¹⁹



D-Cycloserine, a structural analog of D-alanine, is a broad spectrum antibiotic produced by certain strains of *Streptomyces orchidaceus* or *S. garphalus*.¹⁻⁵

D-cycloserine (at 100 to 200 µg/ml) inhibits the synthesis of bacterial cell walls (involving peptidoglycan synthesis) by preventing formation of D-alanine from L-alanine and hence the formation of peptide bonds involving D-alanine.⁴ D-Cycloserine has antibiotic activity *in vitro* against growth phase Gram-negative bacteria including *Escherichia coli* (working concentration of approx. 200 µg/ml)⁴, strains of *Staphylococcus aureus*, *Nocardia* species and *Chlamydia*,³ and some mycobacteria including *Mycobacterium tuberculosis*. The minimum inhibitory concentrations (MIC) *in vitro* for *M. tuberculosis* range from 5-20 µg/ml. Studies *in vitro* show no suppression of growth in cultures made in media containing D-alanine which appears to block the antibacterial action of D-cycloserine.³

D-Cycloserine is an excitatory amino acid and partial agonist at the glycine binding site of the N-methyl-D-aspartate (NMDA) receptor.⁶⁻⁸ At low doses it is a cognitive enhancer that improves learning and memory in several experimental models of disease and cognitive deficit.^{6,7,9-14} At high doses, D-cycloserine is an anti-convulsant.^{15,16} Intermediate doses potentiate the anti-convulsant action of phenytoin but block its long-term memory impairment.¹⁶

The HPLC determination of D-Cycloserine in plasma and urine¹⁷ and a colorimetric method for determination of Cycloserine in biological fluids¹⁸ have been reported. UV, IR, NMR and mass spectra and pharmacokinetics of D-Cycloserine have been reported.²

Applications:

This product is used for preparing diverse selective media for molecular biology and microbiology²¹. Cycloserine is a broad-spectrum antibiotic used to treat tuberculosis. It is used rarely for treating noninfectious diseases.²²

Preparation Instructions:

D-cycloserine is soluble in deionized water up to 100 mg/ml. A solution of 50 mg/ml cycloserine in water is clear and colorless or very faintly yellow. D-cycloserine is also soluble at 1 in 50 parts of 96% ethanol but practically insoluble in chloroform and ether. It is also slightly soluble in methanol or propylene glycol. Stock solutions (e.g. 10 mg/ml) of D-cycloserine may also be prepared immediately before use in 0.1 M sodium phosphate buffer, pH 8.0.



Storage/Stability:

Storage Temperature of the powder is -18°C

It is generally recommended to prepare solutions immediately before use because neutral or acidic solutions are unstable.⁴ However aqueous solutions buffered to pH 10 with sodium carbonate may be stored for up to one week if stored at 2 °C to 8 °C. In addition, aqueous solutions of D-Cycloserine have been stored in aliquots at -20 °C and thawed just prior to use.¹²

References:

1. Stammer, C.H., et al., *J. Am. Chem. Soc.*, 77, 2346 (1955).
2. El-Obeid, H.A. and Al-Badr, A.A, *Analytical Profile of D-Cycloserine in Analytical Profiles of Drug Substances*, Vol. 18, p. 567 (Academic Press, New York, 1989).
3. Goodman and Gilman's *The Pharmacological Basis of Therapeutics*, 9th Ed., Hardman, J.G. et al., (Eds.), p. 1164 (McGraw-Hill, New York, NY, 1995).
4. Raleigh, E.A., et al., *Selected Topics from Classical Bacterial Genetics in Short Protocols in Molecular Biology*, 4th Ed., Ausubel, F.M. et al., (Eds.), Unit 1.4, p. 1-9 (John Wiley & Sons, Inc., New York, 1999).
5. Kuehl, F.A., Jr. et al., *J. Am. Chem. Soc.*, 77, 2344 (1955).
6. Sirvio, J., et al., *Neurosci. Lett.*, 146, 215 (1992).
7. Ohno, M. and Watanabe, S., *Eur. J. Pharmacol.*, 318, 267 (1996).
8. Watson, et al., D-Cycloserine acts as a partial agonist at the glycine modulatory site of the NMDA receptor expressed in *Xenopus oocytes* *Brain Res.* 510, 158-160 (1990)
9. Flood, J.F., et al., *Eur. J. Pharmacol.*, 221, 249 (1992).
10. Schuster, G.M. and Schmidt, W.J., *Eur. J. Pharmacol.*, 224, 97 (1992).
11. Temple, M.D. and Hamm, R.J., *Brain Res.*, 741, 246 (1996).
12. Nakazato, E., et al., Cholinergic and glutamatergic activation reverses working memory failure by hippocampal histamine H1 receptor blockade in rats *Life Sci.* 67, 1139-1147 (2000)
13. Schneider, J.S., et al., Effects of the partial glycine agonist D-cycloserine on cognitive functioning in chronic low dose MPTP-treated monkeys *Brain Res.* 860, 190-194 (2000)
14. Pussinen, R., Sirvio, J., Effects of D-cycloserine, a positive modulator of N-methyl-D-aspartate receptors, and ST 587, a putative alpha-1 adrenergic agonist, individually and in combination, on the non-delayed and delayed foraging behaviour of rats assessed in the radial arm maze *Org. Process R & D* 13, 171 (1999)
15. Chessell, I.P., et al., *Brain Res.*, 565, 345 (1991).
16. Wlaz, P., et al., Influence of D-cycloserine on the anticonvulsant activity of phenytoin and carbamazepine against electroconvulsions in mice *Epilepsia* 37, 610-617 (1996)
17. Musson, D.G. et al., *J. Chromatog.*, 414, 121 (1987).
18. Jones, L.R., *Anal. Chem.*, 28, 39 (1956).
19. Sigma-Aldrich quality control
20. M.K. Jain, *Handbook of Enzyme Inhibitors (1965-1977)* New York 112 (1982)
21. Orth, D.S.: Comparison of sulfite-polymyxin-sulfadiazine medium and tryptose-sulfite-cycloserine medium without egg yolk for recovering *Clostridium perfringens*, *Appl. Envir. Microbiol.*, 33, 986-988 (1977).
22. Hershfield E., *Tuberculosis: 9. Treatment*, Canadian Medical Association, 161(4):405-11. (1999)

Precautions and Disclaimer:

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

The vibrant M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.

Detailed information on trademarks is available via publicly accessible resources.

© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

