



HybriScan®D Beer

Rapid test system for the detection of beer-spoiling bacteria

Product-No.: 62533







Contact information:

HybriScan® - Rapid Test System (R&D)

Dr. Helmut Maucher

Phone: (+49) - 3494 - 6364 15 e-mail: contact@scanbec.de

Sales Organisations

Argentina SIĞMA-ALDRICH DE ARGENTINA S.A.

Free Tel: 0810 888 7446 Tel: (+54) 11 4556 1472 Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD. Free Tel: 1800 800 097 Free Fax: 1800 800 096 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH Tel: (+43) 1 605 81 10 Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA. Free Tel: 0800 14747 Free Fax: 0800 14745 Tel: (+32) 3 899 13 01 Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA. Free Tel: 0800 701 7425 Tel: (+55) 11 3732 3100 Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD. Free Tel: 1800 565 1400 Free Fax: 1800 265 3858 Tel: (+1) 905 829 9500 Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI) TRADING CO. LTD. Free Tel: 800 819 3336 Tel: (+86) 21 6141 5566 Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH S.R.O. Tel: (+420) 246 003 200 Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S

Tel: (+45) 43 56 59 10 Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY Tel: (+358) 9 350 9250 Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.I. Free Tel: 0800 211 408 Free Fax: 0800 031 052 Tel: (+33) 474 82 28 00 Fax: (+33) 474 95 68 08

Germany SIGMA-ALDRICH CHEMIE GmbH Free Tel: 0800 51 55 000 Free Fax: 0800 64 90 000 Tel: (+49) 89 6513 0 Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD. Tel: (+30) 210 994 8010 Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft Ingyenes zöld telefon: 06 80 355

Ingyenes zöld fax: 06 80 344 344 Tel: (+36) 1 235 9055 Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS PRIVATE LIMITED

Telephone

Bangalore: (+91) 80 6621 9600 New Delhi: (+91) 11 4165 4255 Mumbai: (+91) 22 2570 2364 Hyderabad: (+91) 40 6684 5488

Fax

Bangalore: (+91) 80 6621 9650 New Delhi: (+91) 11 4165 4266 Mumbai: (+91) 22 2579 7589 Hyderabad: (+91) 40 6684 5466

Ireland

SIGMA-ALDRICH IRELAND

LTD

Free Tel: 1800 200 888 Free Fax: 1800 600 222 Tel: (+353) 1 404 1900 Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD. Free Tel: 1 800 70 2222 Tel: (+972) 8 948 4100 Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.I. Numero Verde: 800 827018 Tel: (+39) 02 3341 7310 Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K. Tokyo Tel: (+81) 3 5796 7300 Tokyo Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA Free Tel: (+82) 80 023 7111 Free Fax: (+82) 80 023 8111 Tel: (+82) 31 329 9000 Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD Tel: (+60) 3 5635 3321 Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA. SAdeCV Free Tel: 01 800 007 5300 Free Fax: 01 800 712 9920 Tel: 52 722 276 1600 Fax: 52 722 276 1601

The Netherlands SIGMA-ALDRICH CHEMIE BV Free Tel: 0800 022 9088 Free Fax: 0800 022 9089 Tel: (+31) 78 620 5411 Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD. Free Tel: 0800 936 666 Free Fax: 0800 937 777 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS Tel: (+47) 23 17 60 60 Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o. Tel: (+48) 61 829 01 00 Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, Free Tel: 800 202 180

Free Fax: 800 202 178 Tel: (+351) 21 924 2555 Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC Tel: +7 (495) 621 6037 Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD. Tel: (+65) 6779 1200 Fax: (+65) 6779 1822

South Africa

SIGMA-ALDRICH SOUTH AFRICA (PTY) LTD. Free Tel: 0800 1100 75 Free Fax: 0800 1100 79 Tel: (+27) 11 979 1188 Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA. S.A

Free Tel: 900 101 376 Free Fax: 900 102 028 Tel: (+34) 91 661 99 77 Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB Tel: (+46) 8 742 4200 Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH Free Tel: 0800 80 00 80 Free Fax: 0800 80 00 81 Tel: (+41) 81 755 2828 Fax: (+41) 81 755 2815

United Kingdom SIGMA-ALDRICH COMPANY LTD Free Tel: 0800 717 181 Free Fax: 0800 378 785 Tel: (+44) 1747 833 000 Fax: (+44) 1747 833 313 SAFC (UK) Free Tel: 0800 71 71

United States SIGMA-ALDRICH P.O. Box 14508 St. Louis, Missouri 63178 Toll-Free: 800 325 3010 Toll-Free Fax: 800 325 5052 Call Collect: (+1) 314 771 5750 Tel: (+1) 314 771 5765 Fax: (+1) 314 771 5757

Internet

sigma-aldrich.com

Technical Service: flukatec@sial.com

Product Specifications

Cat. No.:

Number of tests: 96 tests, incl. standard series

4 - 8°C, 12 month

Storage: Test duration: approx. 2-2.5 hours (after pre-enrichment)

Pre-enrichment time: 24 - 30 hours

Sensitivity: 1 -10 CFU/L (after pre-enrichment)

Specificity: all beer spoiling bacteria of the genera Lactobacillus, Pediococcus, Pectinatus and Megasphaera





HybriScan®D Beer-Test Protocol

Working Principle

HybriScan® D Beer is a rapid molecular test system for the detection of all beer-spoiling microorganisms in only one test. HybriScan® D Beer is based on the detection of target molecules from the microorganism of interest by means of specific capture and detection probes in a so-called sandwich hybridization. The target molecules of the beer contaminants contained in the sample are captured in a specific microtiter binding plate. All other unbound sample components are removed by several washing steps, so that only beer-spoiling microbes are detected in a highly specific way. In addition to the capture probe, a detection probe is coupled to the target molecule. An enzyme is attached afterwards in a subsequent incubation step. After several washing steps, reaction with a colour substrate gives a blue colouration, which changes to yellow after the addition of a stop solution. The yellow colour enables highly sensitive photometric measurement at 450 nm. Comparison is made with the standard solutions contained in the test kit.

Technical Notes

After starting the test procedure, perform each of the following steps without interruption and within the given time limits:

For each sample use an individual single-use pipette tip to avoid cross-contamination.

Close bottles immediately after use and store them at the temperatures specified on the labels. Do not interchange caps and bottles.

Samples and standards should be tested together for more accurate results.

Do not mix or replace components from test kits of different charges.

Incubation at room temperature refers to a laboratory temperature of 20 to 25°C.

Do not use the test kit after the expiration date listed on the package.

Safety

All reagents contained in the test kit are for in vitro use only.

Test solution D contains formamide. Avoid contact with eyes, skin and the respiratory system. In event of contact with eyes or skin, rinse immediately with plenty of water. If the reagent is inhaled, immediately remove the individual to fresh air and seek medical attention.

Stop solution H contains 1 N sulfuric acid. Avoid contact with eyes and skin. In the event of contact with eyes and skin rinse immediately with plenty of water.

Handling of the kit components and disposal of waste should be performed according to standard laboratory safety guidelines.





Reagents and Storage Conditions

The reagents contained in the test kit are sufficient for at least 96 tests, including 6 standard series. The kit components should be stored between +2 to $+8^{\circ}$ C as indicated on the labels. Do not freeze the test kit components!

Kit components:

1.	Microwell plate, ready to use, 96 wells	1
2.	Binding plate, ready to use, 96 wells	1
3.	Standards 1 – 4 ^{a)} (white screw caps); blank and different concentrations of synthetic RNA as positive control and for preparation of a calibration curve for quantitative measurements	0.2 mL each
4.	Lysis Reagent A (red screw cap), ready to use	1.2 mL
5.	Lysis Buffer B a) (red cap), ready to use	4.5 mL
6.	Lysis Buffer C a) (red cap), ready to use	5.5 mL
7.	Test Solution D (yellow cap), ready to use	4,5 mL
8.	Washing Solution E b) (blue cap), ready to use	90 mL
9.	Enzyme Solution F (green screw cap), dilute a suitable amount 1:100 with Washing Solution E before use	0.120 mL
10. Substrate Solution G b) (green cap), ready to use		10 mL
11. Stop Solution H (green cap) 1 N sulfuric acid, ready to use		5 mL
12. Glass beads (colourless cap), sterile, ready to use		4 mL

^{a)}Components contain SDS, which precipitates at lower temperatures. Equilibrate to room temperature before use.

Additional equipment and materials (required, not supplied with kit)

- Centrifuge for microreaction tubes (1.5 and 2 mL)
- Thermoshaker for microreaction tubes and microwell plates
- Vacuum filtration unit
- 3 Pipettes (2–20 μ L, 20–200 μ L, 200-1000 μ L) with corresponding tips; optional 8-channel pipette (20-200 μ L)
- Microwell plate-photometer
- · Enrichment medium, incubator
- Microreaction tubes (2 mL), cultivation tubes (12 mL), reagent-reservoirs, membrane filter discs (0.45 μm)

b) Equilibrate to room temperature before use.





Specificity:

Following species of beer-spoiling microorganisms can be detected with the HybriScan® D Beer-Kit:

• **Genus Lactobacillus:** Lactobacillus acidophilus

Lactobacillus brevis
Lactobacillus brevisimilis
Lactobacillus buchneri
Lactobacillus casei
Lactobacillus collinoides
Lactobacillus coryniformis
Lactobacillus curvatus
Lactobacillus fermentum
Lactobacillus fructivorans

Lactobacillus malefermentans

Lactobacillus lindneri

Lactobacillus parabuchneri (frigidus)

Lactobacillus paracasei Lactobacillus paraplantarum Lactobacillus plantarum Lactobacillus rhamnosus

• **Genus Pediococcus:** Pediococcus acidilactici

Pediococcus claussenii Pediococcus damnosus Pediococcus inopinatus Pediococcus parvulus Pediococcus pentosaceus

• **Genus Pectinatus:** Pectinatus cerevisiiphilus

Pectinatus frisingensis

Genus Megasphaera: Megasphaera cerevisiae





Test protocol

(1) Sample preparation

Transfer a 2 mL aliquot from the pre-cultivation tube via pipette to a 2 mL microreaction tube that contains a spatula-tip amount of glass beads. Centrifuge the samples for 2 minutes at maximum speed of 13,000 rpm. Remove the supernatant carefully with a pipette.

Note:

Avoid strong shaking after centrifugation to avoid resuspending the bacteria pellet. Centrifuge a second time if necessary. Samples with high yeast content should be centrifuged for 1 minute at 1,000 rpm to allow the yeast to settle. Take 1 mL of the supernatant and proceed with step 1. To increase the sensitivity, start with 2 x 2 mL sample and combine the supernatants.

(2) Cell lysis

Add 40 μ L of **Lysis Buffer B** (bottle with red cap) to the cell pellet and 10 μ L of **Lysis Reagent A*** (microreaction tube with red screw cap), mix well and incubate for 15 minutes at 37°C in a thermoshaker. Next, add 50 μ L of **Lysis Buffer C** (bottle with red cap). Incubate again for 15 minutes at 37°C with shaking at 1,400 rpm in the thermoshaker. Centrifuge the samples for 10 minutes at 13,000 rpm using a microcentrifuge. Use 10 μ L of this supernatant in protocol step 3 (hybridization).

*Note: In the case of a large number of samples prepare a Master Mix of Lysis Reagent A and Lysis Buffer B before use. Pipette 50 µL of the Master Mix to each cell pellet.

Preparation for subsequent steps:

Change the top of the thermoshaker and fix the manifold for microwell plates. Set the temperature to 50°C and shaking speed to 500 rpm. For quantitative analysis of beer-spoiling microorganisms we recommend a repeat determination of the 4 standard solutions (microreaction tubes with white screw cap).

To each well of standards (repeat determinations, e.g. A1-H1) and samples add 45 μ L of Test Solution D (bottle with yellow cap). Incubate the plate at 50°C for 5 minutes in the thermoshaker.

(3) Hybridization

Add 10 μ L of **Standard 1** into the wells A1 and B1; 10 μ L of **Standard 2** into C1 and D1; 10 μ L of **Standard 3** into E1 and F1 and 10 μ L of **Standard 4** into G1 and H1. Apply 10 μ L of each sample (supernatant from step 2) into the respective well position. Cover the microwell plate with a lid and incubate it in the thermoshaker for 10 minutes at 50°C and 500 rpm.

Note:

When adding the standards and samples, to avoid cooling do not remove the microwell plate from the thermoshaker.

The supernatant from step 2 can be stored at -20°C for future use.





(4) Coupling to the binding plate

Transfer 50 μ L of the reaction mixes from each well to the corresponding wells of the binding plate and shake for 10 minutes at 50°C and 500 rpm in the thermoshaker.

Note:

Unused stripes of the plate should be stored in the sealed original packing at 4 to 8 °C.

Preparation for subsequent steps:

The **Enzyme Solution F-Washing Solution E** 1:100 dilution must be prepared immediately before use. It cannot be stored. Prepare only the amount needed for the test, e. g. for 16 reactions combine 1700 μ L **Washing Solution E** and 17 μ L **Enzyme Solution F (microreaction tube with green screw cap)**.

Note

Briefly spin down enzyme solution F prior use to collect the liquid at the bottom of the tube.

(5) Enzymatic reaction

Discard the liquid from each well by inverting and gently tapping of the plate on an absorbent layer. Set the temperature to 25°C. Add 200 μ L **Washing Solution E** (bottle with blue cap) and incubate for 2 minutes at room temperature. Discard the liquid. Pipette 100 μ L of the diluted **Enzyme solution F**, prepared as described above "preparation for subsequent steps", to each well. Cover the binding plate with a lid and incubate it in the thermoshaker for 10 minutes at 25°C and 500 rpm.

(6) Washing

Discard the liquid from each well. Add 200 μ L of **Washing Solution E** (bottle with blue cap) to each well and incubate the microplate (with lid) for 1 minute at 25°C and 500 rpm in the thermoshaker. Repeat washing each well once.

Preparation for subsequent steps:

Switch on computer and the microplate reader.

(7) Substrate Reaction

After discarding the Washing Solution from the second wash step, add 100 μ L of **Substrate Solution G** (bottle with green cap) to each well. Cover the microplate with a lid and incubate it in a thermoshaker for 10 minutes at 25°C and 500 rpm. Stop the reaction by adding 50 μ L of **Stop Solution H** (bottle with green cap) to each well. The addition of acid creates a yellow colour change. Mix shortly (10 sec, 500 rpm) in the thermoshaker and remove air bubbles, if present.

Note:

For qualitative analysis results can be measured by visual inspection. Compared to the blanks (A1, B1), which should be colourless, a blue colour change indicates contamination of the sample.

Data analysis using the HybriScan®-Software

(8) Signal read-out and data analysis using the HybriScan®-Software

Start the reader and open the HybriScan[™]-Software. Insert the microwell plate into the reader, with position A1 rear left. Click button **Start Measurement** and choose the appropriate test. Click **Start**. The instrument measures the absorbance of any well. Results are saved automatically in the database. Alternatively, results of tests saved as text file can be opened by choosing option **File/Open** in the main menu.

The graphical layout of the results corresponds to the position of the sample in the microwell plate. By clicking on each position, you can enter the name of the sample. In the lower pane of the graphical





user interface the user name and other comments can be inserted. For each position you can hide or unhide parameters by choosing the button **View.**

Click the button **Report**. In the upper pane of the window the regression line calculated from the four standards is displayed. Values of standards are represented by red dots. They should be close to the regression curve (blue line). Your results can be printed by choosing the option **Print Report**. For post processing (like in Microsoft Excel) results can be exported by choosing the option **Export Report**.

(9) Data interpretation

Based on the absorbance measurement, HybriScan®-Software will help to evaluate the analyzed samples. The quantitative analysis is performed automatically on the basis of the standard values (A1 to H1). For each well the corresponding number of cells are calculated and indicated by the software numerically and with a colour code. The change of green to red enables a fast visual recognition of contaminated samples.

The cell number calculated by the HybriScan $^{\$}$ -Software reflects the amount of cell equivalents in the 10 μ L sample introduced to the test solution. If you started with 2mL culture, you have to multiply the calculated cell number by 5 to get the amount of cells per mL.

Furthermore, the mean values of optical density of the negative control (standard 1) are used to evaluate, if a sample is positive, questionable or negative.

The HybriScan®-Software enables both, graphical and tabular representation of the measurement data of the respective samples. Optical density and calculated numbers of cells are displayed by the software tool.

Data analysis without the HybriScan®-Software

(10) Signal read-out using VIS-photometer

Start the reader and open the photometer Software. Insert the microwell plate into the reader, with position A1 rear left. Start the measurement. The instrument measures the absorbance of any position at 450 nm.

(11) Internal control

We recommend using all four standards as internal controls for the test procedure. If the analysis is performed correctly your measured data will be near the regression line. We recommend this method especially for inexperienced user.

For linear regression set standard1 = 0, standard2 = 4, standard3 = 8 and standard4 = 16 at the abscissa (x-coordinate). The straight line must subtend the y-axis at S1 = 0.

Regression line offers only verification of analysis-quality; a quantification is not possible.

For experienced user it is sufficient to use only standard S1 and S3. These standards are necessary for a qualitative data analysis.





(12) Qualitative analysis

For the measurement to be valid, the quotient of the mean value of the positive control (S3) divided by the mean value of the negative control (S1) must be greater than 4.0.

Evaluation of the samples is performed using the following formula:

Sample OD% =
$$\frac{\text{OD }_{\text{Sample}}\text{- MV OD}_{\text{NC}}}{\text{MV OD}_{\text{PC}}\text{- MV OD}_{\text{NC}}} \times 72.100\%$$

MV mean value PC positive control (S3) NC negative control (S1)

Sample OD% values are used to evaluate the sample status:

Samples with OD% values under 10 are considered negative.

Samples with OD% values from 10 to < 20 are considered questionable.

Samples with OD% values \geq 20 are considered positive.

(13) Quantitative analysis

For quantitative analysis (semi - quantitative analysis) use following formula to calculate cell count of your samples:

$$Cell count_{Sample} = \frac{OD_{Sample} \times Cell count_{Standard}}{MV OD_{Standard}}$$

MV mean value

For calculating the cell count, choose a standard, which has a OD next to sample - OD.

Cell count - standards:

S1 = 0

S2 = 4000

S3 = 8000

S4 = 16000

Note:

Quantitative analysis is only possible, if sensitivity was reached (see product specifications), your measured data are near the regression line and <u>no</u> pre-enrichment was performed!





Short Protocol

- 1. Filter 100–1,000 mL of sample (vacuum filtration unit; membrane filter disc; 0.45 µm pore size)
- 2. After filtration incubate the filter disc in 5 mL of enrichment medium (e.g. NBB, 28°C, 24-30 hours)
- 3. Remove 2 mL of sample from the enrichment medium, add glass beads, centrifuge (13,000 rpm, 2 minutes) and discard the supernatant
- 4. Add 40 μ L of **Lysis Buffer B** (red cap) to the pellet and add 10 μ L of **Lysis Reagent A** (red screw cap); mix and incubate for 15 minutes at 37°C in a thermoshaker
- 5. Add 50 μ L of **Lysis Buffer C** (red cap) and incubate for 15 minutes at 37°C and 1,400 rpm in the thermoshaker
- 6. Centrifuge for 10 min at 13,000 rpm
- 7. Pipette 45 µL of **Test Solution D** (yellow cap) per sample (including the standards) into the wells of a microplate and pre-incubate for at least 5 minutes at 50°C and 500 rpm in the thermoshaker
- 8. Add 10 μ L of the supernatant from step 6 to each well (row A1–H1 is reserved for the respective standards); cover the microwell plate with a lid and incubate for 10 min at 50°C and 500 rpm in the thermoshaker
- 9. Transfer 50 μ L of reaction mixes to the binding plate and incubate for 10 min at 50°C and 500 rpm in a thermoshaker
- 10. Discard all liquid and wash the plate with 200 μ L **Washing Solution E** (blue cap), discard Washing Solution
- 11. Dilute a suitable amount of **Enzyme Solution F** (green screw cap) $\underline{1:100}$ with **Washing Solution E** (blue cap) and add 100 μ L of the mixture to each well of the microplate; cover the plate with a lid and incubate for 10 minutes at 25°C and 500 rpm in the thermoshaker
- 12. Discard all liquid and add 200 µL of **Washing Solution E** (blue cap) to each well and incubate for 1 minute at room temperature and 500 rpm in the thermoshaker; repeat the washing step once
- 13. Discard all liquid and add 100 μL **Substrate Solution G** (green cap) per sample to the wells of the microplate; cover the plate with a lid and incubate for 5-15 minutes at 25°C and 500 rpm in the thermoshaker
- 14. Add 50 μL **Stop Solution H** (green cap) to each well
- 15. Place the microplate in a microplate reader and measure the optical density in each well at 450 nm; perform data analysis





Overview of the HybriScan®D Beer procedure:



1. Sample filtration (100 – 1000 ml, 15 min)



2. Enrichment culture, optional (NBB-media, 24h)



3. Cell lysis (2 ml sample, 13.000 rpm; 37°C, 45-60 min)



4. HybriScar® test solution (Forming of "sandwich complexes" between specific probes an the sample , 15 min)



5. Immobilisation (Binding of the "sandwiches" to the binding plate, 15 min)




6. Enzyme coupling (Coupling of an enzyme to the "sandwiches", 25 min)



7. Washing (Removal of unbound components, 2x 1 min)



8. Signal read out (Colour reaction and signal read out, 20 min)



9. Test analysis (Analysis of thedata, 10 min)

Advantages

- Simultaneous detection of all beer spoiling bacteria in only one test
- Specific for living cells
- Quantification of cells by standards
- Time saving of 3 to 5 days in comparison to cultivation based assays
- Minimial sample preparation
- High sample throughput using 96-well microplates